SENNHEISER

SENNHEISER ELECTRONIC KG, D-3002 WEDEMARK

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User's Guide

Mode d'Emploi

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RF-Receiver

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
fréquence

EK 1013-6

Bedienelemente
Controls
Éléments de service

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EK 1013-6 RF RECEIVER

Description

The EK 1013 is a single-channel, narrow band, battery-powered RF receiver which is used for wireless audio transmission in combination with the SK 1013-6 Mikroport transmitter.

Both, the circuitry and the housing of the unit, have been designed to meet the needs of people with hearing problems who need an aid to make communication possible, even under difficult acoustic conditions. One application for the system is as an instrument with which one can help integrate such people into all kinds of schools and universities. Equipped with a Sennheiser Mikroport system the user no longer needs to strain to hear the voice of the speaker, distorted and diminished by environmental sounds and the distance between listener and podium; instead, the words spoken are received directly, eliminating the factors of the acoustics of the room and of distance to the speaker. The best way of explaining how the system works is with the help of an example:

The speaker (e.g. a teacher) wears the transmitter, which is equipped with a microphone. The listener, a child or an adult with hearing problems, wears the EK 1013-6 receiver. Depending on the level of the hearing handicap the output of the receiver can be connected either to a pair of headphones or to the listener's own hearing aid. The EK 1013-6 receives the signal and amplifies it so that the listener can hear it loudly and clearly. Even if the listener is sitting in the very back row of a school auditorium with poor acoustics he or she can hear the teacher as clearly as if he were standing next to them and speaking to them directly. The receiver is also fitted with a microphone in order to make it possible for the user to hear the answers of his fellow students to the teacher's questions clearly, as well as the teacher's voice. A special circuit ensures that the teacher's voice receives priority (i.e., it is amplified) during dialogs between students and the teacher. While the teacher is speaking the sounds in the classroom are muted by approximately 20 dB.

And now here is a little information about the technology of the receiver: The EK 1013-6 is a single-channel system, and it can be set to receive a chosen frequency in the range between 30 and 45 MHz (see the section on "Changing Frequencies" for details of the frequencies). The output of the receiver can be used for the connection of headphones, stethoscope headphones, hearing aids, induction plates or the EZT 1011 induction loop. The output voltage/volume is adjustable. Extremely high transmission quality is

guaranteed by the built-in "HiDyn" compander system and an integrated squelch.

The EK 1013-6 is powered by three batteries, which can be either non-rechargeable or rechargeable. Two charging contacts are located on the side of the unit - when rechargeable batteries are being used they can be recharged by connecting the SZL 1011-2 charging unit to these contacts. Two carrying straps are supplied with the unit, a neck strap and a chest strap.

Operation

Inserting the Batteries

Fig. 1:
Opening the Battery Compartment

Fig. 2:
Inserting the Batteries

The transmitter can be operated either with alkali-manganese dry batteries (IEC LR 6) or with rechargeable NiCd batteries (1.2 V/500 mAh sinter accumulators). In the case of dry batteries we strongly advise that you use the alkali-manganese type, as the performance of the transmitter cannot be guaranteed when other types of dry batteries are used. When the transmitter is used frequently, however, it is more economical to use rechargeable NiCd batteries. These batteries can be recharged with the SZL 1011-2 charging unit. For this purpose the EK 1013-6 is fitted with two contacts on the side of its housing, which make it possible to charge the batteries without removing them from the unit. Depending on the types of batteries used the transmitter can be operated for different lengths of time before changing or recharging the batteries:

Rechargeable NiCd batteries: approx. 10 hours

Alkali-manganese batteries: approx. 30 hours

Please note when using rechargeable batteries:

1. Switch the receiver off before charging.
2. Replace the receiver in the charging slot of the charging unit after each use. This makes sure that the batteries are always fully charged.

Connecting Hearing Aids with Audio Inputs

One or two hearing aids equipped with audio inputs can be connected to socket **3**. It is essential that the connection of the units be carried out by a hearing aid specialist. The connection cable required for HdO units with Euro inputs is available as an optional accessory.


Fig. 3: Wiring of the AF Output Socket (switch setting: "3 mV")

Once the connection has been made shift switch **4**, which is located under the sliding cover **8** to the "3 mV" position. Fig. 8 shows the removal of the sliding cover.

Important: Replace the sliding cover once you have changed the switch setting!

Set the desired operating mode with switch **1**.

Operating Mode $\Psi + \text{ID}$

When the switch is in this position the audio signals of the hearing aid and the receiver are mixed. The volume of the receiver signal can be adjusted with knob **6**. The receiver's built-in microphone **2** can also be mixed in with the help of knob **5** making it possible for the user to participate actively in classroom activity.

Operating Mode Ψ / ID

When the switch is in this position an automatic fade-in system is activated which gives priority to the voice of the teacher. While the teacher is speaking the sensitivity of the hearing aid and of the built-in microphone are reduced by approximately 20 dB.

Since the sensitivity reduction of the hearing aid is dependent upon the input impedance of the audio input a hearing aid expert must carry out an adjustment. This adjustment is carried out with a potentiometer which is accessible via an opening on the back of the unit (Fig. 4).

Adjustment Procedure

1. Connect the hearing aid to the EK 1013-6 (see socket wiring diagram, Fig. 3).
 - 1.1 Remove the seal.
2. Switch on the receiver and the transmitter. Take care that the transmitter and receiver frequencies correspond.
3. Speak into the transmitter microphone.
4. Use control **6** to adjust the signal from the transmitter to that of the hearing aid.
5. Speak into the transmitter microphone. At the same time adjust the potentiometer of the EK 1013-6 (Fig. 4) so that the sensitivity of the microphone built into the hearing aid is significantly reduced.
6. Close the opening in the housing with the supplied seal.

Fig. 4

Connecting Hearing Aids without an Audio Input

In the case of hearing aids which do not have an audio input it is possible to transfer the signal to the hearing aid by induction, i.e. with the help of an induction loop (Model EZT 1011, Item No. 1954) or with an induction plate (Model EZT 110, Item No. 2647). Please set switch **4** to the "1.5 V" position when using one of these units.


Fig. 5: Wiring of the AF Output Socket (switch setting: "1.5" V)

The induction loop can be connected directly. In the case of the EZI 110 a connection cable is required, either the HZL 30-6 (for one earpiece) or the HZL 32-6 (for two earpieces). This socket can also be used to drive headphones, so long as their impedance is at least 47Ω . The output voltage, and thus also the volume, can be adjusted with knob **6**. The different operating modes have no effect upon the functioning of the system; switch **1** only has the function of an On/Off switch, and can thus be either in the $\Psi + \square$ or the Ψ / \square position.

Attaching the Neck and Chest Straps

Two straps are supplied with the unit, a neck strap and a chest strap, for attaching the unit to your body. Please slot the fittings on the ends of these straps into the slots in the unit, as shown in Figs. 6 and 7. There should be an audible click as the fittings slot in.

Please take care that the antenna is not pressed against your body and that it does not become entangled in the cable of the HdO unit or the headphones when you are using the receiver, since this has an adverse effect on the reception of the EK 1013-6.

Fig. 6


Fig. 7

Changing Frequencies

It can sometimes be necessary to change the operating frequency, for example when two wireless systems with the same frequency being operated in the same school interfere with one another. In order to enable the operation of two or more systems concurrently in one school it has been made possible to change the frequency of the EK 1013-6. The sliding cover **8** contains the component which determines the operating frequency of the receiver. The EK 1013-6 can be set for reception on another frequency by exchanging this sliding cover (Fig. 8). Your dealer can provide you with information regarding the licenced frequencies in your country (please see the enclosed Service Partner Card).

Fig. 8

Operating permit

The use of wireless microphone systems is generally subject to local licencing or permit arrangements. For further details contact your local Sennheiser distributor.

Technical Specifications

Receiver

Receiving Frequency	1 frequency between 30 and 45 MHz
Rated deviation	± 8 kHz
Sensitivity for 52 dB S/N	≤ 1.5 µV
THD at 1 kHz and U _{IF} = 1 mV	≤ 3 %
Frequency range	40 Hz - 12 kHz
Electronic squelch (internally adjustable)	Factory setting approx. 1 µV

AF output Stage

Output level, switchable	max. 2 V or 6 mV
Internal resistance	
Level switch in "1.5 V" pos.	approx. 1 Ω
Level switch in "3 mV" pos. and Ψ / \square operating mode	approx. 1.5 MΩ
Level switch in "3 mV" pos. and Ψ / \square operating mode	adjustable from 100 Ω to 2.3 kΩ
Minimum load impedance	50 Ω
Output power	50 mW
Output socket	3.5 mm/phone socket

General Specifications

Power supply	3 dry batteries (IEC R 6), or 3 1.2 V rechargeable NiCd batteries
Power consumption at 4.5 V rated voltage	approx. 40 mA
Operating time	approx. 30 hours with alkali- manganese dry batteries approx. 10 hours with rechargeable NiCd batteries
Dimensions in mm (without antenna)	100 x 66 x 37
Weight with batteries	approx. 180 g.
Supplied equipment	1 receiver 1 chest strap (EZU 1011) 1 neck strap (EZU 1011)

The manufacturer reserves the right to make alterations, especially with regard to technical progress.

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