SENNHEISER

SENNHEISER ELECTRONIC KG., D-3002 WEDEMARK

Bedienungsanleitung
User's Guide
Mode d'Emploi

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SK 1013-6 RF-TRANSMITTER

Description

The SK 1013-6 is a battery-powered RF transmitter used for wireless audio transmission in combination with an RF receiver, ideally with the EK 1013-6. One of the possible applications for the SK 1013-6 is as an aid for people with hearing problems in all kinds of schools and universities. Equipped with a Sennheiser Mikroport system consisting of the SK 1013-6 and the EK 1013-6 the user no longer needs to strain to hear the voice of the speaker, distorted and diminished by environmental sounds and the distance between listener and podium; instead, the words spoken are received directly, eliminating the factors of the acoustics of the room and of distance to the speaker. The best way of explaining how the system works is with the help of an example:

The speaker (e.g. a teacher) wears the transmitter, which is equipped with a microphone. The listener, a child or an adult with hearing problems, wears the EK 1013-6 receiver. Depending on the level of the hearing handicap the output of the receiver can be connected either to a pair of headphones or to the listener's own hearing aid. The EK 1013-6 receives the signal and amplifies it so that the listener can hear it loudly and clearly. Even if the listener is sitting in the very back row of a school auditorium with poor acoustics he or she can hear the teacher as clearly as if he were standing next to them and speaking to them directly.

Here is a little information about the technology of this transmitter:

The SK 1013-6 is a narrow band transmitter and is preset to a specific frequency in the range between 30 and 45 MHz. The transmitter's 8-pin input socket can be used to connect microphones, for instance the MKE 4012 or MKE 2012 plug-in electret microphones, and audio signals from tape decks or cassette recorders can be input via the 3.5 mm phono socket. The built-in "HiDyn" compander system ensures extremely high transmission quality.

The SK 1013-6 is powered by three batteries, which can be either rechargeable or non-rechargeable. Two charging contacts are located on the side of the unit when rechargeable batteries are being used they can be recharged by connecting the SZL 1011-2 charging unit to these contacts.

The standard equipment provided with the transmitter consists of a neck strap and the GZK 1011 carrying clip.

Operation

Inserting the Batteries

Fig. 1: Opening the battery compartment cover


Fig. 2: Inserting the batteries

The transmitter can be operated either with alkali-manganese dry batteries (IEC LR 6) or with rechargeable NiCd batteries (1.2 V/500 mAh sinter accumulators). In the case of dry batteries we strongly advise that you use the alkali-manganese type, as the performance of the transmitter cannot be guaranteed when other types of dry batteries are used. When the transmitter is used frequently, however, it is more economical to use rechargeable NiCd batteries. These batteries can be recharged with the SZL 1011-2 charging unit. For this purpose the SK 1013-6 is fitted with two contacts **6** on the side of its housing which make it possible to charge the batteries without removing them from the unit. Depending on the types of batteries used the transmitter can be operated for different lengths of time before changing or recharging the batteries:

Rechargeable NiCd batteries: approx. 10 hours

Alkali-manganese batteries: approx. 30 hours

Please Note when Using Rechargeable Batteries:

1. Switch the transmitter off before charging.
2. Replace the transmitter in the charging unit after each use. This makes sure that the batteries are always fully charged.

Connecting a Microphone

You can connect either the MKE 4012 (cardioid characteristic) or the MKE 2012 (omnidirectional characteristic) to the transmitter's microphone input socket. Of these two plug-in microphones the MKE 4012 is more strongly recommended, as its excellent directional efficiency increases the intelligibility of speech in noisy surroundings. It is recommended that the windscreen supplied with the microphone should always be used, in order to prevent wind and pop noises.

Fig. 3: Socket wiring

Connecting Audio Equipment

Should you wish to input audio signals from a tape deck or cassette recorder to the transmitter, in addition to the spoken word, then you can connect these units to the 3.5 mm phono socket **2**. If a microphone is also connected it is switched off automatically when you insert a jack into the phono socket.

If possible use the record/playback socket when you are connecting your cassette recorder or tape deck to the transmitter. If only a headphones or loudspeaker socket is available, then please note that the phono input is only designed for audio signal voltages between 75 mV and 1.5 V. If this voltage is exceeded for a brief period the overload is levelled out by a high-quality limiting amplifier. Even so, such overloads should be kept to an absolute minimum. You can reduce the input level by turning down the volume control on your cassette recorder or tape deck.

Adjusting the Sensitivity of the Microphone and Phono Inputs

An adjustment knob **4** with which you can adjust the sensitivity of the microphone or phono input to the level of the input signal is located under

the cover **5** (see Fig. 4). In addition to the manual adjustment option the transmitter is also equipped with an automatic level control system, which automatically reduces the input sensitivity when the rated deviation (the value at which the frequency deviation of the transmitter reaches ± 8 kHz) is exceeded, thus preventing overloading.

Please proceed as follows in order to carry out the manual adjustment:

Microphone Input:

1. Connect a pair of headphones to the AF output socket of the receiver and turn the volume control knob of the receiver ("Ψ") all the way to the right.
2. Switch on the transmitter and the receiver. Take care that the transmission and reception frequencies correspond!
3. Turn the adjustment knob **4** all the way to the left.
4. While you are adjusting the level hold the transmitter in the position in which it is to be used.
5. Speak into the microphone at a normal volume. Experience has shown that a long drawn-out French "uuuuu" is particularly effective for level adjustment.

While speaking into the microphone listen to your voice in the headphones. Increase the sensitivity gradually until you can no longer hear any further increase in volume, then shift the adjustment knob **4** one quarter of a turn back to the left.


Fig. 4 a


Fig. 4 b

Phono Input:

Connecting the unit to the record/playback socket of the tape deck or cassette recorder:

1. See points 1 through 3 under "Microphone Input" above.
2. Switch the tape deck or cassette recorder to "play".

Listen to the playback on the headphones. Increase the sensitivity gradually until you can no longer hear an increase in the volume, then shift the adjustment knob **4** one quarter of a turn back to the left.

Sensitivity Adjustment when the Unit is Connected to the Headphones or Loudspeaker Output of the Cassette Recorder or Tape Deck:

1. See points 1 through 3 under "Microphone Input" above.
2. Switch the tape deck or cassette recorder to "play".
3. Turn the adjustment knob **4** to the central position.
4. Listen to the playback with the headphones.

Increase the input signal gradually by turning up the volume control of your tape deck or cassette recorder until you can no longer hear an increase in volume.

Attaching the Neck Strap and the Carrying Clip

A neck strap and a carrying clip for attaching the transmitter to your belt or clothing are supplied with the unit. Please slot the fittings at the ends of the strap into the slots located at the top of the unit, as shown in Fig. 5. There should be an audible click as the strap fittings slot in.

Fig. 5: Attaching the strap

Fig. 6: Attaching the clip


Fig. 7:
Carrying the unit with the Clip

The carrying clip supplied with the transmitter can easily and quickly be attached to your belt or waistband (Fig. 7). When carrying the transmitter in this manner a clip-on microphone such as the MKE 2 R, the MKE 10 R or the MKE 40 R should be used, since the MKE 2012 and MKE 4012 plug-in microphones are not suitable for this purpose.

Please take care that the antenna is not pressed against your body and that it does not become entangled in the cable of the clip-on microphone, since this reduces the operating range of the unit.

Operating permit

The use of wireless microphone systems is generally subject to local licencing or permit arrangements. For further details contact your local Sennheiser distributor.

Technical Specifications

Carrier frequency	1 frequency between 30-45 MHz
Frequency stability at temperatures between -10°C and +55°C and an operating voltage deviation of between +10% and -30%	± 2 kHz
RF output to dummy antenna	10 mW
Radiated interference power	≤ 4 mW
Modulation	FM
Rated deviation	± 8 kHz
Signal to noise ratio, relative to rated deviation (at minimum sensitivity)	68 dB (CCIR 468, peak, typical) 82 dB (DIN 45 500, curve A, eff., typical)
Microphone input	
AF sensitivity (rated deviation)	max. 3.5 mV (Limiter action) min. 70 mV (Limiter action)
Socket wiring	Contact 1: AF Contacts 2 & 3: earth Contact 8: + U _B for electret microphones
Phono Input	
Input impedance	approx. 220 kΩ
AF sensitivity (rated deviation)	max. approx. 70 mV (Limiter action) min. approx. 1.5 V (Limiter action)
Sensitivity adjustment range	approx. 26 dB
Control range of limiting amp. at maximum sensitivity	approx. 40 dB
minimum sensitivity	approx. 20 dB
THD at 1 kHz	< 1%
AF frequency response	40 Hz to 12 kHz
Roll-off above the upper frequency limit	18 dB/oct. (active low-pass filter)
Premphasis	50 μs
Power consumption at 4.5 V	approx. 35 mA
Power supply	3 x 1.5 V batteries (IEC LR 6)
Operating time	approx. 30 hours with alkali-manganese batteries approx. 10 hours with rechargeable NiCd batteries
Dimensions in mm (without carrying accessories)	height 100, 185 with helical antenna width 66, depth 38
Weight (with batteries and plug-in microphone)	approx. 200 g
Supplied equipment	1 transmitter, 1 neck strap (EZO 1011) 1 clip (GZK 1011)
The manufacturer reserves the right to make alterations, especially with regard to technical progress.	

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Warranty conditions

A warranty of 24 months effective from the date of purchase is offered to you on this equipment including the accessories contained in the package, except for the batteries and rechargeable accus. The warranty for accus and accu plugs is 12 month from the date of purchase. The warranty covers – at the discretion of Sennheiser distributor – a free replacement of the parts needed to rectify and fault which may give reason for complaint or a free repair, whichever is considered necessary. To save time, please send your Sennheiser set, all freight charges prepared, directly to the Sennheiser distributor of your country to be found in the enclosed list, include the warranty card and the voucher. Warranty claims due to misuse are excluded. This also applies to ear cushions as they are subject to wear and tear. **Applicable to the UK only:** In the United Kingdom this warranty is in addition to the consumer's statutory rights and does not affect those rights and does not affect those rights in any way.

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Warranty conditions

Applicable to the US. only

Sennheiser Electronic Corporation (N.Y.) warrants that each SENNHEISER product is free from defects in material and workmanship and agrees to repair or replace any part of the product which under normal installation and use discloses such defect, provided the product is delivered through the dealer or by the purchaser to Sennheiser Electronic Corporation (N.Y.), intact, for examination with all transportation charges prepaid to Sennheiser Electronic Corporation (N.Y.) within two years from the date of delivery by the authorized Sennheiser dealer to the purchaser and provided that such examination discloses in the judgment of Sennheiser Electronic Corporation (N.Y.) that it is thus defective.

This warranty does not extend to any parts which have been subjected to misuse, abuse neglect, accident, incorrect wiring, improper installation, or use in violation of instructions furnished by us nor to any Sennheiser product which may have been tampered with, altered or repaired outside of our New York Service Department or our home plant. This warranty is in lieu of all other warranties expressed or implied, and no one is authorized to assume any liability on behalf of Sennheiser Electronic Corporation (N.Y.), or impose any obligation on it in connection with the sale of any product other than as outlined above.

This warranty is conditioned on and becomes effective only after the purchaser returns the detachable card properly filled out, to Sennheiser Electronic Corporation (N.Y.), Service Department, 48 West 38th Street New York N.Y. 10018, within fourteen (14) days after date of delivery to purchaser from an authorized Sennheiser dealer.

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Vi preghiamo di prendere visione delle clausole di garanzia riportate sull'alegato certificato della Exhibo S.p.A.
(AB 10. 3. 87 TEL. 0 51 30/6 00-0)