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
## ELECTRONIC MEDICAL SCALE CLASSE III

Dear Customer,

Thank you for purchasing this product. We hope you will find it satisfactory.

Please read these instructions carefully to take full advantage of your product.

## READ BEFORE USE

The Electronic Medical Scale is designed and manufactured in a facility certified in the following Management Systems: ISO 9001 Quality, ISO14001 Environment and OHSAS18001 Health and Safety.

### APPLICATIONS

The electronic medical scale's intended use is for weighing patients for the purposes of monitoring, diagnosis and medical treatment.

### PRECAUTIONS OF USE

		This symbol indicates general precautions that should be taken when using this product.
	<b>WARNING</b>	This symbol indicates the possibility of serious injury if the product is mishandled or instructions are ignored.
	<b>CAUTION</b>	This symbol indicates the possibility of physical injury or equipment damage if instructions are ignored.

### WARNING

- Use the product on a firm and flat surface.
- Clean after use with a lightly dampened cloth without allowing any water to get into the device. Do not use solvents or immerse the product in water.
- Avoid excessive impact or vibration to the product, such as dropping the product onto the floor.
- Do not jump on the weighing platform.

### CAUTION

- Do not disassemble the scale as incorrect handling may cause injury.


## MAINTENANCE

- The product has been manufactured with precision. The user should not disassemble the product on its own.
- Remove batteries if the product is not in use for an extended period.

## ELECTRONIC MEDICAL SCALE CLASSE III

## PRODUCT SPECIFICATION

- Class I medical device (EC Directive 93/42 EEC).
- Class III OIML compliant.
- Class III Non-Automatic Weighing Instruments (EC Directive 90/384 EEC) compliant.
- Level indicator and adjustable feet.
- On/off function
- Capacity: 220 kg.
- Graduation: 100 g.
- Operates through 1.5 V X 6 (AA/ LR6) or with AC 9V [500 mA] adaptor (adaptor included).



## ITEMS INCLUDED IN THE PACKAGE

- 1 display
- 1 connection cable.
- 1 pillar.
- 1 Base part with platform and adjustable feet.
- 3 screws M5x15.
- 6 screws M5x45.
- 1 screw driver
- Instruction and assembly manual.
- Metrological book.

## ELECTRONIC MEDICAL SCALE CLASSE III

## AFFICHEUR


	<b>Weight \ Off</b>	<ul style="list-style-type: none"> <li>- Enter weighing mode</li> <li>- Power Off</li> </ul>
	<b>Mode \ Hold</b>	<ul style="list-style-type: none"> <li>- Choose mode : BMI, BFA, MEM</li> <li>- Hold weight display</li> </ul>
	<b>Height \ Tare \ Enter</b>	<ul style="list-style-type: none"> <li>- Enter height measurement mode</li> <li>- Automatic zero in weighing mode</li> <li>- Confirm settings in another mode</li> </ul>
	<b>Digits "0" to "9" (on full keyboard)</b>	<ul style="list-style-type: none"> <li>- Enter digits</li> </ul>
	<b>Up \ Print \ Delete</b>	<ul style="list-style-type: none"> <li>- Previous</li> <li>- Print results</li> <li>- Correct input</li> </ul>
	<b>Down \ Memory</b>	<ul style="list-style-type: none"> <li>- Next</li> <li>- Store settings</li> <li>- Memory recall</li> </ul>

## ELECTRONIC MEDICAL SCALE CLASSE III

## PROCEDURE


## SCREEN INSTALLATION AND CONNECTION

1. Insert the connection cable into the pillar


2. Connect the pillar cable to the screen and screw.

3. Using the screwdriver, attach the pillar to the platform by screwing 6 screws M5x45 slightly beneath the platform.

## ELECTRONIC MEDICAL SCALE CLASSE III


4. Insert the screen into the pillar and attach with 2 M5x15 screws.

## SETTING THE BALANCE

- Before you use the scale for the first time or when the scale has been moved to another place, the balance of the scale must be set.


1. Do not step on the scale or/and remove all weight from the platform.
2. Place the scale on a hard and flat surface, use the adjustable feet to adjust the balance until the bubble inside the level indicator is centrally located in the circle.

## PRODUCT USAGE


## INITIALISATION

1. Insert 6pc AA battery into the battery compartment located at the back of the base, or power up the scale through connecting the DC adaptor with the power socket.
2. Remove the batteries if the scale has not been used for an extended period.

## ELECTRONIC MEDICAL SCALE CLASSE III

3. When the scale is powered up, the software version will be displayed for a few seconds (see below), and then all.


The scales will switch off automatically after 2 seconds.

**OPERATION**

Weighing only operation

1. Press the button (🔌/OFF) to power up the scale. The number 0 scrolls on the screen, then the display settles on "0.0".

2. Wait until screen displays "0.0" before stepping on the scale. Then, stand still while your weight is being computed. The screen will display your weight automatically.


3. Step off the scale, the display returns to « 0.0kg ».
4. Press the (🔌/OFF) button to turn the scale off immediately or wait 10 sec, then it will turn off automatically.

## ELECTRONIC MEDICAL SCALE CLASSE III

## WEIGHING WITH BMI FUNCTION AND USER MEMORY

1. Press the **[MODE/HOLD]** key once, after two control screens, the symbol BMI appears on the display screen:


Press **[← / Tare Height]** key to confirm selection.

2. Wait until the scales display «STEP».
3. Step onto the scales.

NB: If you wish to use ID numbers for each user, refer to the setup section.

If you have activated the ID mode, the screen displays «- - -», enter a 4-digit ID then press key **[← / Tare Height]** to confirm..

4. The height display flashes.

#### a - For models without electronic height meter

Input the height with the number keys.

#### b - For models with electronic height meter

- Check that the meter has been returned to its starting position before measuring height.
  - NB: If the height meter has not been returned to its starting position before use, the LCD screen displays «|n|t».
  - Return the height meter to its starting position.
- Extend the height meter and adjust it to the person's height. The person's height will be displayed on the screen. The user can also press key **[▲/Print-CE]** to input data via a keypad rather than using height meter data.

#### c - Validation of the height and recording of user's personal data


- ⇒ If you do not wish to store the user's data
  - Confirm the height immediately by pressing key **[← / Tare Height]**.
- ⇒ If you wish to store the user's personal data

## ELECTRONIC MEDICAL SCALE CLASSE III


Before confirming the height, press and hold key [▼\ **Memory**] for 3 seconds, the screen then displays «M-00».

- On the keypad, enter the selected memory number and confirm by pressing [← / **Tare Height**]
- If the selected number has not been used before, the screen will display ----- . Press key [← / **Tare Height**] to confirm.
- If the number selected for storing the data has been used already, the screen will display USED.
  - \* To delete the previous record, press key [← / **Tare Height**], the screen displays CLEAR, press key [← / **Tare Height**] again to confirm the record number and the height.
  - \* If you do not wish to delete the previous record, select another number to store the data, press key [← / **Tare Height**], the screen displays -----, press key [← / **Tare Height**] again to confirm the record number and the height.

5. After confirming the height, stay still on the scales and wait for the weight to stabilise: STA.

6. The display will toggle between Height, Weight and BMI. Each result will be displayed for more than 4 seconds. To speed up the display, use the [▲] or [▼] keys.

7. To print the results, press and hold key [▲/Print-CE] for 3 seconds.

NB: If the user presses any key, after displaying the results, the unit will exit the current mode.

### RECALLING THE USER'S PERSONAL STORED DATA

1) Press key [MODE/HOLD]. Two control screens appear followed by the BMI screen.

2) Press key [MODE/HOLD] a second time, the screen then displays «BFA»:


3) Press key [MODE/HOLD] a third time, the screen then displays «MEM»:

## ELECTRONIC MEDICAL SCALE CLASSE III


- 4) Confirm the selection by pressing [  $\leftarrow$  / Tare Height].
- 5) Use the arrow or number keys to select the memory to be recalled and confirm the selection by pressing key [  $\leftarrow$  / Tare Height].

- 6) The user ID, if previously recorded, and the height will appear on the screen.  
Confirm the memory selection again by pressing key [  $\leftarrow$  / Tare Height].
- 7) The screen then scrolls «0», then «0.0 Kg».
- 8) Repeat from point 5 of the section above: Validation of height and storing users' personal data.

**RECALLING RESULTS FROM THE LAST RECORDED WEIGHING**

- 1) Press key [  $\nabla$  \ Memory] to enter the «recall results» mode.
- 2) Press on the arrow or number keys to select the memory to be displayed on the screen.  
Confirm the selection by pressing key [  $\leftarrow$  / Tare Height].

- 3) The values are displayed in the following sequence: ID number (if any), height, weight and BMI.  
If the recalled memory is empty, the screen will display «- - -».
- 4) To print the results, press and hold key [  $\blacktriangle$  / Print-CE] for 3 seconds.

**WEIGHING WITH IMPEDANCE (BODY FAT) ANALYSIS WITH USER MEMORY**

- 1) Press key [MODE/HOLD] once, after two control screens, the symbol BMI appears on the display:

## ELECTRONIC MEDICAL SCALE CLASSE III

2) Press key **[MODE/HOLD]** again, the screen will display:

Press key **[← / Tare Height]** to confirm selection.


3) Wait for the scales to display «STEP».

4) Step onto the scales with bare feet.

NB: If you wish to use ID numbers for each user, refer to the setup section.


If you have activated the ID mode, the screen will display «- - -», enter a 4-digit ID, then press key **[← / Tare Height]** to confirm.

5) Enter your personal data. Firstly, the screen displays: AGE


⇒ Enter age and confirm by pressing key **[← / Tare Height]**

6) Then, the screen displays: 

⇒ Select the gender of the user using the arrow keys and confirm by pressing key **[← / Tare Height]**

7) Finally, the screen displays: **N**

⇒ Select the user mode using the arrows:  
N = non-athlete - A = athlete. Confirm by pressing key **[← / Tare Height]**

8) Enter or measure the height

**a - For models without electronic height meter**

Input the height with the number keys.

## ELECTRONIC MEDICAL SCALE CLASSE III

**b - For models with electronic height meter**

- Check that the meter has been returned to its starting position before measuring height.

NB: If the height meter has not been returned to its starting position before use, the LCD screen displays «InIt». Return the height meter to its starting position.

- Extend the height meter and adjust it to the person's height. The person's height will be displayed on the screen. The user may also press key [▲/Print-CE] to enter data via a keypad rather than using height meter data.

**c - Validation of the height and recording of user's personal data**


⇒ If you do not wish to store the user's data. Confirm the height immediately by pressing key [←/Tare Height]

⇒ If you wish to store the user's personal data


Before confirming the height, press and hold key [▼\Memory] for 3 seconds, the screen will then display «M-00».

- On the keypad, enter the selected memory number and confirm by pressing [←/Tare Height]
- If the selected number has not been used before, the screen will display «-----». Press key [←/Tare Height] to confirm.
- If the number selected for storing the data has been used already, the screen will display USED.
  - \* To delete the previous record, press key [←/Tare Height], the screen displays CLEAR, press key [←/Tare Height], again to confirm the record number and height.
  - \* If you do not wish to delete the previous record, select another number to store the data, press key [←/Tare Height], the screen displays «-----», press key again to confirm the record number and height.

9) After validating the height, stay still on the scales and wait for the weight to stabilise: STA

10) Stay still on the scales, the body fat calculation is in progress, the screen displays:


11) A first screen appears showing personal data:


## ELECTRONIC MEDICAL SCALE CLASSE III

12) Scroll through the results. Each result will be displayed for more than 4 seconds. To speed up the display, use the [▲] or [▼] keys.


- Screen 1: Weight

- Screen 2: Body fat percentage

- Screen 3: Fat mass


- Screen 4: Fat-free mass

- Screen 5: BODY WATER percentage

- Screen 6: Muscle mass

- Screen 7: Bone mass

13) To print the results, press and hold key [▲/Print-CE] for 3 seconds.

NB: If the user presses any key, after displaying the data, the unit will exit the current mode.

## ELECTRONIC MEDICAL SCALE CLASSE III


**Recalling the user's personal stored data**

- 1) Press key **[MODE/HOLD]**. Two control screens appear followed by the BMI screen.
- 2) Press key **[MODE/HOLD]** again, the BFA screen is displayed.
- 3) Press key **[MODE/HOLD]** a third time, the MEM screen is displayed.

Refer to the above section: Validation of the height and storing the user's personal data

**Recalling results from the last recorded weighing**

- 1) Press key **[▼ Memory]** to enter recall results mode.
- 2) Press on the arrow or number keys to select the memory to be displayed on the screen.  
Confirm the selection by pressing key **[← / Tare Height]**.

- 3) The values are displayed in the following sequence: ID no. (if any), Height, Weight, Age, Gender and mode:
- 4) Secondly, the results will appear.  
If the recalled memory is empty, the LCD screen displays «-----».
- 5) To print the results, press and hold key **[▲/Print-CE]** for 3 seconds.

**SETUP MODE****DEACTIVATION OF THE PRINT MODE**

In the factory setup, the print mode is activated. If you wish to deactivate it:

- Press and hold **[MODE]** for more than 10 seconds to go into the setup mode.
- After a control screen, the screen displays «Print».
- Press the **[▲]** or **[▼]** keys to select «OFF» then press key **[← / Tare Height]** to confirm.

**ACTIVATION OF THE ID MODE**

In the factory setup, the ID mode is deactivated. If you wish to activate it:

- Press and hold key **[MODE]** for more than 10 seconds to go into the setup mode.
- Go from the PRINT function setup by pressing key.
- The screen will display «Id no». Press the **[▲]** or **[▼]** keys to select «ON» and press key **[← / Tare Height]** to confirm.

## ELECTRONIC MEDICAL SCALE CLASSE III

**PERSONALISATION of the body fat analysis results display (BFA)**

In the factory setup, the following data is displayed: body fat, body water mass, bone mass and muscular mass. You may choose to remove particular data from the display.

- Press and hold key **[MODE]** for more than 10 seconds to go into the setup mode.
- Go from the PRINT function and ID NO setup by pressing [**←**] / **Tare Height**].
- The screen then displays «BFA». For each data, use keys [**▲**] or [**▼**] to select ON or OFF and thus deactivate or activate its display. Confirm your selection by pressing key [**←**] / **Tare Height**].

**PRINT FORMAT FROM THE BMI BALANCE**

MEASUREMENT RESULTS			
<b>ID</b>	0001	<b>Body fat</b>	10 %
<b>Age:</b>	30 yrs	<b>Body fat</b>	5,0 kg
<b>Gender:</b>	Female	<b>Fat-free mass</b>	60,0 kg
<b>Activity</b>	Normal	<b>Body water mass</b>	50 %
<b>Height</b>	175,0 cm	<b>Bone mass</b>	4,5 kg
<b>Weight</b>	68,0 kg	<b>Muscle mass</b>	20,5 kg
<b>BMI</b>	21,5		

**SPECIAL INDICATIONS**

Err0	Initial zero error. Get off the scales, wait for the display to show 0.0 and restart the measurement.	AdErr	Hardware problem. Contact your retailer.
Err1	Unstable weight: restart the weighing.	EEP	Hardware problem. Contact your retailer.
Err2	Overload: remove the weight immediately.		
Err3	Body impedance beyond the limit (100 ohm to 1kohm): restart the measurement. If Err3 appears often and depending on your skin type (thin or thick), it may be necessary to moisten your feet to ensure better contact with the electrodes.		
Err4	Body fat percentage out of range 4-50%: restart the measurement.		
Err5	Water volume out of range. Moisten your feet and restart.		
bAtt	Replace the battery.	Ustb	Unstable zero. Wait for the display to show 0.0 and restart the measurement

## ELECTRONIC MEDICAL SCALE CLASSE III

## PROTECT THE ENVIRONMENT

Dispose of used batteries in a designed container so that they can be collected and recycled. Do not mix different types of batteries. Do not mix new and used batteries.

After life, dispose of this product at a designated waste recycling point.

## GUARANTEE

This product is guaranteed against defects in materials or manufacturing. During the guarantee period, any such defects will be repaired free of charge (proof of purchase must be shown in the case of a complaint under guarantee). The guarantee does not cover damage resulting from improper use or negligence. In the event of a claim, first contact the shop where you purchased the product.

## AFTER SALES SERVICE

For calibration and after sales services, please contact Terrailon consumer service at 0 826 88 1789 or at the following address: [serviceconsommateurs@terraillon.fr](mailto:serviceconsommateurs@terraillon.fr)

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### **Terrailon UK Ltd**
### **Terrailon Deutschland GmbH**
### **Terrailon Asia Pacific Ltd**
### **Terrailon Corp**