TOPCOM[®]

HBM Watch 1010 Sport

**USER GUIDE
HANDLEIDING
MANUEL D'UTILISATEUR
BEDIENUNGSANLEITUNG
MANUAL DE USUARIO
BRUKSANVISNING
BRUGERVEJLEDNING
ANVÄNDARHANDBOK
KÄYTTÖOHJE**

V1.1

- UK** The features described in this manual are published with reservation to modifications.
- NL** De in deze handleiding beschreven mogelijkheden worden gepubliceerd onder voorbehoud van wijzigingen.
- FR** Les possibilités décrites dans ce manuel sont publiées sous réserve de modifications.
- DE** Die in dieser Bedienungsanleitung umschriebenen Möglichkeiten, werden vorbehaltlich Änderungen publiziert.
- ES** Las características descritas en este manual pueden ser objeto de futuras modificaciones.
- SE** Funktionerna i denna bruksanvisning publiceras med reservation för ändringar.
- DK** Vi forbeholder os retten til ændringer af de specifikationer, der er beskrevet i denne brugsanvisning.
- NO** Funksjoner beskrevet i denne manualen kan endres uten nærmere informasjon.
- FI** Tässä ohjekirjassa julkaissut tiedot voivat muuttua ilman ennakoilmoitusta.

This product is in compliance with the essential requirements and other relevant provisions of the R&TTE directive 1999/5/EC.
The Declaration of conformity can be found on :

<http://www.topcom.net/support/cedclarations.php>

1 INTRODUCTION

Congratulations on your purchase of the HBM Watch 1010. We hope you will enjoy using the Heart Meter and the benefits it offers. In exercise or sports, your heart naturally speeds up in pumping blood to the body in accordance with the increased energy level.

The HBM Watch 1010 features a heart measuring and data recording functions in addition to the normal time functions of a watch. The styling will allow you to wear it comfortably during exercise as well as everyday use.

2 SAFETY NOTES

Using a Heart Monitor is NOT a substitute for professional medical attention. If you suffer from heart disorder, illness or have recently taken up fitness training, it is advisable that you seek a doctor or talk to your trainer and work out a fitness plan that suits you.

The HBM Watch 1010 should only be used as a guide in maintaining a safe heart beat rate for exercises. It should NOT be relied upon as a medical heart rate monitor. If you are uncertain about how your exercise program will impact on your heart, please consult a physician.

To ensure your safety, please use the HBM Watch 1010 under a doctor or coach's direction if you have one of the following conditions:

1. Cardiopulmonary disease
2. Obesity
3. No exercise for an extended period of time

3 CARE AND MAINTENANCE

Never use force to open the case or backside cover. DO NOT expose the watch to extreme weather conditions or cleaning agents and solvents. Avoid dropping or rough handling of the watch. DO NOT use the watch when diving or under water. Clean the watch with a dry soft cloth or a soft cloth moistened with water and mild detergent. Store the watch in a dry place when not in use.

Take care of you chest belt. Wash the chest belt with with water. Let it dry naturally, avoiding putting the chest belt under the environment of high temperature or touching the corrosive as strong acid or alkalis material.

4 LIMITS AND EXECISE ZONE

Exercise zones are established by setting Upper and Lower Heart Rate Limits. These limits constitute a certain percentage of your Maximum Heart Rate (MHR).

The following formula will help you to make an educated guess :

$$\text{MHR} = 220 - (\text{Your age})$$

The HBM Watch 1010 will automatically calculate 3 Target Zones based on your age and weight.

One extra target zone is freely programmable.

5 BATTERY INSTALLATION

If you are not confident with electronics, let change the battery by a professional. Do not open the watch yourself.


5.1 WATCH

Remove the back cover of the watch with a small screwdriver. Install the 3V(CR2016) battery with the positive (+) pole facing up and replace the cover. Remove the battery and reinsert it if the display shows irregular figures.


5.2 CHEST BELT

Remove the battery cover from the bottom of the belt using a small coin. Install the 3V battery with the positive (+) pole facing the battery cover and replace the cover.


6 WEARING THE BELT

Fasten the fastener at one end. Try the chest belt on your chest and loosen the stretch band if necessary.

Adjust the length of the stretch band until it feels conformable. The stretch band must cling to the chest. Fasten the fastener at the other end. Adjust the sensor to the center of your chest and be sure the backside of the sensor clings to your chest and touches the skin.


The chest belt has to be worn while using the heart rate function.

Moisten the chest belt conductive pads with saliva or ECG-gel (available at your local pharmacy) in order to ensure good contact with the skin at all time.

7 HOW TO USE THE HBM WATCH 1010

The HBM watch 1010 can be worn like a watch on your wrist or it can be fixed on the handle bar of your bike with the included optional mounting system. The maximum distance transmission should be less than 60cm.

8 PRECAUTIONS

- Avoid fastening the strap too tightly.
- Do not bend or fold the chest belt for preventing damage.
- Wetting the skin, where there is contact with the conductive area of the chest belt will improve the conduction and get a more stable signal.
- The physical condition of individual might effect the intensity of measured signal.
- Avoid using the heart rate monitor close to a trolley car, a tram stop, transformer, electric substation and high-tension distribution lines etc. The radio signal will be affected under the environment from high voltage and strong magnetic fields.

9 BUTTONS

A. MODE BUTTON

- Press to toggle between **Clock** and **Pulse** mode

B. SELECT BUTTON

- In **Pulse** mode, press to select the different heart rate monitor functions.
- In **Clock** mode, press to enter Alarm mode.
- In **Setting** mode, press to select the digits.

C. SET BUTTON

- In **Pulse** mode, hold to clear the date or to enter the setting mode
- In **Clock** mode, hold to enter setting mode or press to confirm the setting.

D. LIGHT BUTTON

- Press to turn on the back light.

10 GETTING STARTED

When you use the watch for the first time, the display is blank. Press a button to wake up, the default age '23' is displayed.


When the watch was already running you have to **RESET** the watch. Press and hold both buttons **B** and **C** until all segments are displayed.

10.1 SETTING THE EXERCISE ZONE

After wake up, the HBM Watch 1010 will automatically calculate 3 Target Zones after setting your age and weight. During heart rate measurement, these settings will also be used for the Calorie Counter and Fat Loss Counter.


1. Setting your age

- Press button **B** to adjust the age. Hold button **B** for fast advance.
- Press button **C** to confirm your setting and go to weight setting mode.


2. Setting your weight.

- Press and hold button **B** to select the unit of weight 'Kilogram (KG) or Pounds (LB)'.
- Press button **C** to confirm
- Press button **B** to adjust the weight. Hold button **B** for fast advance.
- Press button **C** to confirm your setting and go to exercise zone setting mode.

3. Exercise zone setting

After setting your age and weight, the HBM 1010 will automatically calculate exercise target zones as a percentage of your Maximum Heart Rate (MHR). (Zone 1: 50-65%, Zone 2: 65-80% and Zone 3: 80-95%) Example when your age is 30.

Zone 1 : 50% - 65% Zone 2 : 65% - 80% Zone 3 : 80% - 95%

Press button **B** to select the different zones. You will see the corresponding Upper and Lower limit.

If you want to use one of the Target Zone, press button **C** to confirm.

TOPCOM HBM Watch 1010 Sport

When you do not want to use an automatic defined Target Zone, you can freely program one extra Target Zone.


- Press button **B** to select the programmable Target Zone indicated with **U**.
- Press button **C** to confirm

Lower limit setting

- Press button **B** to adjust the digits of the lower limit. Hold button **B** for fast advance.
- Press button **C** to confirm

Upper limit setting

- Press button **B** to adjust the digits of the lower limit. Hold button **B** for fast advance.
- Press button **C** to confirm


After confirmation you will enter the **Clock mode**.


11 CLOCK SETTING

In **clock mode** :

- Hold button **C** to enter the **Clock Setting mode**.
- Press button **B** to select 12 hour or 24 hour format.
- Press button **C** to confirm.
- Press button **B** to adjust the digits and press button **C** to confirm each setting of the time and date as follow:

1. Hour
2. Minutes
3. Year
4. Month
5. Day

After confirming the day you will enter clock mode.

12 DAILY ALARM

12.1 DAILY ALARM SETTING

- In **Clock Mode**, press button **B** to enter **Alarm mode**.
- Hold button **C** to enter **Alarm Time Setting mode**. The hour will blink on the display.
- Press button **B** to adjust the digits and press button **C** to confirm.
- Repeat the same steps for setting the minutes.

12.2 DAILY ALARM ON/OFF


- In **Clock mode**, press button **B** to enter Alarm mode.
- In **Alarm mode**, press button **C** briefly to turn the alarm On (🔔) or Off.
- Press button **B** to go back to **Clock mode**.

13 PULSE MODE

13.1 PULSE MODE DISPLAY

1. Current Heart Rate
2. Average Heart Rate
3. Target Zone graphical bar
 - 3a. Upper Zone Range
 - 3b. In Zone Range
 - 3c. Below Zone Range
4. Upper/Lower Zone limit
5. Percentage of the Maximum Heart Rate. (See §4)


13.2 CURRENT HEART RATE

Make sure you wear the Heart Rate Belt as described in §6.

Press button **A** to enter **Pulse mode**. The current heart rate is displayed on the upper line of the display.


The  icon is flashing in the right corner of the display if a signal is received from the belt. The icon will freeze and the last recorded heart rate will be displayed if no signal is received.

For power saving, the HBM Watch 1010 will automatically return to clock mode if no signal is received.

13.3 % OF MAXIMUM HEART RATE

This is the data showing the percentage of the Current Heart Rate compared to the Maximum Heart Rate (MHR by age). It is useful to monitor your heart performance safely and effectively during exercise.

13.4 TARGET ZONE GRAPHICAL BAR

This graphical bar is the visual relationship between your Current Heart Rate and the Maximum Heart Rate.

The bar moves up and down from Below-Zone, In-Zone to Above-Zone (see §13.1). It is a very important and effective tool to monitor your workout both safely and effectively.


13.5 TARGET ZONE AUDIBLE ALARM

13.5.1 ALARM SOUND ON/OFF

In pulse mode:

Press button **B** until you enter **Target Zone mode (ZN)**.

Press button **C** to select alarm sound on () or off.

13.5.2 EXERCISE ZONE ALARM ACTIVATION

In pulse mode.

- Press button **B** until you enter the **Target Zone mode (ZN)**
- Hold button **C** until ∇/\blacktriangle flashes on the display.
- Press button **B** to select the target zone for the alarm.

\blacktriangle - The alarm will be activated only when your heart rate exceeds the upper limit.

∇ - The alarm will be activated only when your heart rate drops below the lower limit.

∇/\blacktriangle - The alarm will be activated both when your heart rate exceeds the upper or drops below the lower limit.

- Press button **C** to confirm.

13.5.3 EXERCISE ZONE SELECTION


Right after the Zone Alarm Activation.

- Press button **B** to select Zone 1, 2, 3 or U.
- Press button **C** to confirm.

13.6 AVERAGE HEART RATE

The average heart rate is displayed on the second row of the LCD.

13.7 AUTO EXERCISE TIMER (TMR)


In pulse mode:

Press button **B** until you enter **Timer mode (TMR)**.


The timer is activated automatically whenever the pulse signal from the chest belt is received.

Press button **C** to read the different times from your zone.


1. Below Zone - time


2. In Zone - time

3. Above Zone - time

4. Total exercise time


To reset the average heart rate and the timer. Press button **B** to enter the timer mode (TMR) and then hold button **C** to clear all data.

13.8 STOPWATCH (SW)

In pulse mode :


- Press button **B** until you enter the **Stopwatch mode (SW)**
- Press button **C** to start or to stop.
- Press and hold button **C** to reset.

13.9 MAXIMUM AND MINIMUM HEART RATE (MEM)

In pulse mode:


- Press button **B** until you enter the **Memory mode (MEM)**
- Press button **C** to check the maximum heart (HI) and minimum heart rate (LO)
- Press and hold button **C** to reset the data.

13.10 AUTO CALORIE COUNTER (CAL)

In pulse mode:


- Press button **B** to enter the **Calorie mode (CAL)**
- The approximate energy consumption in CALORIE, based on your age and weight, will be counted automatically whenever a pulse is received from your chestbelt.
- Press and hold button **C** to reset.

13.11 FAT LOSS COUNTER (FAT) IN GRAMS

In pulse mode:

- Press button **B** until you enter the **Fat Loss mode (FAT)**
- The approximate fat loss in gram, based on your age and weight, will be counted whenever your heart beat reaches 65% of your MHR.
- Press and hold button **C** to reset the data.

14 TOPCOM WARRANTY

14.1 WARRANTY PERIOD

The Topcom units have a 24-month warranty period. The warranty period starts on the day the new unit is purchased. Consumables or defects causing a negligible effect on operation or value of the equipment are not covered.

The warranty has to be proven by presentation of the original purchase receipt, on which the date of purchase and the unit-model are indicated.

14.2 WARRANTY HANDLING

A faulty unit needs to be returned to a Topcom service centre including a valid purchase note.

If the unit develops a fault during the warranty period, Topcom or its officially appointed service centre will repair any defects caused by material or manufacturing faults free of charge.

Topcom will at its discretion fulfil its warranty obligations by either repairing or exchanging the faulty units or parts of the faulty units. In case of replacement, colour and model can be different from the original purchased unit.

The initial purchase date shall determine the start of the warranty period. The warranty period is not extended if the unit is exchanged or repaired by Topcom or its appointed service centres.

14.3 WARRANTY EXCLUSIONS

Damage or defects caused by incorrect treatment or operation and damage resulting from use of non-original parts or accessories not recommended by Topcom are not covered by the warranty.

The warranty does not cover damage caused by outside factors, such as lightning, water and fire, nor any damage caused during transportation. No warranty can be claimed if the serial number on the units has been changed, removed or rendered illegible.

<p>The CE symbol indicates that the unit complies with the essential requirements of the R&TTE directive.</p>

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