# Instructions For Use


## **REVOLUTION<sup>®</sup>** **45MHz ROTATIONAL** **IMAGING CATHETER** **REF 89000**

<b>English</b>	<b>Page 1</b>
<b>Nederlands</b>	<b>Page 4</b>
<b>Français</b>	<b>Page 7</b>
<b>Deutsch</b>	<b>Page 10</b>
<b>Italiano</b>	<b>Page 13</b>
<b>Español</b>	<b>Page 16</b>

# REVOLUTION 45MHz ROTATIONAL IMAGING CATHETER

ENGLISH

## **CAUTION:**

1. U.S. Federal Law restricts this device to sale by or on the order of a physician.
2. Prior to use, read this entire Instructions For Use.

## **INTENDED USE:**

The Revolution<sup>®</sup> catheter is intended for the intravascular ultrasound examination of coronary arteries. Intravascular ultrasound imaging is indicated in patients who are candidates for transluminal interventional procedures.

## **DESCRIPTION:**

The Revolution 45MHz rotational IVUS imaging catheter consists of two main assemblies; the imaging core and the catheter body. The catheter body comprises of three sections; distal section w/ .014" compatible F/X port, proximal section (single lumen), telescope section.

The distal section and proximal (single lumen) sections comprise the "working length" of the catheter, the telescoping section remains outside of the guiding catheter. The telescoping shaft (section) allows the imaging core to be advanced and retracted for up to 150 mm of linear movement. The corresponding movement of the transducer occurs from the proximal end of the guidewire exit port to the proximal end of the window portion of the distal section.

The imaging core is composed of a hi-torque, flexible, rotating drive cable with a distal outward looking 45MHz ultrasonic transducer. An electromechanical connector interface at the proximal end makes the connection to the patient interface module (PIM). The PIM-catheter interface consists of an integrated mechanical drive assembly and electrical connection.

A flushing port with a one-way valve (Fig. 1) is used to displace the air initially present within the catheter. The catheter must be flushed with heparinized saline prior to use, as this provides the acoustic coupling media required for ultrasonic imaging. The one-way valve helps retain saline in the catheter during use.

The catheter body has a distal guidewire lumen with a proximal exit port located 2 cm from the distal end (Fig 1). A radiopaque (RO) marker is embedded in the catheter body at 0.5 cm from the tip. In addition, an insertion depth indicator is located on the catheter body at 100 cm, corresponding to femoral insertions.

The catheter is for use with the In-Vision Gold imaging system with software V5.0 or higher or the Volcano s5 and Volcano s5i imaging system. Consult your System Operator's Manual.

## **CONTRAINDICATIONS:**

This device is not currently indicated for use in cerebral or peripheral vessels. Use of IVUS Imaging Catheters is contraindicated where introduction of any catheter would constitute a threat to patient safety. Contraindications include: bacteremia or sepsis, major coagulation system abnormalities, patients disqualified for CABG surgery, patients disqualified for PTCA, severe hemodynamic instability or shock, patients diagnosed with coronary artery spasm, and total occlusion.

## **ADVERSE EFFECTS:**

Bleeding at the entry puncture site, injury to the vascular wall, thrombosis of the vessel, and peripheral embolization has occurred with the use of percutaneous intravascular catheter devices.

## **WARNINGS:**

- Use of the Revolution catheters is restricted to specialists who are familiar with, and have been trained to perform, the procedures for which this device is intended.
- DO NOT advance the catheter if resistance is encountered. The catheter should never be forcibly inserted into lumens narrower than the catheter body or forced through a tight stenosis.
- Care should be taken when utilizing devices that comprise a short monorail; in such instances advancement of the device distal to a deployed stent can result in exposure of the guidewire to the stent struts.
- In instances where the device has crossed a deployed stent, care should be taken when retracting the device to ensure that entanglement does not occur. Fluoroscopy should be used to monitor guidewire position with the respect to the imaging catheter and the stent; at no time should the imaging catheter be retracted if there is evidence of guidewire prolapse or if significant resistance to withdrawal is experienced. If either of these events occur, advance the imaging catheter distal of the stent and then carefully remove the whole system under the guidance of fluoroscopy.
- Care should be taken when re-advancing a guide wire after stent deployment. A guide wire may exit between stent struts when re-crossing a stent that is not fully apposed to the vessel wall. Subsequent advancement of the catheter could cause entanglement between the catheter and the stent. Care should be taken to slowly remove the catheter from a stented vessel.

## **PRECAUTIONS:**

The Revolution device is a delicate scientific instrument and should be treated as such. Always observe the following precautions:

- Contents supplied STERILE using an EtO (ethylene oxide) process. Do not use if sterile barrier is damaged. If damage is found call your Volcano Corporation representative.
- To maintain optimal patient safety, inspect the product prior to use. Do not use if saline leaks from any location other than the vent port in the monorail section.
- For single use only. Do not re-use, reprocess or re-sterilize. Re-use, reprocessing or re-sterilization may compromise the structural integrity of the device and/or lead to device failure which in turn may result in patient injury, illness or death.
- Re-use, reprocessing or re-sterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patients.
- The catheter has no user serviceable parts. Do not attempt to repair or to alter any component of the catheter assembly.
- Do not attempt to connect the catheter to electronic equipment other than the designated systems.
- Never attempt to attach or detach the catheter while the PIM motor is running. To do so may damage the connector.
- Avoid any sharp bends, pinching, or crushing of the catheter.
- Do not kink or sharply bend the catheter at any time. This can cause drive cable failure. An insertion angle greater than 45° is considered excessive.
- Turn the PIM "OFF" before withdrawing the imaging catheter.

## INSTRUCTIONS FOR USE:

### Materials and Equipment

Revolution Catheter  
Sterile PIM cover  
10 in. extension tubing  
3 cc and 10 cc syringes  
3-way stop cock  
Pre-formed guide catheter [0.064 in. (1.63 mm) I.D. min.] with Y-adapter assembly\*  
In-Vision Gold imaging system with software V5.0 or higher, Volcano s5 and Volcano s5i imaging system\*  
Heparinized, physiologic saline solution\*  
Guide wire, 0.014 in. (0.36 mm) max. diameter\*  
\*not packaged with catheter

### Inspection Prior to Use

Carefully inspect the package prior to use for any breach of the sterile barrier or damage to the contents. If the sterile barrier integrity is compromised or the contents damaged, contact your Volcano Corporation representative.

### Preparation For Use

Refer to the operator's manual or user's guide for instrumentation and PIM Setup.

Using sterile technique, remove the catheter from its sterile packaging. Remove the Packaging Coil protecting the catheter. Retract the movable Imaging Core completely to the proximal position via the telescope shaft.

Connect the 3 cc and 10 cc syringes to the 3-way stopcock, then connect the assembly to the extension tube and fill both syringes with heparinized saline. Ensure that all air is expelled from the system. Do not use if saline leaks from any location other than the vent port in the monorail section.

Connect the extension tube to the one-way valve on the catheter hub. The 10 cc syringe is to be used as a reservoir for refilling the 3 cc flushing syringe.

Flush the imaging catheter TWICE continuously with 3 cc volume each time. DO NOT USE EXCESSIVE PRESSURE. Advance the imaging core to its fully distal position; via the telescope shaft.

Connect the imaging catheter to the PIM by inserting the proximal end of the connector through the opening in the Sterile PIM cover, gently twisting the connector until it locks into place. To ensure that the hub is fully seated in the PIM, gently tug on the catheter hub.

Begin imaging by pressing the IMAGE button on the PIM long enough to ensure proper function of the catheter by observing a pattern of partial bright concentric rings on the monitor. Refill the 10 cc syringe as needed and reattach to the stopcock without introducing air into the line.

### Place Guide Catheter

Prepare the entry site with a sheath introducer according to the standard practice. Prior to inserting the imaging catheter, ensure the patient has been prepared using standard procedure for interventional treatment.

Backload the guide wire into the distal end of the catheter. Advance the guidewire into the catheter until the guidewire exits from the wire exit port. Place the guide catheter and Y-adapter. Introduce the guidewire and advance it to the region of interest. Introduce imaging catheter into guide catheter.

Note: Guidewires that supply more stiffness near the distal tips are recommended.

Note: Always wipe down the guidewire with heparinized saline prior to loading the catheter onto the guidewire.

CAUTION: Never advance the imaging catheter without guidewire support.

CAUTION: Never advance or withdraw the imaging catheter without the imaging core assembly in the most distal position.

CAUTION: Never advance or withdraw the imaging catheter without direct, fluoroscopic visualization.

CAUTION: Never advance the distal tip of the imaging catheter near the very floppy end of the guidewire. This part of the guidewire will not adequately support the catheter. A catheter advanced to this portion may not follow the guidewire when it is retracted and cause the guide wire to buckle into a loop. The catheter may then drag along the inside of vessel and catch on the guide catheter tip. If this occurs, remove the catheter assembly, guidewire and the guide catheter together. If the catheter is advanced too near the end of the guidewire, advance the guidewire while holding the imaging catheter steady. If this fails, withdraw the catheter and guidewire together.

Continue to advance the imaging catheter into the guide catheter, up to the femoral marker. Tighten the hemostasis valve on the guide catheter's Y-adapter. Tighten only enough to prevent fluid/blood leakage.

Note: AN EXCESSIVELY TIGHTENED HEMOSTASIS VALVE MAY DISTORT THE IMAGE DUE TO BINDING OF THE ROTATING DRIVE CABLE.

### Catheter Placement and Imaging

With the PIM image "OFF" and using fluoroscopy, advance the imaging catheter over the guide wire until the distal marker crosses a minimum of 3 cm beyond the region of interest in the vessel/lesion.

Keeping the catheter body and guidewire fixed, turn PIM image "ON" and retract the imaging core slowly along its 150 mm travel, imaging any region of interest.

Note: Always turn the PIM image "OFF" before advancing the imaging core within the catheter.

When finished, stop imaging by pressing the IMAGE button on the PIM, in manual mode, advance the imaging core to its most distal position. Maintain the position of the wire and remove the catheter.

### Troubleshooting

If your system menu does not include "Revolution Catheters", contact your Volcano Corporation representative before proceeding. If the images fade during use, flush the catheter with heparinized saline. If shadowed areas persist after flushing in situ, the distal lumen or catheter body may contain air bubbles.

### STORAGE AND HANDLING:

Products should be stored in a dry place with the temperature not exceeding 54 degrees Celsius (54° C) in their original cardboard box.

### PRODUCT SPECIFICATIONS:

Model	Revolution
Catalog number	89000
Crossing profile at transducer	3.2F
Maximum guide wire	0.014" (0.36 mm)
Minimum guide catheter	6F (1.63 mm)
Usable length	135 cm

Acoustic Output Parameter	B-Mode
I <sub>SPTA.3</sub> (mW/cm <sup>2</sup> )*	70.778
I <sub>SPPA.3</sub> (W/cm <sup>2</sup> )*	95.533
Pr.3 (MPa)**	1.901
PD (μs)	0.048
PRF (Hz)	15360
Center Freq (MHz)	42.3
MI**	0.281
TI*	0.010

Uncertainty

\* +/- 29.1%

\*\* +/- 14.6%

TI: Thermal Index defined as  $TI = \frac{W_{01x1}f_c}{210}$

W<sub>01x1</sub>: Bounded-square Output (mW)

f<sub>c</sub>: Center Frequency (MHz)

MI: Mechanical Index defined as  $MI = Pr.3/(f_c^{1/2})$

I<sub>SPPA.3</sub>: Derated Intensity, Spatial Peak Pulse Average (W/cm<sup>2</sup>)

I<sub>SPTA.3</sub>: Derated Intensity, Spatial Peak Temporal Average (mW/cm<sup>2</sup>)

Pr.3: Derated Peak Negative Pressure at a location of the maximum derated pulse intensity integral (MPa)

W<sub>0</sub>: Total Power (mW)

PD: Pulse Duration (μs)

PRF: Pulse Repetition Frequency (Hz)

### **LIMITED WARRANTY:**

Subject to the conditions and limitations on liability stated herein, Volcano Corporation ("VOLCANO") warrants that the Revolution catheter (the "Catheter"), as so delivered, shall materially conform to VOLCANO'S then current specification for the Catheter upon receipt for a period of one year from the date of delivery. ANY LIABILITY OF VOLCANO WITH RESPECT TO THE CATHETER OR THE PERFORMANCE THEREOF UNDER ANY WARRANTY, NEGLIGENCE, STRICT LIABILITY OR OTHER THEORY WILL BE LIMITED EXCLUSIVELY TO CATHETER REPLACEMENT OR, IF REPLACEMENT IS INADEQUATE AS A REMEDY OR, IN VOLCANO'S OPINION, IMPRACTICAL, TO REFUND OF THE FEE PAID FOR THE CATHETER. EXCEPT FOR THE FOREGOING, THE CATHETER IS PROVIDED "AS IS" WITHOUT WARRANTY OF ANY KIND, EXPRESSED OR IMPLIED, INCLUDING WITHOUT LIMITATION, ANY WARRANTY OF FITNESS, MERCHANTABILITY, AND FITNESS FOR A PARTICULAR PURPOSE OF NONINFRINGEMENT. FURTHER, VOLCANO DOES NOT WARRANT, GUARANTEE, OR MAKE ANY REPRESENTATIONS REGARDING THE USE, OR THE RESULTS OF THE USE, OF THE CATHETER OR WRITTEN MATERIALS IN TERMS OF CORRECTNESS, ACCURACY, RELIABILITY, OR OTHERWISE. Licensee understands that VOLCANO is not responsible for and will have no liability for any items or any services provided by any persons other than VOLCANO. VOLCANO shall have no liability for delays or failures beyond its reasonable control. Additionally, this warranty does not apply if:

1. The Catheter is used in a manner other than described by VOLCANO in the Instructions For Use supplied with the Catheter.
2. The Catheter is used in a manner that is not in conformance with purchase specifications or specifications contained in the Instructions For Use.
3. The Catheter is re-used or re-sterilized.
4. The Catheter is repaired, altered, or modified by other than VOLCANO authorized personnel or without VOLCANO authorization.

If claims under this warranty become necessary, contact VOLCANO for instructions and issuance of a Return Material Authorization number if the Catheter is to be returned. Equipment will not be accepted for warranty purposes unless the return has been authorized by VOLCANO.

This product, and the use thereof, may be covered by one or more of the following U.S. and international patents: 5081993; 5183048; 5257629; 5453575; 5601082; 5921931; 6036650; 6165128; 6200268; 6213950; 6381350; DE69031685.2; EP0386058; EP0473665; EP0611291; EP0707453; EP0871043; GB2233094; GB2246632; JP2905489; JP3188470; JP3194582. Other U.S. and international patents pending.

This product is licensed to the customer for single use only. Any re-sterilization or subsequent re-use is an unlicensed use and therefore constitutes patent infringement.

Revolution is a trademark of Volcano Corporation and is registered in the United States and other countries.

Volcano and the Volcano logo are trademarks of Volcano Corporation and are registered in the United States and other countries.


ADDITIONAL QUESTIONS REGARDING THIS PRODUCT SHOULD BE DIRECTED TO:

#### **Manufactured By:**

Volcano Corporation  
2870 Kilgore Road  
Rancho Cordova, CA 95670 USA  
(800) 228-4728 (916) 638-8008 (916) 638-8112 fax


#### **Authorized European Representative:**

Volcano Europe SA/NV  
Excelsiorlaan 41  
B-1930 Zaventem Belgium  
+32.2.679.1076 +32.2.679.1079 fax

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Revolution es una marca comercial de Volcano Corporation, y está registrada en los Estados Unidos y otros países.