

# DCMS

## Display Consistency Management Software

### User Manual

**K5960008-03**

---

The logo for BARCO, consisting of the word "BARCO" in a bold, sans-serif font with a small circle to the right of the "O".

Visibly yours

---

*Intentionally left blank*

---

# 1 About

## 1.1 Contents of the user manual

This guide consists of the following chapters:

**Chapter 1: About**

**Chapter 2: Important notice**

**Chapter 3: About DCMS**

**Chapter 4: Installing DCMS**

**Chapter 5: Connecting the components**

**Chapter 6: Running the DCMS application**

**Chapter 7: Display connection and calibration options**

**Chapter 8: Calibration process**

**Chapter 9: Verification**

**Chapter 10: The log file**

**Chapter 11: List of abbreviations**

**Chapter 12: Appendix A: Open Source Software**

**Chapter 13: Appendix B: Specifications**

## 1.2 Change record

Revision	Date	Description
00	March-08	Initial Release
01	May-08	Added/changed RS232 connection issues
02	July-09	Added NSL-4601
03	April-10	Software update V1.03

## 1.3 Notation convention

Following notations are applicable to this manual and should be respected through the manual.

---



**WARNING:**

Warnings – presented in this manual, provide information, which if not adhered to, may result in personal injury or death.

---



**CAUTION:**

Cautions – presented in this manual, provide information, which if not adhered to, may result in damage to the equipment.

---



**NOTE:**

Notes – presented in this manual, provide information, which emphasize points, significant to understand and operate the unit.

---



**IMPORTANT:**

Important – presented in this manual, provide information, which is important to highlight.

---

## 2 Important notice

### 2.1 Commercial In Confidence

The information contained herein is Barco confidential information. No part of the information contained herein may be disclosed outside of the organization of the recipient, its sub-contractors, and customers in any form or by any means and/or stored in a database or retrieval system without the prior written consent of Barco.

### 2.2 Copyright

This document is copyrighted. All rights are reserved. Nor this document, nor any part of it, may be reproduced or copied in any form or by any means - graphical, electronic, or mechanical including photocopying, taping or information storage and retrieval systems - without written permission of Barco.

© Barco N.V. All rights reserved.

### 2.3 Open Source Software

The RXTXComm library used to communicate over RS232 is Open Source. The library carries the RXTX License (<http://users.frii.com/jarvi/rxtx/license.html>), which is a modified LGPL 2.1. LGPL can be found at <http://www.gnu.org/licenses/lgpl.txt>. The RXTX library is included in its whole (including sources) to the DCMS CD-ROM.

The RXTX extension to the LGPL is for completeness listed in appendix A.

### 2.4 Embedded pictures

**NOTE:**

The pictures embedded in this document are for illustration purposes only and may look different depending on the device used to visualize this document (printer, display, ...).

---

### 2.5 Product safety notice

**CAUTION:**

Check the documentation of your display before physically connecting the monitor. Check that the cable you are about to use is fit for use with your monitor.

---

## 2.6 Software License Terms

**1** This Software License is between you and BARCO NV, a corporation organized and existing under the laws of Belgium registered under number BE 0473.191.041, Commercial Companies' Register of Kortrijk, having its registered office President Kennedypark, 35 at B-8500 Kortrijk, Belgium ("Barco") for the use of the Software.

**2** The Display Consistency Management Software and related documentation which are supplied herewith (together "Software") and the relating USB Sensor ("Materials"), used to operate the equipment identified in the documentation ("Equipment"), is licensed to you under this software license agreement ("Software License"), for use by any persons employed by or otherwise authorized by you ("Authorized Users"), subject to the use limit stated herein and accepted by you when agreeing to enter into this License.

You must read this Software License and agree to their terms before you may install and run the Software that has been supplied to you. You indicate YOUR ACCEPTANCE of this Software License by breaking the seal on the CD wallet and/or arranging for installation of the Software, and/or by proceeding to run the Software. If you do not accept the terms set out in this Software License, please return the software to Barco NV ("Barco") or the appointed Barco reseller who sold the Software to you and, as appropriate, discontinue the use of the Software.

**3** The Software is intended to be used to calibrate the displays identified in the documentation and you agree to limit the use of the Software for such intended purpose. You hereby undertake to inform all Authorized Users of the terms of this Software License and to bind all Authorized Users to accept all such terms of this Software License as applies to them. In consideration of payment of the agreed upon license fees, Barco grants you a limited, non-exclusive, non-assignable, non-transferable license (without the right to grant sublicenses) and right to use the Software by you and/or any Authorized User(s) to operate the Equipment.

**4** Barco (and Barco's licensors, as appropriate) retain ownership of all intellectual property rights in the Software and any copies you or any Authorized User may make of such Software. The Software is protected by national copyright laws, international copyright treaties and conventions, and other applicable laws. All rights not expressly licensed to you in this Software License are reserved to Barco and Barco's licensors, as appropriate. The Software contains certain other licensed materials (and including open source software) and Barco's licensors may protect their rights in the event of any violation of this Software License. Neither you nor any Authorized User may, whether in whole or in part, copy, translate, reverse engineer, derive source code from, modify, disassemble, decompile, create derivative works based on the Software, or remove any proprietary notices or labels on the Software, save as may be permitted by law, under the open source software license terms or this Software License, without the prior consent, in writing, of Barco.

**5** You and Authorized Users are entitled to use the Software for the purposes and in the manner set out herein, but neither you nor any Authorized User are entitled to: (i) sell or grant a security interest in the Software to other parties in any way, or to rent, lease or sub-license the Software to others without the express prior written consent of Barco; or (ii) exploit the Software or any of its component parts for any commercial purpose, other than use by you and/or Authorized Users of the Software as an ancillary visualization platform for computer integrated manufacturing data.

**6** The duration of this Software License will be from the date of your acceptance (as set forth above) of the Software, with no termination date, unless otherwise specified. You may terminate this Software License at any time by destroying all

copies of the Software then in your possession and returning all associated materials and documentation, to Barco or the appointed Barco reseller that sold or provided these to you. Barco may terminate this Software License forthwith by informing you at any time if you and/or any Authorized User are in breach of any of the Software License's terms.

**7** Barco warrants that (a) the Software will perform substantially in accordance with the accompanying written materials for a period of three (3) months from the date of your receipt of the Software (the "Warranty Period"), and (b) the Materials which may accompany the Software will be free from defects in materials and workmanship under normal use and service for a period of twelve (12) months from the date of your receipt of the Software and Materials.

You acknowledge that the Software is a complex computer software application, and that the performance thereof may vary depending hardware platform, software interactions and configuration. You acknowledge that the Software is not designed and produced specifically to meet your specific requirements and expectations and the selection of the Software by you is entirely your own choice and decision.

**8** Except for the open source software incorporated in or delivered with the Software, Barco shall hold you harmless and indemnify you from and against direct damages, losses and expenses arising from infringement or alleged infringement of any patent, trademark or copyright of such third party by the license and the right to install the Software as permitted by this Software License and defend and settle at its sole expense any claim, action, suit or proceeding brought against you, provided that (i) you promptly notify Barco in writing after a claim has been asserted against you or the commencement of any claim, action, suit or proceeding, and (ii) Barco shall assume sole control of the defense and any settlement negotiations related to any claim, action, suit or proceeding, and (iii) you shall not negotiate, settle or compromise any claim, action, suit or proceeding without the prior written consent of Barco and (iv) you, at your cost, shall cooperate with Barco and provide assistance and support, as may reasonably required by Barco, in connection with the defense and any settlement negotiations related to any claim, action, suit or proceeding. Barco shall have no indemnity obligation for any Software, or any portion thereof, (i) that is based on specifications, drawings, models or other data furnished by you or, (ii) that is not provided by Barco or, (iii) that is modified, in spite of the prohibition for you to modify the software or, (iv) to the extent that you continue allegedly infringing activity after having been provided modifications that avoid the alleged infringement, or (v) where the use of the Software, or the combination or thereof with other products, processes or materials or the distribution thereof rather than the Software itself is the primary cause of an alleged infringement. In case it has been determined by a finally awarded judgment that Barco infringed or misappropriated such third party rights or earlier, at Barco's discretion, it may, at its option and cost, (i) modify the Software in such a way that it shall not infringe upon or misappropriate the rights of the third party or (ii) obtain for you a license or other right to use the rights allegedly infringed or (iii) replace the Software in question with non-infringing Software. The remedies set forth in this paragraph shall constitute your sole and exclusive remedy and Barco's sole and exclusive liability for a third party claim that the Product infringes or misappropriates any intellectual property right of a third party.

**9** YOU UNDERSTAND THAT THE SOFTWARE IS BEING PROVIDED TO YOU "AS IS" WITH ONLY A LIMITED WARRANTY AS SET OUT IN CLAUSES 7 AND 8 HEREOF. BARCO DOES NOT MAKE NOR INTENDS TO MAKE ANY OTHER WARRANTIES OR REPRESENTATIONS, EXPRESS OR IMPLIED AND SPECIFICALLY DISCLAIMS ALL IMPLIED WARRANTIES OF MERCHANTABILITY FITNESS, FOR A PARTICULAR PURPOSE AND NON-INFRINGEMENT OF INTELLECTUAL PROPERTY AND DOES NOT WARRANT THAT THE SOFTWARE WILL BE FREE FROM ERRORS OR THAT

SUCH ERRORS WILL BE CORRECTED BY BARCO AND YOU ARE SOLELY RESPONSIBLE FOR ALL COSTS AND EXPENSES ASSOCIATED WITH RECTIFICATION, REPAIR OR DAMAGE CAUSED BY SUCH ERRORS.

**YOU ALSO ACKNOWLEDGE AND AGREE THAT:**

(i) BARCO ACCEPTS NO LIABILITY FOR ANY DAMAGES, LOSSES OR CLAIMS YOU OR ANY THIRD PARTY MAY SUFFER AS A RESULT OF YOUR USE OF THE SOFTWARE, AND YOU HEREBY AGREE TO INDEMNIFY, KEEP INDEMNIFIED, DEFEND AND HOLD HARMLESS BARCO AND BARCO'S AFFILIATES AND SUBSIDIARIES FROM AND AGAINST ANY AND ALL ACTIONS, PROCEEDINGS, LIABILITY, LOSS, DAMAGES, FEES AND COSTS (INCLUDING ATTORNEYS' FEES), AND OTHER EXPENSES INCURRED OR SUFFERED BY BARCO ARISING OUT OF OR IN CONNECTION WITH ANY BREACH BY YOU OF THE TERMS OF THIS SOFTWARE LICENSE. TO THE MAXIMUM EXTENT PERMITTED BY LAW, IN NO EVENT WILL BARCO BE LIABLE FOR ANY INDIRECT, SPECIAL, PUNITIVE, INCIDENTAL OR CONSEQUENTIAL LOSS OR DAMAGES OF ANY KIND WHICH MAY ARISE OUT OF OR IN CONNECTION WITH THE SOFTWARE, THIS SOFTWARE LICENSE OR THE PERFORMANCE OR PURPORTED PERFORMANCE OF OR FAILURE IN THE PERFORMANCE OF BARCO'S OBLIGATIONS UNDER THIS SOFTWARE LICENSE OR FOR ANY ECONOMIC LOSS, LOSS OF BUSINESS, CONTRACTS, DATA, GOODWILL, PROFITS, TURNOVER, REVENUE, REPUTATION OR ANY LOSS ARISING FROM WORK STOPPAGE, COMPUTER FAILURE OR MALFUNCTION OF THE SOFTWARE AND ANY AND ALL OTHER COMMERCIAL DAMAGES OR LOSSES WHICH MAY ARISE IN RESPECT OF USE OF THE SOFTWARE, EVEN IF BARCO HAS BEEN ADVISED OF THE POSSIBILITY OF THEIR OCCURRENCE; AND

(ii) YOU ACKNOWLEDGE THAT NO WARRANTY IS GIVEN BY BARCO ON THE OPEN SOURCE SOFTWARE INCLUDED IN THE SOFTWARE; AND

(iii) BARCO'S SOLE LIABILITY AND YOUR SOLE REMEDY UNDER THIS SOFTWARE LICENSE SHALL BE THE REPLACEMENT OF THE SOFTWARE IF ANY CD-ROM CARRYING THE SOFTWARE OR THE MATERIAL IS PROVEN TO BE DEFECTIVE BY BARCO'S PERSONAL INSPECTION, IN EACH CASE WITHIN THE PERIOD SET OUT IN CLAUSE 7, FOLLOWING INSTALLATION.

**10** Each party shall treat as confidential all information obtained from the other pursuant to this Software License which is marked 'confidential' or the equivalent or has the necessary quality of confidence about it and shall not divulge such information to any persons (except to such party's own employees or agents and then only to those employees or agents who need to know the same and are obligated in writing to confidentiality at least to the same extent as provided herein) without the other party's prior written consent provided that this Paragraph 6 shall not extend to information which was rightfully in the possession of such party prior to the commencement of the negotiations leading to this Software License, which is already public knowledge or becomes so at a future date (otherwise than as a result of a breach of this paragraph 10), is required to be disclosed by law or which is trivial or obvious. Each party shall ensure that its employees are aware of and comply with the provisions of this paragraph 10. The foregoing obligations as to confidentiality shall survive any termination of this Software License.

**11** This Software License (and documentation which is hereby integrated herewith) is the only understanding and agreement between you and Barco for use of the Software by you and/or Authorized Users. The Software License supersedes all other communications, understandings or agreements we had prior to this Software License (with the exception of any continuing confidentiality agreement) although nothing in this Software License purports to exclude liability for fraudulent misrepresentation. You may not export or re-export the Software or any copy or adaptation in violation of any applicable laws or regulations. Neither this Software License shall be altered, amended or varied save by an instrument in writing executed by duly authorized representatives of



both parties. Except as otherwise expressly provided for in this Software License, no forbearance, delay or indulgence by either party in enforcing the provisions of this Software License shall prejudice or restrict the rights of that party, nor shall any waiver of such parties rights operate as a waiver of any subsequent breach, and no right, power or remedy herein conferred upon or reserved for either party, is exclusive of any other right, power or remedy available to that party and each such right, power or remedy shall be cumulative. If any provision of this Software License is determined by the parties to be illegal, void or unenforceable, or if any court of competent jurisdiction in any final decision so determines, this Software License shall continue in full force save that such provision shall be deemed to be deleted with effect from the date of such decision, or such earlier date as the parties may agree.

**12** The construction, validity and performance of this Software License shall be governed in all respects by the laws of Belgium without recourse to its conflict of law principles. All disputes arising in any way out of or affecting this Software License shall be subject to the exclusive jurisdiction of the courts of Kortrijk, to which the parties agree to submit, without prejudice to enforcement of any judgment or order thereof in any other jurisdiction.

**13** The United Nations Convention on Contracts for the International Sale of Goods (the "Convention") shall not apply to this Software License, however, if the Convention is deemed by a court of competent jurisdiction to apply to this Software License, Barco shall not be liable for any claimed non-conformance of the Software under Article 35(2) of the Convention.

## 3 About DCMS

### 3.1 What is DCMS?

The Barco Display Consistency Management Software (DCMS) is a software tool to calibrate displays in terms of luminance, gamma and gray tracking\*. DCMS is able to calibrate displays to a predefined set of parameters or match displays to each other. The software uses a USB color sensor to measure color and luminance.

The term Display Consistency applies to the consistency of color throughout the gray scale and the consistency of color and luminance between multiple displays of the same or different type.

*\* Gray tracking means that the color of every gray level from black to white is constant, and equal to the white point target. In practice, because of LCD limitations, this holds for gray levels higher than 20% drive.*

### 3.2 Why using DCMS?

Because of spread in backlight, color filter and liquid crystal properties, all LCDs have different native white points and gamma. The same user settings will effectively produce different images, visible when the displays are put next to each other. Appropriate calibration can solve this by matching the luminances, gammas and white points.

This calibration requirement is general for all markets that use tiled LCD displays. More specifically, broadcasters must be able to trust the color they are seeing, so as good as possible absolute calibration (to an absolute, defined target) is also necessary. This is a matter of cost – used sensor, bit-depth of display and electronics, cost of color correction software and cross-talk removal, spread of color filters etc.

DCMS consists of software and sensor performing the given calibration.

## 3.3 Key components

### 3.3.1 Calibration system

#### **Calibration with internal test pattern generation (without disconnecting the display from its source)**

The calibration can be done without disconnecting the display. In this case, a laptop or desktop PC is needed to run the DCMS software. This laptop/PC needs to be connected to the RS232 service port of the display. The USB color sensor should be connected to the laptop/PC. The internal test pattern generator of the display is used to create the colored fields needed for the calibration.



Figure 1: External calibration system

#### **Calibration with external computer graphics (by using the input source as calibration system)**

Another way of setting up the calibration system is by using the input source as calibration system. The computer is connected to the display via RS232 and the color sensor is connected via USB. In this scenario, the colored test fields are generated as computer graphics.

For this, the display needs to be connected through DVI and driven at native resolution at 24 bits per pixel color depth.




Figure 2: The input source as part of the calibration system

There is no difference in calibration results between the two calibration systems.



**IMPORTANT:**

When calibration is done with external computer graphics, please restore all color settings of the host's graphics controller to their default value and disable any host services that may affect the color of the output picture.

---

### 3.3.2 Display

Each pixel of the display has 3 sub-pixels. These sub-pixels each have a color that corresponds with the three wavelengths for which the human eye is the most sensitive. These colors are red, green and blue (RGB). Each of these sub-pixels typically has 256 discrete levels. Colors are reproduced by combinations of these three primary colors.

Gamma is the perceptual linearity of light. Gamma is expressed by a power function.

Color temperature (or white point) is the 'color' of white. We say that the display has a warm white point or color temperature if the white of the display is reddish. The display has a cold or cool white point or color temperature if the white of the display is bluish. The color temperature or white point is expressed as coordinates in a color space (for example: xy values in CIE1931xyY space) or the temperature in degrees Kelvin at which a black body gets this color. Contradictory, we speak of warm white when the color temperature is low (for example: 3200K) and cold white when the color temperature is high (for example: 9300K). Daylight has a color temperature of 6500K.

Gray tracking is keeping the color temperature constant for all gray levels between black and white.

The brightness (light intensity of full white) is also calibrated as the power applied to the backlights. If the brightness target is set lower than the dimming capabilities of the backlight driver then additional dimming can be applied in the video if the display supports this function.

### 3.3.3 Sensor

The sensor used for calibration is the Barco USB color sensor (not sold separately). This sensor has a USB interface to the host PC running the DCMS application. Accuracy of the sensor is  $< 0.005$  xy in CIE1931xyY color space. Luminance accuracy is  $< 5\%$ . The sensor is shipped with a suction cap to attach the sensor to the monitor and a pod for dark current calibration of the sensor.

## 3.4 Package content

- USB color sensor (Marked 'Configured for DCMS')
- RS232 cable male-female, 3-wire straight, 1.8m
- RS232 cable male-female, 3-wire crossed, 1.8m
- Calibration sheet of the sensor
- Sensor pod for sensor calibration
- Contra weight for sensor
- Felt cup to prevent damage to displays without front glass
- Suction cup for attaching sensor to front glass of display
- CD-ROM containing DCMS software
- CD-ROM containing DCMS User Manual
- Installation manual
- DCMS license certificate

## 4 Installing DCMS

### 4.1 System requirements

- Host: Pentium II or later
- Memory: 256MB
- Hard disk space: 80MB (including Java Run-Time Environment)
- Drives: CD-ROM + HDD
- Operating system: Windows XP
- Run-Time Environments: Java SE 1.5 minimum (installed with DCMS)
- Boards & drivers: Driver for USB Optisense color sensor is available in the installation folder.
- Display: DVI capable of driving the display at its native resolution in 24bpp (optional)
- Other peripherals: USB, RS232

### 4.2 Before installing DCMS

- You must have Administrator privileges to install the software
- If there are older versions of DCMS installed, remove them first. Refer to the paragraph "*Uninstalling DCMS*" further in this chapter.


### 4.3 Installing DCMS

The DCMS software installer is available as an executable '**install.exe**' in the following folder on the installation CDROM: Disk1\InstData\VM.

*The DCMS installer is a standard application installer which guides you through the installation and gives the option to go to the next window, go back to the previous window or cancel the installation at any time.*


To install the DCMS software, follow these steps:

- **Double click the 'install.exe'**



The installer unpacks the necessary files. A dialog with progress bar is displayed.


- **Introduction**



Read the introduction and press '**Next**' to continue.

It is strongly recommended to quit all programs before continuing with the installation. If you want to change something on a previous screen, click the 'Previous' button. The installation can be cancelled by clicking the 'Cancel' button.


- **License Agreement**



Read and accept the Software License Terms.


As long as the licence agreement is not accepted, the installation process cannot be continued. Press '**Next**' to continue.

- **Enter License Key**



To continue the installation, enter your personal License Key, supplied on the DCMS license certificate. Press **'Next'** to continue.

- **Choose Install Folder**




Type or choose (with the folder browser) a location to install DCMS. The default location is 'C:\Program Files\Barco\DCMS'. By pressing the 'Restore Default Folder' button, you can always reset the default location.

Press **'Next'** to continue.



- **Choose Shortcut Folder**




Choose the location where the shortcuts will be located. There is a shortcut available for the DCMS application in the Barco DCMS program Group.

It is also possible to skip the creation of shortcuts.

Press '**Next**' to continue.

- **Pre-Installation Summary**



This window gives a summary of the selected/specified options. This panel gives the user the chance to review the settings before starting the actual installation.


To start the actual installation, press '**Install**'.

- **Installing...**



During the installation, a dialog with progress bar is displayed.

- **Install Complete**




The Install Complete panel shows if the installation was successful or not. If not successful, the installer will refer to the installation log for more details.


## 4.4 Uninstalling DCMS

To remove DCMS from your system:

- Be sure you have Administrator privileges to uninstall the software
- Exit the DCMS program
- From the Start menu, select Settings > Control Panel
- Double-click the Add/Remove Programs icon
- Select DCMS from the list box
- Click the Add/Remove... button



- Click the Uninstall button to remove the DCMS program.



**NOTE:**

Log files and files containing calibration parameters are not automatically removed.

---

## 5 Connecting the components

### 5.1 Optical sensor connection

**CAUTION:**

The USB sensor is delivered with a suction cup head that can be attached to the sensor for easy attachment to a display with front glass. The use of this suction cup on displays without front glass is not advised. Therefore a felt cup head is also delivered.

If you decide to use the suction cup head despite this warning, be sure to turn the sensor when removing the sensor from the LCD and do not pull it off perpendicular to the LCD.

To change the sensor head:

- Remove the suction cup head by turning the head adapter in counterclockwise sense as illustrated below.
- Place the felt cup sensor head on the sensor and turn the head adapter in clockwise sense.




Figure 3: Sensor head removal

To connect the sensor:


- Plug in the sensor in one of the USB connectors of your PC. See also 'Calibration system connections' paragraph below.
- Windows pops up a message box 'Found new hardware'



- Select 'No, not this time' and click 'Next' to continue.




- Select 'Install the software automatically' and click 'Next' to continue.



- Browse to the path where the file 'seqcal.sys' is located and click 'OK'. If the DCMS is successfully installed, the file should be available in the sub-folder 'drivers' of the installation folder on your hard disk.



## 5.2 Calibration system connection



- 1 | USB connection color sensor
- 2 | RS232 connection from display to PC
- 3 | DVI connection from display to PC

Figure 4: Calibration system connection

- Plug the sensor (delivered together with the DCMS software) in one of the USB connections of your PC. See also 'Optical sensor connection' paragraph. The USB color sensor measures color and luminance.
- Connect the display to be calibrated to the PC by using the correct RS232 cable. The RS232 is used to control the display and to upload the calibration parameters. The RS232 can also be used to set RGB test patterns in stead of DVI or when no DVI connection to the monitor is available.
- Connect a DVI cable (not delivered) between PC and display to be calibrated. The DVI connection is used to display color patches.



**NOTE:**

If the display is already connected to a different DVI source, the serial connection can be used to control the display's internal test pattern generator for colored test fields. In this case, the input source is part of the calibration system and the colored test fields are generated as computer graphics.

---



**IMPORTANT:**

The DVI of the PC must be capable of driving the display at its native resolution with 24 bit color resolution.

---





**CAUTION:**

Barco displays are connected with a 3-wire, crossed or straight serial cable. DCMS provides 2 serial cables: one 3-wire crossed cable and one 3-wire straight cable. The 3-wire crossed cable (Z3498894) is to be used for ISIS and MDPW displays; The 3-wire straight cable (K3494481) is to be used for ADP, NSL-4601 and LC displays.


---

# 6 Running the DCMS application

## 6.1 User Interface

### 6.1.1 DCMS User Interface

The DCMS User Interface window has 3 main panels:



- 1 | Connection panel
- 2 | Sensor panel
- 3 | Calibration target panel

Figure 5: DCMS User Interface

## 6.1.2 Connection panel

The user can select the display type and display communication details in the Connection Panel. The displays that are supported are the Barco ISIS, MDPW476, ADP361, NSL-4601 and the LCx4x series displays. The specifics for each of these displays are described later in this document. Display communication is done via serial communication. Every serial port of the host is listed in the 'COM port' combo box. The Baud rate of the selected serial port can be chosen in the 'Baud rate' combo box. In the 'Monitor ID' spinner, the user can select the serial communication ID of the monitor. This is useful for displays that are chained after each other on the same serial port. How to change the Monitor ID of your display is described in the User Manual of the display.

Once the correct display type and communication details are set, press 'Connect' to start communication with the display. On successful communication with the display, the button will change text to 'Disconnect' and the Status field will color green and show 'CONNECTED'. Press the button that is now labeled 'Disconnect' to stop communication with the selected monitor. The Status field will color red again and show 'DISCONNECTED'. Changing display type or communication details is not possible as long as the connection is active.

## 6.1.3 Sensor panel

### Calibration of the sensor

The Sensor Panel contains the USB sensor related functions. Press the 'Calibrate' button to open a connection to the sensor and calibrate it. After pressing the 'Calibrate' button, a dialog window will appear that explains the actions to be taken: Connect the USB color sensor to the PC and Place the USB color sensor on a dark surface or in the sensor pod and press 'Calibrate Sensor'. Perform this step in dark room conditions for optimal results.



Figure 6: Sensor Calibration Dialog

Press the 'Calibrate Sensor' button to start calibration. This process takes about 10 seconds to complete.

After successful connection and calibration of the sensor, the progress bar of the Sensor Calibration dialog will color green and show 'Calibration passed'. The Status field of the Sensor panel will color green and show 'CALIBRATED'. If connection or calibration failed, both progress bar and Status field will color red and show 'CALIBRATION FAILED' and 'NOT CALIBRATED'. The calibration will fail if the sensor is not plugged in or if the sensor is placed on a bright surface.

Press the 'Close' button to return to the main application.



**IMPORTANT:**

For optimal performance of the sensor it is necessary to place the sensor on a black and darkened surface or place the sensor in the sensor pod. An incorrect calibration of the sensor leads to invalid results. If the sensor is taken from a cold storage room, allow the sensor to acclimatize to the ambient temperature before use.

---

A recalibration of the sensor is required in the following conditions:

- the sensor was disconnected
- the ambient temperature has changed
- the last calibration of the sensor is more than a few hours old

**Single measurement**

For various reasons, the user may wish to know the color coordinates and luminance of a test patch on a display. This function is provided under 'Single Measurement'. This button hides a dialog box. Place the sensor on the sample and press 'Measure'. When the measurement is done, the progress bar will stop running and the color coordinates and luminance are shown. Press 'Close' to return to the main application window. The measurement is done in CIE1931xyY color space.




Figure 7: Single measurement

## 6.1.4 Calibration target panel

The Calibration Target Panel contains 8 sections. Each of these sections are listed and described below. Depending on the type of display that is selected, some of these sections may be disabled. A disabled section means that this functionality is not supported by the display.

### Test patch

The user can select the preferred way of measuring. To acquire calibration parameters (the 'Get' function that is described below) or to calibrate a display, the DCMS software needs to measure colored test patches. These test patches can be generated as a colored field in the application itself or use the internal built-in-test pattern generator that is available in each display.

Select 'External computer graphics' to generate the test patches as a colored field in the application itself. Choose this option if your host is the video source for the display that you wish to calibrate. Note that the calibration dialog window will always be on top of all applications to prevent obfuscation of the test patches by other applications or dialog windows during the calibration.




Figure 8: External test patch

Select 'Internal display pattern generator' to let DCMS instruct the display to generate the test patches. Choose this option if your host is not the video source for the display that you wish to calibrate. The size of the test patch depends on the display's implementation of this function.




Figure 9: Internal test patch

Note that the input source will not be visible for the duration of the calibration when the display pattern generator is selected.

### White point

The White point of a display is the color of full white. A color display's white is the combination of all three primary colors red, green and blue. The proportion red, green and blue defines the color of white.




Figure 10: Additive colors

The White point can be expressed in Color Temperature (the temperature in degrees Kelvin at which a heated black-body radiator matches the color of the light source's color). A bluish white has a high color temperature (commonly referred to as cold white) while a reddish white has a low color temperature (commonly referred to as warm white).



Example of 5000K



Example of 7500K

Figure 11: Color temperature examples




Figure 12: Color temperature scale\*

*\*This picture is licensed under Creative Commons Attribution ShareAlike 2.5 Poland License.*

Select the 'Color temperature' radio button to define the white point in color temperature and use the slider to set the desired color temperature between 3200 and 9900K. Note that daylight is typically 6500K. The selected color temperature can be read in the text field on right of the slider.

For expert use, you can also give the desired white point in CIE1931xyY coordinates. Figure below shows the CIE1931xyY color space. In this chart, all color temperatures between 2000 and 10000K are also shown.




Figure 13: CIE1931xyY color space

Select the 'CIE1931xyY' radio button if you wish to set the desired white point in CIE1931xyY color coordinates in the two text fields below, labeled 'x' and 'y'. For your convenience, the color temperature slider may still be used to set an initial x and y value to start from. In the text field on the right, labeled 'Tc(K)' you can read the correlated color temperature. Below this text field, a second text field labeled 'u'v'' shows how far the entered x and y coordinates deviate from the Planckian locus (this is the path that a black-body follows as it is heated). This deviation is expressed in distance in the CIE1976UCS color space (which is a color space that is proportional to human perception of color differences). Both these text fields will color red if no correlated color temperature can be calculated for the entered x and y coordinates.

### Mid gray

For some applications (like color matching of displays) the mid\* gray point can also be defined. This setting is intended for expert use. The desired color coordinates and luminance (in CIE1931xyY) of mid gray can be changed here.

*\*Mid gray in 3x8 bit RGB color space is defined here as (128, 128, 128)*

The three input text fields 'x', 'y' and 'Y' are disabled as long as the check box 'Follow white point and gamma' is selected. This means that changing white point, gamma or brightness settings will influence the mid gray point. The color coordinates will match those of the white point and the luminance is calculated from the gamma and the brightness.

Uncheck the 'Follow white point and gamma' will allow you to set the desired mid gray point if this should differ from the white point.



Figure 14: Match of white and mid gray color



Figure 15: Mismatch of white and gray color

### Gamma

Human vision has a non-linear perception of changes in luminance (light intensity). This non-linearity can be characterized by a power-law expression. This power-law is used in display technologies to map the number of possible gradations in luminance to a more perceptually uniform space. The Greek letter gamma is the power in this power-law expression. Typical gamma values for displays range from 1.8 to 2.8. The pictures below illustrate changes in gamma.




Figure 16: The gamma power-law






Figure 17: Gamma 4




Figure 18: Gamma 2



Figure 19: Gamma 1

The gamma slider can be used to set the desired gamma of the display. This value can also be seen in the text field on the right of the gamma slider. The luminance of the gray point changes by definition when the gamma changes.

Note that the gamma slider is disabled when the mid grey point is defined by the user (when the 'Follow white point and gamma' is unchecked). The user defined luminance of the mid grey point now defines by definition also the gamma. The text field of the gamma value will color red if the resulting gamma is out of the slider's range. So, the mid grey luminance is in function of the desired gamma or the gamma is in function of the desired mid grey luminance.

### Brightness

The brightness of a display is the maximum luminance of full white. The slider allows the user to set the desired brightness. Values between 25 and 250 cd/m<sup>2</sup> in steps of 5 can be chosen for almost all displays except for NSL-4601 displays, where the brightness can be set from 100 to maximum 650 cd/m<sup>2</sup> in steps of 10. The text field on the right shows the desired brightness. If the desired brightness cannot be reached because of the limitations of the display, the calibration routine will continue with the closest match but will report a failed calibration at the end of the routine.

### Status field

The Status field has 4 possible states:

	No display has been connected
Unknown Calibration History	Display
CALIBRATION SUCCESSFUL	Display was calibrated successfully
CALIBRATION FAILED	Display calibration failed, target settings were out of range or calibration was interrupted

### Get, Load, Save and Reset

Besides using the controls above, the user has the option to 'Load' a previously defined set of parameters from file or 'Save' the current set to file.

Use the 'Get' function to acquire these settings by measurement. This is useful for display matching purposes. How to match displays is explained in the section *"Using DCMS, Match displays"* further in this document.

Use the 'Reset' function to perform a reset of all calibration parameters (gamma, brightness,...) to their factory default values. Since not all parameters are supported by every display, unsupported parameters will be skipped.

### Calibrate

Once the parameters are set to the desired values, press the 'Calibrate now...' button to start the calibration routine. See the section on *"Using DCMS, Calibrate a display"* for a complete explanation of the functionality behind this button.




Figure 20: Barco DCMS – Calibration window

## 6.2 Using DCMS

### 6.2.1 Proper use and general remarks

- Store the sensor at the same ambient temperature as the display, facing down in the sensor pod.
- When placing the sensor on the test patch, lead the cable over the top of the display. Use the contra weight on the USB cable of the sensor to balance the weight of the sensor. This prevents damage to the sensor if the suction cup would release its grip.
- If the suction cup does not hold, slightly moisten the suction cup.
- Allow displays a warm up time of a few hours to stabilize before calibrating them.

### 6.2.2 Calibrate a display

In this section we will calibrate a display to a set of target values. Calibrating several displays is just a matter of repeating these steps for each of the displays that need calibration.




Figure 21: Display calibration

#### Step 1: Connection and startup

- Physically connect the display via serial connection to the host. Refer to the documentation of the display for pin-out of the serial port.
- Connect the USB sensor to the host.
- Start the DCMS application.

#### Step 2: Calibration and serial communication

- Select the display type, serial port number, baud rate and monitor ID (Default monitor ID is 1). Press 'Connect' and verify in the connection status field that the connection to the display says 'CONNECTED'.
- Start the sensor calibration routine by pressing 'Calibrate...'. Follow the instructions of the calibration dialog as mentioned in the *"User Interface,*

*Calibration target panel*" paragraph above. Verify that the sensor status field states 'CALIBRATED'.

### Step 3: Calibration target

- Select the type of test patches to be used. Select 'External computer graphics' if the host is the input source of the display and you want the test patterns to be created in the calibration dialog. Select 'Internal display pattern generator' when the display has an input source different from the host. See the *"User Interface, Sensor panel"* paragraph for more details on the selection of the test patch.
- Set the desired white point, mid gray, gamma and brightness by manually changing the settings or load a previous set from file. See the *"Sensor panel"* paragraph for details on calibration target settings and the paragraph *"Display connection and calibration"* for specifics of the display.

### Step 4: The calibration routine

- Press 'Calibrate now...' to open the calibration dialog.
- Follow the instructions as given in the calibration dialog and press 'Start' to start the calibration routine.
- During the calibration routine, the user may interrupt at any time by pressing the 'Close' button to return to the main application.



#### NOTE:

Note that the display will not return to a previously known state but will leave the display in the state that it was at the time of pressing the 'Close' button.

---

- The user is informed of each step of the calibration in the text area of the calibration dialog.
- During calibration the progress bar is colored blue.
- When the calibration is finished, the progress bar will display if the calibration was successful or not and color green when successful, red if failed. The same message will appear in the calibration status field of the main application window. As soon as this message appears, you may close this dialog box.
- If the dialog box is not closed after calibration, a verification process will start. This process measures all gray levels in steps of 16 from 0 to 255 and displays the measured results together with the desired values. See section 6 for details on the verification table.
- After the verification routine, the test patch will show a green color if the calibration was successful or red if it failed. Text in the text area of the calibration window can be copied and pasted into your favorite text editor and saved to file.
- Press 'Close' to return to the main application. Pressing 'Start' will redo the calibration routine all over again. If 'Internal display pattern generator' was selected, the display will switch back to the previously selected input when the calibration dialog is closed.

### 6.2.3 Match displays

In this section we will calibrate a display (target display) to match another display (reference display). The displays do not need to be of the same type.




Figure 22: Matching of displays

#### Step 1: Connection and startup

- Physically connect the reference display via serial connection to the host. Refer to the documentation of the display for pin-out of the serial port.
- Connect the USB sensor to the host.
- Start the DCMS application.

#### Step 2: Calibration and the reference display

- Select the display type, serial port number, baud rate and monitor ID (Default monitor ID is 1). Press 'Connect' and verify in the connection status field that the connection to the display says 'CONNECTED'.
- Start the sensor calibration routine by pressing 'Calibrate...'. Follow the instructions of the calibration dialog as mentioned in the 'Sensor Panel' section described above. Verify that the sensor status field says 'CALIBRATED'.

#### Step 3: Set parameters

- Select the type of test patches to be used. Select 'External computer graphics' if the host is the input source of the display and you want the test patterns to be created in the calibration dialog. Select 'Internal display pattern generator' if the display has an input source different from the host.

See the *"User Interface, Sensor panel"* paragraph above for more details on the selection of the test patch.

- There is no need to set other parameters as it is the objective of this process to acquire these from the reference display.

#### **Step 4: The acquisition**

- Press 'Get...' to open the acquisition dialog.
- Follow the instructions as given in the calibration dialog and press 'Start' to start the acquisition routine.
- During the acquisition routine, the user may interrupt at any time by pressing the 'Close' button to return to the main application.
- The user is informed of each step of the acquisition in the text area of the acquisition dialog.
- During acquisition the progress bar is colored blue.
- When the acquisition is finished, the progress bar will display if the acquisition was successful or not and color green when successful, red if failed. As soon as this message appears, you may close this dialog box.
- After the routine, the test patch will show a green color if the routine was successful or red if it failed. Text in the text area of the calibration window can be copied and pasted into your favorite text editor and saved to file. White point, mid gray, gamma\* and brightness are automatically copied to the main application window.
- Press 'Close' to return to the main application. Pressing 'Start' will redo the routine all over again. If 'Internal display pattern generator' was selected, the display will switch back to the previously selected input when the acquisition dialog is closed.

\* *Gamma and mid gray luminance may differ slightly if recalculated as gamma is calculated as a best fit from 8 gray levels between 0 and 255.*

#### **Step 5: The target display**

- Press 'Disconnect' to stop communication with the reference display.
- Physically connect the target display via serial connection to the host. Refer to the documentation of the display for the serial port pinning.
- Select the display type, serial port number, baud rate and monitor ID (Default monitor ID is 1) of the target display. Press 'Connect' and verify in the connection status field that the connection to the display says 'CONNECTED'.

#### **Step 6: Set parameters**

- Select the type of test patches to be used. Select 'External computer graphics' if the host is the input source of the display and you want the test patterns to be created in the calibration dialog. Select 'Internal display pattern generator' when the display has an input source different from the host. See the *"User Interface, Calibration Target Panel"* paragraph above for more details on the selection of the test patch.
- All calibration target values are automatically set during the acquisition routine on the reference monitor. If needed, changes can be still be made here.

**Step 7: The calibration routine**

- Press 'Calibrate now...!' to open the calibration dialog.
- Follow the instructions as given in the calibration dialog and press 'Start' to start the calibration routine.
- During the calibration routine, the user may interrupt at any time by pressing the 'Close' button to return to the main application.

**NOTE:**

Note that the display will not return to a previously known state but will leave the display in the state that it was at the time of pressing the 'Close' button.

- The user is informed of each step of the calibration in the text area of the calibration dialog.
- During calibration the progress bar is colored blue.
- When the calibration is finished, the progress bar will display if the calibration was successful or not and color green when successful, red if failed. The same message will appear in the calibration status field of the main application window. As soon as this message appears, you may close this dialog box.
- If the dialog box is not closed after calibration, a verification process will start. This process measures all gray levels in steps of 16 from 0 to 255 and displays the measured results together with the desired values. See section 6 for details on the verification table.
- After the verification routine, the test patch will show a green color if the calibration was successful or red if it failed. Text in the text area of the calibration window can be copied and pasted into your favorite text editor and saved to file.
- Press 'Close' to return to the main application. Pressing 'Start' will redo the calibration routine all over again. If 'Internal display pattern generator' was selected, the display will switch back to the previously selected input when the calibration dialog is closed.

## 6.2.4 Single measurement

The single measurement function allows the user to use DCMS as a simple luminance and chromaticity meter. Note that a serial connection to a display is not required for this function.

**Step 1: Connection and startup**

- Connect the USB sensor to the host.
- Start the DCMS application.

**Step 2: Calibration**

- Start the sensor calibration routine by pressing 'Calibrate...!'. Follow the instructions of the calibration dialog as mentioned in the 'Sensor Panel' section above. Verify that the sensor status field says 'CALIBRATED'.

### **Step 3: Measurement**

- Press the 'Single measurement' button to open the measurement dialog.
- Place the sensor on a color patch on the display to be measured and press 'Measure'. The luminance and xy color coordinates will appear in the text field as soon as the measurement is done. Press 'Close' to return to the main window or press 'Measure' for another measurement.



# 7 Display connection and calibration options

## 7.1 ADP361

### Connection

The Barco ADP361 can be connected to the host by means of a 3-wire straight (Rx-Rx, Tx-Tx, GND-GND) male-female cable. This cable is provided in the DCMS package.

### Calibration options




Figure 23: DCMS with ADP361

The ADP display supports all the functionality of DCMS. However, for calibration of the gray levels, only the mid gray (RGB 128, 128, 128) is calibrated. The color of the mid gray point can be set to a value different from the white point as this display is typically used in combination with the Barco MDP471 which has a mid gray point that is slightly different from its white point.

## 7.2 LCx-4x

### Connection

The Barco LCx-4x series displays (LCN/LCS-42/47) can be connected to the host by means of a 3-wire straight (Rx-Rx, Tx-Tx, GND-GND) male-female cable. This cable is provided in the DCMS package.

### Calibration options




Figure 24: DCMS with LCx-4x

The LCx-4x display supports all the functionality of DCMS. Note that the option to set the mid gray color is disabled. Mid gray calibration is not needed as a full gray tracking is done. Gray tracking is done with the help of high precision look-up tables embedded in the display and programmed by DCMS.

Note that a valid input source is required on this type of display for the internal test pattern generator to work as the generator uses the applied input timing for the generation of the test patches.

## 7.3 MDP471 (ISIS)

### Connection

The Barco MDP471 (code name ISIS) display can be connected to the host by means of a 3-wire crossed male-female cable (Rx-Tx, Tx-Rx, GND-GND). This cable is provided in the DCMS package.

### Calibration options





Figure 25: DCMS with MDP471

The ISIS display currently only supports the calibration of luminance. All other controls are disabled. In typical use, this display is often the reference display for the matching functionality of DCMS in combination with the Barco ADP361.

After calibration following dialog box pops up:



Here you can choose whether the calibration value should be stored or not. Pressing the Cal pushbutton on the MDP471 control unit will take the stored brightness value and set the brightness of your display to that value.

---



**NOTE:**

This can only be done when no OSD menu is active or when OSD level 1 is active. More info about the OSD and its levels can be found in the user manual of the unit.

---

## 7.4 MDPW476

### Connection

The Barco MDPW476 display can be connected to the host by means of a 3-wire crossed male-female cable (Rx-Tx, Tx-Rx, GND-GND). This cable is provided in the DCMS package.

### Calibration options




Figure 26: DCMS with MDPW476

The MDPW476 display supports all the functionality of DCMS. Mid gray calibration is not needed as a full gray tracking is done. Gray tracking is done with the help of high precision look-up tables (LUTs) embedded in the display and programmed by DCMS. The MDPW is able to store up to 4 color calibration settings that are also selectable via OSD menu (on-screen-display menu). Before the calibration dialog is displayed, an MDPW specifics dialog is displayed first. The progress bar below is colored blue while DCMS is retrieving this information from the display and will color green as soon as this information is available.




Figure 27: MDPW LUT selection

The user can now select which of the 3 look-up-table slots the user wishes to overwrite with the new calibration (slot 0 is reserved for a native 1-to-1 LUT). Select one of the 3 available slots. Each slot can also be given a name for easy identification of the calibration settings that is stored in the slots. Do not use special characters for the naming of these slots.

A second MDPW specific function is the 'Extended dimming'. In case of low ambient light, the user may wish to dim the maximum luminance beyond the maximum dimming point of the display's backlight system. This function uses the high-precision look-up-tables to dim the display beyond this point.

Press 'Cancel' to return to the main window application or 'Continue' to advance to the calibration dialog.

## 7.5 NSL-4601

### Connection

The Barco NSL-4601 displays can be connected to the host by means of a 3-wire straight (Rx-Rx, Tx-Tx, GND-GND) male-female cable. This cable is provided in the DCMS package.

For a single display calibration:

- Connect the RS cable (delivered with DCMS) and the signal cable (RGB/DVI) to the NSL-4601 display. Disconnect all other cables (e.g. DVI out).
- Make sure the display is set as master display.

For the calibration of a display in a wall structure:

- Connect the RS cable (delivered with DCMS) to the master display
- Make sure that you can display the BCM dialog on the to be calibrated display in the wall.  
See BCM manual for more detailed information about the BCM software package.

### Calibration options




Figure 28: DCMS with NSL-4601

The NSL-4601 displays support all the functionality of DCMS. Note that the option to set the mid gray color is disabled. Mid gray calibration is not needed as a full gray tracking is done. Gray tracking is done with the help of high precision look-up tables embedded in the display and programmed by DCMS.

The monitor ID is the address of a monitor. The monitor address is a double (x, y) where x represents the horizontal position (column) and y the vertical position (row) of the monitor in the wall. The top left monitor (seen from front) has the address (0,0): Numbering starts with zero and goes from left to right and from top to bottom. The address is set in binary code via two 4-dip switches at the rear side of the monitor, one for x and one for y. The following table shows how the dip-switch settings are mapped to monitor ID's:

		Horizontal Dip Switches															
		0	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15
Vertical Dip Switches	0	0	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15
	1	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31
	2	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47
	3	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63
	4	64	65	66	67	68	69	70	71	72	73	74	75	76	77	78	79
	5	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95
	6	96	97	98	99	100	101	102	103	104	105	106	107	108	109	110	111
	7	112	113	114	115	116	117	118	119	120	121	122	123	124	125	126	127
	8	128	129	130	131	132	133	134	135	136	137	138	139	140	141	142	143
	9	144	145	146	147	148	149	150	151	152	153	154	155	156	157	158	159
	10	160	161	162	163	164	165	166	167	168	169	170	171	172	173	174	175
	11	176	177	178	179	180	181	182	183	184	185	186	187	188	189	190	191
	12	192	193	194	195	196	197	198	199	200	201	202	203	204	205	206	207
	13	208	209	210	211	212	213	214	215	216	217	218	219	220	221	222	223
	14	224	225	226	227	228	229	230	231	232	233	234	235	236	237	238	239
	15	240	241	242	243	244	245	246	247	248	249	250	251	252	253	254	255

Figure 29: Monitor ID mapping table

Note that a valid input source is required on this type of display for the internal test pattern generator to work as the generator, as it uses the applied input timing for the generation of the test patches.

Once the calibration is completed, a dialog box appears to ask if you want to continue. When pressing continue, the calibration data will be written to the flash. When finished, the display will be automatically restart.




Figure 30: Calibration NSL-4601 completed



---

## Color versus brightness

---

**NOTE:**

Note that reducing the color temperature reduces the overall achievable brightness. The table below provides an indication of the relative brightness compared to the brightness at native color temperature, which is 10000K.


<b>Color temperature</b>	<b>Brightness (relative)</b>
3200K	between 50% and 55%
5000K	between 70% and 75%
5500K	between 75% and 80%
6500K	approx. 85%
7500K	between 90% and 95%
9300K	approx. 95%
10000K	100% as this is the native color temperature

Table 1: Color temperature versus brightness

---

## 8 Calibration

Calibration is done in 3 major steps:



- 1 | Calibration reset
- 2 | Measurement of native characteristics
- 3 | Iterative calibration

Figure 31: Calibration steps flowchart

Each of these steps is explained in detail in the following subparagraphs.

### Calibration reset

The first step in the calibration process is resetting all parameters that are to be calibrated or may influence the calibration. Not all parameters are supported by every display. Unsupported parameters are skipped.

Resetting these parameters to their default, neutral or disabled value allows us to measure the native characteristics of the display in the next step.

### Native characteristics

In the second step of the calibration process, the native characteristics are measured and calculated. First, the black is measured. This will later be used in the algorithms to compensate for black. In this context, black is not the total absence of light but is the darkest possible luminance that can be achieved. Since this black has a luminance that is not equal to zero, it also has a color and therefore needs to be reckoned with when calibrating the display.

Next, we loop over the three primary colors red, green and blue. For each of these colors, we measure luminance, color and native gamma.

Finally, we measure the luminance and color of the maximum native white. This is our starting point from where we calibrate to the required white point.

**Calibration**

The third step in the calibration process is the actual calibration. Everything is done by gray tracking calibration as full white can be seen as just a special kind of gray. This gray tracking is done for a reduced set of gray levels to reduce calibration time without compromising the end result.

## 9 Verification

The verification table below lists an example output of a verification table. After the calibration process, the user can start a verification of the results by pressing the 'Start verification' button in the calibration dialog. The layout of this table is always the same but the content may depend on the type of display that is calibrated. The example below is taken after calibration of a Barco MDPW476.

```

Verifying results ...
Target white point: (Y: 150.00, x: 0.313, y: 0.329)
  DAC  Y      Yth      dY(%)  dEuv*,L x      y      du'v'  dEuv*,C  dEuv*  Tc
  255  151.28  150.00    0.85   0.33  0.312  0.329  0.000    0.40   0.52  6522
  240  128.80  130.11   -1.02   0.37  0.312  0.330  0.001    1.06   1.12  6553
  224  110.63  110.66   -0.03   0.01  0.313  0.330  0.001    0.65   0.65  6499
  208  93.79   93.01    0.84   0.28  0.312  0.329  0.000    0.35   0.45  6515
  192  77.06   77.09   -0.05   0.01  0.312  0.329  0.000    0.48   0.48  6526
  176  62.49   62.87   -0.61   0.18  0.313  0.331  0.001    0.95   0.97  6468
  160  50.32   50.29    0.05   0.01  0.313  0.331  0.001    0.93   0.93  6483
  144  39.67   39.30    0.91   0.23  0.313  0.330  0.000    0.21   0.31  6486
  128  30.65   29.85    2.62   0.60  0.312  0.329  0.001    0.54   0.81  6562
  112  22.06   21.86    0.91   0.19  0.313  0.328  0.001    0.35   0.40  6497
   96  15.97   15.28    4.33   0.81  0.318  0.332  0.003    1.63   1.82  6223
   80  10.22   10.02    1.93   0.31  0.313  0.328  0.001    0.39   0.49  6521
   64   6.02    6.01    0.17   0.02  0.311  0.329  0.001    0.30   0.30  6586
   48   3.15    3.15   -0.13   0.01  0.305  0.320  0.006    1.34   1.34  7026
   32   1.35    1.33    1.57   0.13  0.300  0.309  0.013    1.37   1.38  7511
   16   0.54    0.42   22.34   0.75  0.299  0.294  0.023    0.95   1.19  7950
    0   0.19    0.19    0.17   0.00  0.313  0.280  0.036    0.53   0.53  7082

Gamma:
Target: 2.35
Measured (least square log-log): 2.34
Correlation coefficient      : 0.99905
VERIFICATION COMPLETED

```

Table 2: Verification table

Verification is done by measuring the color and luminance of a set of gray levels, including white. The measured results are compared to the calculated theoretical values and printed on screen.

The report shows:

Target white point:	Shows the target luminance and color coordinates for full white.
Target mid gray point:	Shows the target luminance and color coordinates for mid gray.
DAC:	Gray levels from 255 to 0 in steps of 16.
Y:	Measured luminance in cd/m <sup>2</sup> .
Y <sub>th</sub> :	Theoretical target luminance in cd/m <sup>2</sup> .
dY(%):	Luminance error in percentage from theoretical target.
dE <sub>uv*</sub> ,L:	Luminance error in $\Delta E_{uv^*}$ from theoretical target.
X:	Color coordinate x in CIE1931xyY color space.
Y:	Color coordinate y in CIE1931xyY color space.
du'v':	Color deviation from target in CIE1976UCS color space.
dE <sub>uv*</sub> ,C:	Chromaticity error (color error) in $\Delta E_{uv^*}$ .
dE <sub>uv*</sub> :	Total error of luminance and chromaticity in $\Delta E_{uv^*}$ .
T <sub>c</sub> :	Correlated color temperature in degrees Kelvin.
Gamma target:	Target gamma.
Gamma measured:	Gamma calculated as best fit over the measured gray levels. (Calculation is done with least square error linear regression on logarithmic levels and luminance)
Gamma correlation coefficient:	This coefficient gives an indication of the reliability of the calculated gamma. The result is more reliable as it approaches 1.

Table 3: Verification table explanation

## 10 The log file

Each time the DCMS application is started, a new log file is created. This file is named 'dcmsxxx.log' where xxx is the timestamp of creation. The log file contains all text that is displayed in the instructions panel of the 'calibration dialog' and 'get calibration dialog'. It also lists the measurement results of the verification routine and single measurements. The file is stored in the installation folder of DCMS.

# 11 List of abbreviations

**A**

ADP Auxiliary Display Panel  
ALC Ambient Light Controller

**B**

BCM Barco wall Control Manager  
BIT Built In Test  
Bkl Backlight  
BLOS Backlight Optical Stabilisation  
Bri Brightness

**C**

Cal Calibrate  
CCFL Cold Cathode Fluorescent Lamp  
cd candela  
COMM Communication  
Conn. Connector  
CRT Cathode Ray Tube

**D**

DCMS Display Consistency Management Software  
Diff Differential  
DVI Digital Video Interface

**L**

LCD Liquid Crystal Display  
LUT Look-Up Table

**M**

MDPW Main Display Panel Wide screen  
MISC Miscellaneous (keypad)  
ms Millisecond

**N**

N/A Not Applicable  
Nr Number  
NSL Near Seamless LCD

**O**

Opt Optical  
OSD On Screen Display

**P**

PWR Power

**R**

RGB Red Green Blue  
RX Receiver

**S**

SMPTE	Society of Motion Picture and Television Engineers
STD	Standard
SVGA	Super Video Graphic Array
SW	Software
SXGA	Super extended Graphics Array
S/N	Serial Number
<b>T</b>	
Temp	Temperature
TX	Transmitter
<b>U</b>	
UMAN	User Manual
USB	Universal Serial Bus
UXGA	Ultra eXtended Graphic Array
<b>V</b>	
V	Volt
VGA	Video Graphic Array
VHDL	Is the VHSIC Hardware Description Language
<b>W</b>	
WUXGA	Widescreen Ultra eXtended Graphics Array
<b>X</b>	
XGA	EXtended Graphic Array
<b>Y</b>	
Y	Luminance



# 12 Appendix A: Open Source Software

## DCMS copyright

DCMS is copyrighted by Barco nv with exception of the RXTXComm library. RXTXComm is open source software, licensed under the 'RXTX extension to LGPL'. The full RXTXComm and LGPL license can be found on the DCMS CD-ROM.

## RXTX extension to LGPL

RXTX License v 2.1 - LGPL v 2.1 + Linking Over Controlled Interface. RXTX is a native interface to serial ports in java. Copyright 1997-2007 by Trent Jarvi [tjarvi@qbang.org](mailto:tjarvi@qbang.org) and others who actually wrote it. See individual source files for more information.

A copy of the LGPL v 2.1 may be found at <http://www.gnu.org/licenses/lgpl.txt> on March 4th 2007. A copy is here for your convenience.

This library is free software; you can redistribute it and/or modify it under the terms of the GNU Lesser General Public License as published by the Free Software Foundation; either version 2.1 of the License, or (at your option) any later version.

This library is distributed in the hope that it will be useful, but WITHOUT ANY WARRANTY; without even the implied warranty of MERCHANTABILITY or FITNESS FOR A PARTICULAR PURPOSE. See the GNU Lesser General Public License for more details.

An executable that contains no derivative of any portion of RXTX, but is designed to work with RXTX by being dynamically linked with it, is considered a "work that uses the Library" subject to the terms and conditions of the GNU Lesser General Public License.

The following has been added to the RXTX License to remove any confusion about linking to RXTX. We want to allow in part what section 5, paragraph 2 of the LGPL does not permit in the special case of linking over a controlled interface. The intent is to add a Java Specification Request or standards body defined interface in the future as another exception but one is not currently available.

<http://www.fsf.org/licenses/gpl-faq.html#LinkingOverControlledInterface>

As a special exception, the copyright holders of RXTX give you permission to link RXTX with independent modules that communicate with RXTX solely through the Sun Microsystems CommAPI interface version 2, regardless of the license terms of these independent modules, and to copy and distribute the resulting combined work under terms of your choice, provided that every copy of the combined work is accompanied by a complete copy of the source code of RXTX (the version of RXTX used to produce the combined work), being distributed under the terms of the GNU Lesser General Public License plus this exception. An independent module is a module which is not derived from or based on RXTX.

Note that people who make modified versions of RXTX are not obligated to grant this special exception for their modified versions; it is their choice whether to do so. The GNU Lesser General Public License gives permission to release a modified version without this exception; this exception also makes it possible to release a modified version which carries forward this exception.

You should have received a copy of the GNU Lesser General Public License along with this library; if not, write to the Free Foundation, Inc., 59 Temple Place, Suite 330, Boston, MA 02111-1307 USA. All trademarks belong to their respective owners.

# 13 Appendix B: Specifications

## 13.1 Specifications of the USB sensor

- Connection type: USB
- Power supply: through USB.
- Accuracy:  $\bar{y}(\lambda) = v(\lambda)$
- Maximum deviation from reference sensor (Minolta CA 210) on reference color and grayscale LCD displays:
  - ✓  $\Delta Y \leq 5\%$  from 0.2 to 1000 cd/m<sup>2</sup> (typical)
  - ✓  $\Delta x \leq 0.005$  from 1 to 1000 cd/m<sup>2</sup> (typical)
  - ✓  $\Delta y \leq 0.005$  from 1 to 1000 cd/m<sup>2</sup> (typical)
  - ✓  $\Delta x \leq 0.01$  from 0.2 to 1 cd/m<sup>2</sup> (typical)
  - ✓  $\Delta y \leq 0.01$  from 0.2 to 1 cd/m<sup>2</sup> (typical)
- Repeatability over 5 measurements:
  - Y:
    - ✓ 0.01 cd/m<sup>2</sup> in the range of 0.2 to 1 cd/m<sup>2</sup> (typical)
    - ✓ 0.05 cd/m<sup>2</sup> in the range of 1 cd/m<sup>2</sup> to 10 cd/m<sup>2</sup> (typical)
    - ✓ 0.1 cd/m<sup>2</sup> in the range of 10 cd/m<sup>2</sup> to 100 cd/m<sup>2</sup> (typical)
    - ✓ 1 cd/m<sup>2</sup> in the range of 100 cd/m<sup>2</sup> to 1000 cd/m<sup>2</sup> (typical)
  - x, y :
    - ✓ 0.005 in the range of 0.2 cd/m<sup>2</sup> to 0.5 cd/m<sup>2</sup> (typical)
    - ✓ 0.003 in the range of 0.50 cd/m<sup>2</sup> to 1000 cd/m<sup>2</sup> (typical)

## 13.2 DCMS specifications

For a typical target of 6500 K, gamma 2.35 and luminance 150 cd/m<sup>2</sup>, the user may expect an accuracy of

- ✓  $\Delta E_{uv} \leq 2$  for MDPW476
- ✓  $\Delta E_{uv} \leq 3$  for ADP361 and LCx-4x
- ✓ gamma 2.35  $\pm$  0.025

for gray levels down to 20% drive level.

These values are not binding as results may vary depending specified targets, display and display type.

# 14 Table of contents

1	About .....	3
1.1	Contents of the user manual .....	3
1.2	Change record .....	3
1.3	Notation convention .....	4
2	Important notice .....	5
2.1	Commercial In Confidence .....	5
2.2	Copyright .....	5
2.3	Open Source Software .....	5
2.4	Embedded pictures .....	5
2.5	Product safety notice .....	5
2.6	Software License Terms .....	6
3	About DCMS .....	10
3.1	What is DCMS? .....	10
3.2	Why using DCMS? .....	10
3.3	Key components .....	11
3.3.1	Calibration system .....	11
3.3.2	Display .....	12
3.3.3	Sensor .....	12
3.4	Package content .....	13
4	Installing DCMS .....	14
4.1	System requirements .....	14
4.2	Before installing DCMS .....	14
4.3	Installing DCMS .....	14
4.4	Uninstalling DCMS .....	19
5	Connecting the components .....	21
5.1	Optical sensor connection .....	21
5.2	Calibration system connection .....	24
6	Running the DCMS application .....	26
6.1	User Interface .....	26
6.1.1	DCMS User Interface .....	26
6.1.2	Connection panel .....	27
6.1.3	Sensor panel .....	27
6.1.4	Calibration target panel .....	29
6.2	Using DCMS .....	35
6.2.1	Proper use and general remarks .....	35
6.2.2	Calibrate a display .....	35
6.2.3	Match displays .....	37
6.2.4	Single measurement .....	39
7	Display connection and calibration options .....	41
7.1	ADP361 .....	41
7.2	LCx-4x .....	42
7.3	MDP471 (ISIS) .....	43
7.4	MDPW476 .....	45
7.5	NSL-4601 .....	47
8	Calibration .....	50
9	Verification .....	52
10	The log file .....	54

11	List of abbreviations .....	55
12	Appendix A: Open Source Software.....	57
13	Appendix B: Specifications .....	58
13.1	Specifications of the USB sensor.....	58
13.2	DCMS specifications.....	58

## 14.1 List of figures

Figure 1:	External calibration system .....	11
Figure 2:	The input source as part of the calibration system .....	11
Figure 3:	Sensor head removal .....	21
Figure 4:	Calibration system connection .....	24
Figure 5:	DCMS User Interface .....	26
Figure 6:	Sensor Calibration Dialog .....	27
Figure 7:	Single measurement .....	28
Figure 8:	External test patch .....	29
Figure 9:	Internal test patch .....	29
Figure 10:	Additive colors .....	30
Figure 11:	Color temperature examples .....	30
Figure 12:	Color temperature scale* .....	30
Figure 13:	CIE1931xyY color space .....	31
Figure 14:	Match of white and mid gray color .....	32
Figure 15:	Mismatch of white and gray color .....	32
Figure 16:	The gamma power-law .....	32
Figure 17:	Gamma 4 .....	33
Figure 18:	Gamma 2 .....	33
Figure 19:	Gamma 1 .....	33
Figure 20:	Barco DCMS – Calibration window .....	34
Figure 21:	Display calibration .....	35
Figure 22:	Matching of displays .....	37
Figure 23:	DCMS with ADP361 .....	41
Figure 24:	DCMS with LCx-4x .....	42
Figure 25:	DCMS with MDP471 .....	43
Figure 26:	DCMS with MDPW476 .....	45
Figure 27:	MDPW LUT selection .....	46
Figure 28:	DCMS with NSL-4601 .....	47
Figure 29:	Monitor ID mapping table .....	48
Figure 30:	Calibration NSL-4601 completed .....	48
Figure 31:	Calibration steps flowchart .....	50

## 14.2 List of tables

Table 1:	Color temperature versus brightness .....	49
Table 2:	Verification table .....	52
Table 3:	Verification table explanation .....	53

*End of document*

## Free Manuals Download Website

<http://myh66.com>

<http://usermanuals.us>

<http://www.somanuals.com>

<http://www.4manuals.cc>

<http://www.manual-lib.com>

<http://www.404manual.com>

<http://www.luxmanual.com>

<http://aubethermostatmanual.com>

Golf course search by state

<http://golfingnear.com>

Email search by domain

<http://emailbydomain.com>

Auto manuals search

<http://auto.somanuals.com>

TV manuals search

<http://tv.somanuals.com>