



# Ampli500<sup>®</sup>



**Caution:** This telephone product is specifically designed for people with hearing difficulties and is provided with a handset earpiece volume control for individual requirements. Due care must be taken by all users that the handset earpiece volume control is set to the lowest level acceptable by each respective user. Care should therefore be exercised to ensure that any new users are familiar with this requirement.

English

Español Français

# CONTENTS

CONTENTS	2	VOLUME AND TONE	21
INTRODUCTION		Sending Level Adjustment	21
Unpacking the Ampli500™	3	Receiving Level Adjustment	21
DESCRIPTION	4-8	Handset Amplifier	21
INSTALLATION		PHONEBOOK	22
Setting Up	9	Create a Phonebook Entry	22
SETTINGS	10	Create and Entry from the	
Ringer Volume	10	Caller Display Log	23
Ringer Tones	10	Make a Call from a Phonebook Entry	23
Alarm & Ringer Settings	10	Delete a Phonebook entry	24
Dialing Mode	11	Modify a Phonebook Entry	24
Date & Time Setting	11	CALLER DISPLAY	25
Language Setting	12	Use a Number from the	
Screen Contrast Adjustment	13	Caller Display log	25
Auto On Hook	13	Delete a number from the	
Auto screen Sizing	14	Caller Display Log	26
Area & LDS Code Setting	15	Delete all Numbers from the	
USING THE TELEPHONE	16	Caller Display Log	26
Answering a Call	16	Call Waiting	26
Making a Call	16	NEW CALL/MESSAGE LED	27
Switching between Speakerphone		Message Waiting LED	28
and Handset	17	TROUBLESHOOTING	29
Call Time Counter	17	SAFETY INFORMATION	31
Last Number Redial	17	GUARANTEE	32
Recalling Recent Numbers	17		
Deleting Recent Numbers	18		
Mute function	18		
Flash	18		
One Touch Memory Keys	18		
OPTIONAL ACCESSORIES	20		
Headset	20		
Neckloop	20		
Shaker	20		

# INTRODUCTION

Congratulations on purchasing your Geemarc Ampli500™ telephone.

This is a multifunction telephone which offers features such as hands free use, phonebook and caller display\*. It also has the added benefits of outgoing voice amplification and receiving amplification.

It is important that you read the instructions below in order to use your Geemarc Ampli500™ to its full potential. Keep this user guide in a safe place for future reference. With any questions please contact customer service at 1-888-864-2446 or [RMA@sonicalert.info](mailto:RMA@sonicalert.info).

\*To use the caller display facility, you must subscribe to your network provider's caller display service.

**CAUTION: THIS TELEPHONE PRODUCT IS SPECIFICALLY DESIGNED FOR PEOPLE WITH HEARING DIFFICULTIES AND IS PROVIDED WITH AN AMPLIFIED HANDSET.**

## Unpacking the Ampli500™



When unpacking the telephone, you should find the following in the box:

- 1 Telephone base unit
- 1 Telephone handset
- 1 Handset curly cord
- 1 8 inch short telephone line
- 1 6ft telephone line cord
- 1 AC Adapter
- 1 Screwdriver
- 1 Wall plate adapter
- 1 User guide

# DESCRIPTION




# DESCRIPTION



# DESCRIPTION

## Screen example (On hook position)



OK button is used to confirm action or a choice on the screen.

▲ ▼ “UP / DOWN arrow” buttons are used to highlight or scroll through the menu.

C button is used to cancel or confirm an action/entry or clear a stored Caller ID call.

## DESCRIPTION

Your Amp500™ Phone has a screen with different icons, which help you to use the telephone



The telephone is supplied by battery + battery indicator



The telephone is powered by AC adapter.



Memory storing is required.



Screen contrast level.



“Mute” is activated.



Speakerphone is activated.



Precede the CID call record number.



To show new records.

# DESCRIPTION

## Connection



(\* ) Classified TNV-3 according to EN60950 standard.

(\*\* ) Classified „hazardous voltage“ according to EN60950 standard.

**Important: Only Shaker, Neckloop, Headset and Adaptor with Geemarc reference can be connected to the sockets.**



# INSTALLATION

## Setting Up

Connect the curly cord and the handset.

Remove the screw of the battery compartment, using the screwdriver provided, and open the battery compartment. Insert 4 good quality alkaline batteries (AAA size, 1.5V type, not supplied) into the battery compartment located on the base of the telephone.

**Note :** Do not use rechargeable batteries. Do not put batteries into the household rubbish, pay attention to battery regulations regarding battery disposal.  
Replace the battery cover.

## Desk Installation

Connect the 6ft line cord to the socket located at rear of telephone then connect the plug to the telephone wall socket.


Connect the power supply to the socket located on the back of the telephone then plug into the wall socket. Please choose a socket near your telephone socket to enable you to unplug the mains power quickly in case of a problem.

Lift the handset, if there is a dial tone, the telephone is ready for use. If there is no dial tone, check that everything is plugged in as described above.

**Note: Mains power is required to power the multifunction features of the telephone. Fit 4 x AAA sized batteries in order to save and use stored numbers in case of power failure. Without mains or battery power, the telephone can still be used to answer calls or to manually dial out calls but all stored numbers will be lost.**

## Wall Mount Installation

1. Slide the handset hanger bar out and rotate it 180 degrees. Slide it back into place. The hook will now point up. This ensures the handset will not fall out of the cradle when it is mounted on the wall.
2. Plug the 8in short phone cord into the telephone jack on the Ampli500™, run the phone cord underneath the wall plate adapter and through the opening.
3. Insert the Wall Plate adapter. Hold the wall plate adapter so that the raised mounting pins line up with the phone's keyholes on the Ampli500™. Once the holes are lined up, push the wall plate adapter in towards the Ampli500™ and then slide the plate upwards until it is securely seated on the phone.
4. Once the Wall Plate adapter is installed, plug the 8 inch short cord into the jack on the wall.
5. Holding the phone slightly above the mounting screws on the wall jack, push the phone against the mounting screws so they are hooked into upper and lower key-hole slots on the back of the phone. Slowly slide the phone down until it snaps into place.
6. Connect the handset curly cord to the handset and the other end to the telephone.
7. Lift the handset, if there is a dial tone, the telephone is ready for use. If there is no dial tone, check that everything is plugged in as described above.




# SETTINGS

## Ringer Volume

The ringer volume can be adjusted as high as 95 +dB.

A switch located on right side allows the ringer to be switched to OFF, LOW, HIGH. Your telephone comes with the ringer switched to 'High'.

## Ringer Tones



## Alarm and Ringer Settings


A 3-position switch located on the right side of the phone allows the settings below:



Ringer and strobe flash are activated (Normal Setting)

The other two settings are only effective when an optional vibrator pad has been added.



 Ringer, shaker and strobe flash are activated



Ringer and shaker are activated

# SETTINGS

## Dialing Mode

The Tone/Pulse switch is found on the base of the phone near the battery compartment. Set the switch to (T) if you have Touch Tone service. The Tone/Pulse switch is pre-set to Tone position (T).


If your phone does not dial out, it is possibly being used from an older private switchboard. In this case, the Tone/Pulse switch must be moved to Pulse (P)

## Date & Time Setting

If you subscribe to Caller I.D. service via your telephone service provider, these parameters will usually be set automatically when you receive your first call.

### **Date Format**


You can choose between Day-Month or Month-Day.



# SETTINGS

## Hour Format


You can choose between 12 hours or 24 hours.



## Language Setting

You can choose between 16 different languages including English and French. The default is American English (English 1).


**Note:** Do not change to a language you cannot read or understand as it will be difficult to return to a language on screen that you do understand!



# SETTINGS


## Screen Contrast Adjustment

Your screen has a backlight that automatically switches off after 15 seconds of inactivity. To decrease or increase the contrast of the backlight:



## Auto On Hook


When this feature is ON, if you leave the telephone off the hook by mistake, the Ampli500™ will automatically disconnect the line after 5-8 minutes and your telephone will be ready to take another incoming call.



# SETTINGS

## Auto Screen Sizing –DIM:X-->x


The phone can show a maximum of 8 digits at full size. You can choose to see the complete number (up to a maximum of 18 digits) by selecting 'Auto Screen Sizing'. When activated, the phone will display at full size but automatically switch to display small numbers when there is a number with more than 8 digits. The default on arrival is activated. However, if you ever need to change this setting, see below:



# SETTINGS

## Area and Long Distance Code Setting

When you use your phone for the first time or if you have moved to a different area code, you need to enter your area code into the phone. This assures that when dialing a call out of memory or your phonebook, the number is formatted properly. To set your area code & LDS Code:



# USING THE TELEPHONE

## Answering A Call

During an incoming call, your phone will ring and then the number and/or the name of the caller will be displayed on the screen, if you have subscribed to a Caller I.D. service via your telephone provider.

### Using the handset

- Lift the handset and speak.
- To end the call, replace the handset.

### SPEAKER PHONE (hands free)

- Your Ampli500™ allows you to speak hands free. When the phone rings, press **SPKR**. The caller's voice will be heard through the loudspeaker and they will hear your voice via the microphone near the **SPKR** button. The speakerphone symbol will appear on the display.
- To end the call, press **SPKR** again.

### Using the inductive coupler for hearing aids

- Turn your hearing aid to the 'T' setting. (Check with your audiologist or hearing aid manufacture to determine whether or not your hearing aid/s is equipped with T-Coils for the T-Setting).
- Use the phone handset as usual but the inductive coupler will feed the incoming call directly into your hearing aid.

**Note:** users of digital hearing aids must make sure their aids have been programmed for use with the 'T' setting.(check with your audiologist or hearing aid manufacturer to determine whether or not your hearing aid/s are equipped with T-Coils for the T-Setting)

## Making A Call

### Using the handset

- Lift the handset and when you hear the dialling tone, dial the number on the keypad.
- To end the call, replace the handset.

### Pre-dialling

- Key in the required number without lifting the handset.
- Use the **C** key if necessary to cancel and correct a wrong digit.
- Lift the handset and the number will be automatically dialled.
- To end the call, replace the handset.

### SPEAKER PHONE (hands free)

- Key in the required number without lifting the handset.
- Use the **C** key if necessary to cancel and correct a wrong digit.
- Press **SPKR** and the number will be automatically dialled. The speakerphone symbol will appear on the display.
- When the call is answered, simply speak into the microphone located on the base
- You can adjust the speaker level with the slide button located on the left side of your telephone.
- To end the call, press the "SPKR" **SPEAKER PHONE (hands free)** button again..



# USING THE TELEPHONE

## Switching between Speakerphone and Handset

To stop using speakerphone during a call, simply lift the handset.

To activate speakerphone during a call, press **SPKR** and replace the headset on the cradle at the same time

Hands free calls can also be made using an optional Geemarc headset –see Optional Accessories chapter

## Call Time Counter

The call time counter will be displayed a few seconds after the phone number is connected and starts ringing.

To stop the time counter, press **C**

If you press **C** twice, the time counter will be reset to 0.

## Last Number Redial

Lift the handset or press **SPKR** or press Headset  
(if headset plugged in)

Await the dial tone and press **REDIAL**

The previously dialed number will be automatically redialled

## Recalling Recent Numbers

You can call back any of the last 30 numbers dialed

Press **OK**


Screen will display '**DIALLED**', the last number dialed, total number of numbers dialed plus time and date of displayed number

Use **▲** or **▼** to to select the number you wish to dial

Lift the handset or press **SPKR** or press Headset  
(if headset plugged in), the number will be automatically dialed


# USING THE TELEPHONE

## Deleting Recent Numbers



## Mute Function

If you wish, you can speak privately to someone else in the room without your caller overhearing. You will still be able to hear the caller but they cannot hear you during the mute operation.



## Flash

This is used with the new services provided on digital exchanges e.g. the flash button can be used to disconnect a call and re-establish a dial tone or switch to another caller when using services such as Call Waiting.

For details please contact your network operator.

## One Touch Memory Keys

Your phone has a total of 11 one-touch memory keys.


The 3 priority P1, P2, P3(CS) memories allow you to store priority numbers in case of emergencies, such as:

- Local police station,
- Your doctor
- CS customer Service

There are also 8 one-touch memories for frequently used numbers.


# USING THE TELEPHONE

## To store a number



**Note:** Numbers will be lost after a few hours if you disconnect your phone.

## Dialling stored Numbers




There is an index card located under the base of your phone for you to make a list of which numbers are stored under which one touch memory key.

# OPTIONAL ACCESSORIES

## Headset


You can also use your Ampli500™ in hands free mode with an optional Geemarc headset, this is ideal for use with a computer screen.

### To use the headset



You can adjust the volume and tone level of calls received with the VOLUME and TONE slide buttons located on the top of the phone (if function is active) and the sending level with the 3 position switch located on the right side.

## Neckloop



**Note:** The Ampli500™ must be connected to the mains in order to use a neck loop.

## Shaker

If you would like to be alerted to incoming phone calls by a vibrating shaker, you can use an optional Sonic Alert bed shaker model ss12v.

Connect the shaker into the shaker socket located at the back of your phone. When an incoming call is received, the shaker will be activated and vibrate while the phone rings.

One of the shaker settings must be selected, using the 3 position switch located on the right side of the phone.

(See Settings chapter)

Important: only connect a Geemarc shaker, any other shaker model can damage your telephone.

# VOLUME AND TONE

## Sending Level Adjustment

If your voice is quiet or the person you have called has difficulty hearing, you can increase the sending level volume with the switch located on the right side. You can also decrease the sending volume if required. There are three levels, Min, Normal and Max (outgoing volume up to 12dB).

## Receiving Volume and Tone Adjustment

You can increase the receiving volume and adjust the receiving tone according to your own hearing requirements. The AMPLI button allows extra receiving and tone functions to be activated during the conversation in an ad hoc fashion. When the function is activated, the red light to the left of the button is lit. Use the TONE and VOLUME slider buttons to adjust the level according to your own hearing requirements. When **AMPLI** is on, the receiving volume adjustment is 15- 30dB. The tone control allows you to either boost low or high frequency range. When **AMPLI** is off, the receiving volume adjustment is 0 - 15dB. The tone control cannot be used.

When the call is finished, the telephone will reset to standard mode i.e. **AMPLI** turned off.

### **Amplification over ride switch**

In order to override the standard setting i.e. **AMPLI** switched off, turn the amplification over ride switch located on the rear panel of the telephone, to On. (use the tip of a pencil or small screwdriver if necessary).

**ON:** When the amplification over ride switch is ON, the extra receiveing volume and tone functions are automatically activated each time you use the phone. You do not need to use the **AMPLI** button. (during any call, you can temporarily switch off and on the amplification override by pressing the **AMPLI** button.)

**OFF:** When the switch is switched OFF, you will need to activate the function by pressing the **AMPLI** button in order to get use of the extra volume and tone functions.

This override is useful if the person using the telephone will always need the amplification setting.

## Handset Amplifier

In both cases, you can also boost the reception volume by an extra 12 dB by holding down the handset amplifier button located on the handset.


**Note:** The button must be released while talking or it will act as a mute button.

**Caution: Exposure to high amplification levels to people without hearing loss may cause damage to hearing.**

# PHONEBOOK

## Create a Phonebook Entry

You can store 99 entries into the phone book with the name and phone numbers of your correspondent.



You can use this table or the letters printed on the keypad to help you create a phone book entry.

Key	Characters
1	1. ? - " ' < 1
2	ABCabc2ÆØ
3	DEFdef3 4 GHghi4
5	J K L j k l 5
6	MNOmno6ÖØ
7	P Q R S p q r s7
8	T U V t u v 8 U
9	WXYZwxyz9
10	0+&@/\$%

Example: To input the letter 'b' press the **2** key 5 times. The letter 'b' will show in the display

To insert a 'space', press on the **OK** key twice.

If the next letter is the same as the previous letter (or located on the same key), press first on the **OK** key to move and select the desired letter.


# PHONEBOOK

Use the **C** key to delete a wrong number or letter.


**Note** : If you wish to insert a pause in the stored number, press the **REDIAL** key whilst entering (a P will be displayed on the screen).

## Create an Entry from the Caller I.D. Display Log

You can store 99 entries into the phone book with the name and phone numbers of your correspondent.




## Make a call from a Phonebook Entry




# PHONEBOOK

## Delete a Phonebook Entry



## Modify a Phonebook Entry





# CALLER I.D. DISPLAY LOG

The features described in this section are only available if you subscribe to a Caller Display service from your network provider.

Caller Display means you can see who is calling on your LCD display (provided the number is not withheld, unavailable or is an international call). When your telephone rings, the time, date, the caller's phone number and name and a count number (1-99) is displayed on the LCD display. The most recent call will have the highest count.

Every received call that is not answered is stored in the call log and the new calls count together with total calls count displayed on the standby screen increases by one. The LCD display will return to standby 10 seconds after the last ringing burst if the call is not answered.

This new call will also activate the new call indicator to blink. This will continue to blink until the new calls are reviewed and saved/deleted in the call log.

**Note:** Area code and LDS code must be set correctly (see Settings chapter) for the telephone to recognize and dial calls from the Caller Display log.


If the caller does not wish to show his phone number or his name, the screen will display '**PRIVATE**'.

If the number is unknown from the telephone centre (international incoming call for example), the screen will display '**OUT OF AREA**'.


To review the new calls, use the **←** or **→** keys. When all new calls have been reviewed, the new call LED will disappear and the screen will come back to idle state after a few seconds. The total number of incoming call entries on the incoming call log will still be displayed on the screen. When the incoming call log is empty (if records have been deleted for example) the screen will display '**00**'.

## Use a Number from the Caller I.D. Display LOG

There are two ways to use numbers from the call log. Firstly, **↵** the number and then lift the handset:




Alternatively you can lift the handset and then **↵** the number:




Note: You cannot use either of these methods to recall a number that displays **PRIVATE** or **OUT OF AREA**.

# CALLER I.D. DISPLAY

## Delete a Number from the Caller I.D. Display Log




## Delete all Numbers from the Caller I.D. Display Log



## Call Waiting

With this service, you are made aware that a new call is coming through when you are in conversation. You can take this call and switch between callers. To use this feature, you must have subscribed to this service with your network provider.



# CALLER I.D. DISPLAY

## **Call Waiting in Practice**

When you are in conversation with caller no.1, their name and number will be displayed on the screen.

A beep tone will be heard when a second call is incoming and their name and number will be displayed on the screen. If you chose to answer the call, press FLASH or the phone hook button to transfer.

Note : If you don't answer the second call, the number will be recorded into the incoming call log and you will be able to call them back later.


For more information about subscribing to Call Waiting, contact your network provider.

# NEW CALL/MESSAGE LED

The red LED on the front of the phone will flash to indicate missed calls.


To cancel the flashing light, review your new calls. Simply use ▲ or ▼ to access the log and scroll through all the new calls.

Alternatively, you can delete the log of unreviewed calls



If you do not subscribe to caller display or message waiting service simply press ▲ or ▼ once.

To turn off this feature, set the 'MSG' switch on the base of the phone to 'OFF'. New calls will still be indicated with 'N' on the screen but the LED will not flash.



## Message Waiting LED

If you are connected to certain type of switchboards, you may get a message waiting/voicemail service. 'MSG WAITING' will be displayed on screen and the red LED will be lit continuously. You must retrieve your voicemail in order to turn the red light off.

To turn the message waiting feature off there is a switch on the bottom of the phone with a symbol of an envelope, move this switch to the OFF position.

# TROUBLESHOOTING

## No display

Make sure that the power adaptor is plugged in correctly and is not damaged

Make sure the batteries are inserted correctly

Make sure that the batteries are fully charged

## Telephone does not ring

Make sure that the telephone line cord is plugged in correctly and is not damaged

Connect another telephone to the telephone PTT socket or move your telephone to another socket to track if the failure is coming from your telephone or your PTT socket

There may be too many communication devices connected to a single telephone line e.g. another telephone, modem or facsimile machine. Contact your network provider for help in calculating the limit for your home or business line

Check that the button ringer on the right side is not set to OFF position.

Check that the phone is in the ON HOOK position

## No dial tone

Make sure that the telephone line cord is plugged in correctly and is not damaged

Connect another telephone to the telephone PTT socket or move your telephone to another socket to track if the failure is coming from your telephone or your PTT socket

There may be too many communication devices connected to a single telephone line e.g. another telephone, modem or facsimile machine. Contact your network provider for help in calculating the limit for your home or business line

## No number is displayed when the telephone rings:

Please contact your telephone network provider to obtain the caller display service or ensure that it has been activated if you have subscribed already

Even if you have subscribed to caller display, if the caller does not want to show their number, the display will show Private number or if the call is from an area not providing caller display, Unavailable number will be displayed

It could be a call from a PBX. Caller display may not work properly if the telephone is connected to a PBX

## No response on LCD display

If you have experienced a power failure the batteries may become weak. If the power is restored and the LCD display still does not respond then replace the batteries with new ones.

### NOTE:

Contact customer service Monday-Friday 10:00AM EST-5:00PM EST. at 1-888-864-2446 or via Email at [RMA@sonicalert.info](mailto:RMA@sonicalert.info). Allow up to 1 business day to receive a response via Email.

# SAFETY INFORMATION

## General

Only use the correct Geemarc authorized power supply with the product.

Do not open the unit (except to replace the batteries). Contact the helpline for all repairs.

Never dispose of batteries in a fire. There is a serious risk of explosion and/ or the release of toxic chemicals.

## Cleaning

Clean the telephone with a soft cloth. Keep the base and handset charging contacts clean. Never use polishes or cleaning agents - they could damage the finish or the electrics inside.

## Environmental

Do not expose to direct sunlight.

Always ensure there is a free flow of air over the surfaces of the telephone.

Do not place any part of your product in water and do not use it in damp or humid conditions e.g. bathrooms.

Do not expose your product to fire or other hazardous conditions.

Unplug your phone from the telephone socket during thunderstorms because it could damage your telephone. Your guarantee does not cover thunderstorm damage.

Do not use the telephone if there is a gas leak. Any spark that might be created by the electrical contacts found in the telephone could ignite a large amount of leaking gas.

The telephone is designed for working within a temperature range of 5°C to 45°C.

# REGULATORY COMPLIANCE

## Part 68 of FCC rules information

- This equipment complies with Part 68 of the FCC rules and the requirements adopted by the ACTA. On the bottom of this equipment is a label that contains, among other information, a product identifier in the format US:AAAEQ##TXXXX. If requested, this number must be provided to the telephone company.
- All applicable certification jacks Universal Service Order Codes (USOC) for the equipment is provided (i.e. RJ11C) in the package with each piece of approved terminal equipment.
- A plug and jack used to connect this equipment to the premises wiring and telephone network must comply with the applicable FCC Part 68 rules and requirement adopted by the ACTA. A compliant telephone cord and modular plug is provided with this product. It is designed to be connected to a compatible modular jack that is also compliant. See installation instructions for details.
- The REN is used to determine the number of devices that may be connected to the telephone line. Excessive RENs on a telephone line may result in the devices not ringing in response to an incoming call. In most but not all areas, the sum of RENs should not exceed five (5.00). To be certain of the number of devices that may be connected to a line, as determined by the total RENs, contact the local telephone company service provider. For products approved after July 23, 2001, the REN for this product is part of the product identifier that has the format US:AAAEQ##TXXXX. The digits represented by ## are the REN without a decimal point (e.g., 03 is a REN of 0.3). For earlier products, the REN is separately shown on the label.
- If this equipment causes harm to the telephone network, the telephone company will notify you in advance that temporary discontinuance of service may be required. But if advance notice isn't practical, the telephone company will notify the customer as soon as possible. Also, you will be advised of your rights to file a complaint with the FCC if you believe it is necessary.
- The telephone company may make changes in its facilities, equipment, operation or procedures that could affect the operation of the equipment. If this happens the telephone company will provide advance notice in order for you to make necessary modifications to maintain uninterrupted service.
- Contact customer service Monday-Friday 10:00AM EST-5:00PM EST. at **1-888-864-2446** or via Email at [RMA@sonicalert.info](mailto:RMA@sonicalert.info). Allow up to 1 business day to receive a response via Email.

If the equipment is causing harm to the telephone network, the telephone company may request that you disconnect the equipment until the problem is resolved.

- Please follow instructions for repairing if any; otherwise do not alter or repair any part of device except as specified.

Contact customer service Monday-Friday 10:00AM EST-5:00PM EST. at **1-888-864-2446** or via Email at [RMA@sonicalert.info](mailto:RMA@sonicalert.info). Allow up to 1 business day to receive a response via Email.

(To setup the auto dial feature for Sonic Alert / Geemarc customer service the Green CS key must be programmed with the following phone number **1-888-864-2446**. Follow the directions for the emergency key programming on page **18** to program the CS customer Service Key.)

**Note:** For California Program Customer Service please program the phone with the following numbers:

**English:** 1-800-806-1191

**Spanish:** 1-800-949-5650

**TTY:** 1-800-806-4474

**Mandarin:** 1-866-324-8747

**Cantonese:** 1-866-324-8754

**Hmong:** 1-866-880-339

- Connection to party line service is subject to state tariffs. Contact the state public utility commission, public service commission or corporation for information.
- If your home has specially wired alarm equipment connected to the telephone line, ensure the installation of this telephone equipment does not disable your alarm equipment. If you have questions about what will disable alarm equipment, consult your telephone company or a qualified installer.
- This equipment is hearing aid compatible.

# GUARANTEE

Sonic Alert warrants the phone against any defect in materials or workmanship for the period of one year from the date of purchase.

Should you experience a problem, contact our customer service department

Be sure to save your sales receipt as proof of purchase date should you need warranty service.

Within a period of one year from purchase date, Sonic Alert will repair or replace (our discretion) your phone at no cost, if a defect in materials or workmanship is found. If we elect to replace your phone, we may replace it with a new or reconditioned product of the same or similar design.

Repair or replacement will be warranted for a period of 90 days or the original time on the original warranty, whichever is longer.

The warranty does not cover accidents, negligence or breakage to any parts. This includes shipping damage, failure to follow instructions, misuse, fire, floods, use of incompatible accessories, Acts of God or failure in your phone service carrier's line service. The product must not be tampered with or taken apart by anyone who is not an authorized

Sonic Alert representative. Tampering with the phone will void ant written or implied warranties

Sonic Alert shall not be responsible for loss of time, inconvenience, property damage caused by your phone or any other accidental or consequential damages

Warranty service is available only with proof of purchase.

Simply send the Ampli500™ (postpaid) and a copy of your sales slip as proof of purchase to:  
Sonic Alert Inc., 1050 E. Maple Road, Troy MI 48083

**IMPORTANT : YOUR RECEIPT IS PART OF YOUR WARRANTY AND MUST BE RETAINED AND PRODUCED IN THE EVENT OF A WARRANTY CLAIM**



## SonicAlert

1050 East Maple Rd.

Troy, MI 48083

Sales: 248-577-5400 Fax: 248-577-5433

Customer Service & Warranty: 1-888-864-2446

[www.sonicalert.com](http://www.sonicalert.com)





## Free Manuals Download Website

<http://myh66.com>

<http://usermanuals.us>

<http://www.somanuals.com>

<http://www.4manuals.cc>

<http://www.manual-lib.com>

<http://www.404manual.com>

<http://www.luxmanual.com>

<http://aubethermostatmanual.com>

Golf course search by state

<http://golfingnear.com>

Email search by domain

<http://emailbydomain.com>

Auto manuals search

<http://auto.somanuals.com>

TV manuals search

<http://tv.somanuals.com>