


PRO PLUS® STRETCHER SCALE
MODEL 2000KL



OPERATION MANUAL

AA0339-0

PRO PLUS® STRETCHER SCALE MODEL 2000KL

Thank you for your purchase of this product.

Please read this manual carefully and keep it handy for ready reference.

TABLE OF CONTENTS

CAUTION AND WARNING.....	2
SPECIFICATIONS.....	3
ASSEMBLY INSTRUCTION.....	4
SET UP.....	10
OPERATION INSTRUCTIONS.....	11
MAINTENANCE.....	23
TROUBLESHOOTING.....	24
CALIBRATION PATH.....	25
EXPLODED VIEW OF SCALE.....	25
PARTS LIST.....	27
WARRANTY.....	28

CAUTION AND WARNING

To prevent injury and damage to your scale, please follow these instructions very carefully.

- To ensure proper operation of the ProPlus® Stretcher Scale, this operation manual should be reviewed carefully before operation. Keep this manual available for future reference, and also for use in the orientation of new personnel.
- If the scale becomes damaged, it should not be operated until properly serviced. All repairs should only be performed by authorized Healthometer® Professional service personnel.
- This scale is designed for static weighing of patients only. **No scale should be used for patient transfer. Healthometer® Professional / Pelstar LLC assumes no responsibility for patient injury or scale damage should this caution be ignored.**
- **Do not leave patients unattended while on the scale.**
- For accurate weighing the scale must be placed on a flat and stable surface with the locking screw locked properly and the weighing mechanism not touching the floor.
- For accurate weighing, verify proper operation according to the procedure described in this manual before each use.
- Do not use in the presence of flammable materials.
- Operating at other voltages and frequencies than specified could damage the equipment.
- If the “LOW BAT” indicator activates, for accurate weighing, replace the batteries or connect the scale to an AC power source as soon as possible.
- Failure to comply with the following will void scale warranty:
 - Do not gas sterilize or autoclave the scale.
 - Do not place liquids on top of the scale display console.
 - Do not exceed recommended weigh limit of 400 lb/181kg for this scale.
 - This scale contains delicate sensors. Do not bang into doors, elevators, etc., as damage may result.

SPECIFICATIONS

GENERAL

Health o meter's ProPlus® Stretcher Scale Model 2000KL uses highly sophisticated microprocessor technology. This scale is designed to provide accurate, reliable and repeatable weight measurements and features that make the weighing process simple, fast and convenient. The scale is set up to use motion-sensing weighing technology, to determine the actual weight of a moving patient. The scale may be changed to measure live weight; see page 19 for instructions on changing the scale setting.

The weight can be displayed in pounds (decimals, fractions of a lb or lb/oz) or in kilograms. The scale features a folding stretcher for easy mobility and a display assembly stand. The unit can be operated using its AC adapter or by 6-D cell batteries (not included).

SCALE SPECIFICATIONS

Capacity and Resolution	400 Lb x 0.2 Lb / ¼ lb / 4 oz (181 Kg x 0.1 Kg)
Power Requirements	Adapter 120VAC - 6VDC - 60Hz or 6 x D size batteries
Environmental	Operating temperatures: 50°F to 95°F (10°C to 35°C) Storage temperatures: 30°F to 125°F (0°C to 50°C) Humidity: 85%
Stretcher Dimensions	Closed: Length: 70 7/8" (180 cm) Width: 6" (15 cm) Weight: 16.5 Lb (7.5 Kg) Fully extended: Length: 70 7/8" (180 cm) Width: 29 1/2" (75 cm)

ASSEMBLY INSTRUCTIONS

Before Assembly

The Model 2000KL ProPlus® Stretcher Scale is shipped disassembled in two (2) separate cartons. Carefully inspect all cartons for shipping damage before unpacking. If damage is found, contact your shipper or Healthometer Professional immediately. Claims must be filed with the shipper as soon as possible after receipt of the damaged package.

Remove each assembly from the carton and unwrap packing material carefully to prevent scratching the unit's parts. It is recommended that two people perform the assembly.

The scale is shipped in two (2) separate cartons as follows:

A. Large carton containing:

- (1) Hydraulic pump assembly
- (1) Lifter mast
- (1) Base adjustment lever
- (1) Lifter base with casters

B. Small carton containing:

- | | | |
|------------------------|-------------------------------|---------------------------------|
| (1) Lifter boom | (1) Stretcher storage tray | (4) Stretcher bars |
| (1) Display support | (1) Adapter holder | (2) Stretcher bar joint |
| (1) Display assembly | (1) AC adapter | (1) Hexagon screw with lock nut |
| (1) Load cell assembly | (1) Display rear cover | (2) Plastic washer |
| (2) Stretcher hoops | (6) Screws for display | (1) Velcro strap |
| (1) Stretcher sheet | (2) Shoulder screw and washer | |

Tools required: Phillips head screwdriver, adjustable wrenches, 7/16" Allen wrench.

Parts of Stretcher Scale




Figure 1

ASSEMBLY INSTRUCTIONS

Assembling the Patient Lifter

1. Place the Stretcher Storage Tray on the Lifter Base, as shown below, before attaching the Lifter Mast.





Figure 2

2. Remove black plastic protective caps from Lifter Mast. Insert the bottom of the Lifter Mast and Hydraulic Pump assembly into the mast sleeve (through the Stretcher Tray Holder) on the Lifter Base. The notch on the end of the mast fits over the round tube at the bottom of the Mast Sleeve. The mast must lock into position in the mast sleeve, making it impossible to rotate. Check that the mast is locked into position and that the hole in the mast is aligned with the threaded hole in the mast sleeve.
3. Insert the threaded locking device (the plastic knob and stud chained to the base) into the threaded hole in the mast sleeve. Tighten to lock the mast to the base. Keep mast and base locked at all times except when removing mast from base for storage.

ASSEMBLY INSTRUCTIONS

Attaching the Display and Scale to the Patient Lifter



ASSEMBLY INSTRUCTIONS

4. Connect the Lifter Boom and the Display Support to the Lifter Mast using the NC $\frac{3}{8}$ -20x2 $\frac{1}{2}$ " hex screw and tighten the self lock nut on its other side.

Note: This screw is the axis of the Lifter Boom, so it is imperative that the three parts – Lifter Boom, Lifter Mast and Display Support – be assembled in the proper order.

Place the Plastic Washers on either side of the Lifter Boom screw hole to prevent friction. With the Boom hooks facing down, place the Lifter Boom and Washers inside the two sides of the Lifter Mast and align the holes. Place the Display Support over the Lifter Mast and Lifter Boom and align holes. Thread the screw through all three parts and tighten the self lock nut on the other side well.

5. Insert the shaft end of the Pump into the eared bracket on the Lifter Boom and align the holes. Slip the 1-3/4" nut through the bracket and the Pump shaft hole. Tighten with locknut.
6. Insert the Base Adjustment Lever into its socket at the back of the Lifter Base. Check to see that the bump on the outside of the socket fits securely in the notched adjusting plate. To adjust the base legs, pull back on the adjusting lever to unlock. Move the lever to the right to widen the base.
7. A wheel lock is provided on the Lifter Base for parking only. To lock caster, step on cam lever on side of caster. To unlock, step on higher cam lever. NOTE: Never lock wheel lock or block wheels when lifting someone. Wheel lock is provided for parking the Lifter only,
8. Check the operation of the Hydraulic Pump by pumping the handle to elevate the Boom. Make certain the control valve knob – located on the pump near the handle) is fully tightened (clockwise).
9. To lower the Boom, slowly turn the control valve knob counterclockwise. The rate of descent can be controlled by how much the knob is turned. The Boom may not lower readily with weight. To test, simply apply downward force to the Boom.

Note: If bolts, mounting pump or Lifter Mast are too tight, the Boom may not lower properly.

ASSEMBLY INSTRUCTIONS




Figure 4

10. Position the Display Assembly close to the Display Support. Connect the load cell cable connector to the load cell port on the rear of the Display Assembly and insert and secure the cable into the center channel in the display.

Note: Do not attach the power cord at this point.

11. Slide the display assembly onto the display support by inserting the two metal supports into the two slots in the rear display assembly. Tuck the excess load cell cable in the Lifter Boom tube.

Note: The load cell cable should now be secured in the channel between the display assembly and the display support.

12. Insert 4 screws into the display support and tighten.

13. Connect power cord to the power port in the display assembly and insert and secure the cable into the exposed portion of the 2nd channel from the left in the display.

14. Insert the power cord into the retaining slot on the rear cover. Align the rear cover tabs with the slots on the rear of the display assembly. Apply pressure to snap the rear cover into place and secure with 2 screws.

ASSEMBLY INSTRUCTIONS




Figure 5

15. Hang the Load Cell Assembly on the two hooks at the end of the Lifter Boom, using the loops on top of the Load Cell Assembly.
16. Connect the end of the Load Cell Cable coming out of the center of the Load Cell Assembly to its port in the Lifter Boom.
17. Use the two washers and shoulder screws to connect the two Stretcher Hoops to the Load Cell Assembly and tighten well.
18. Connect the two right stretcher bars by screwing them onto one Stretcher Bar Connector, in order to make one long Stretcher Bar.
19. Repeat step 9 with the left Stretcher Bar.
20. Insert the two bars into the open side of the channel-pockets in the Stretcher Sheet and through the two hanging loops sewn into both sides of the Stretcher. The Stretcher Bars should not protrude from the Stretcher Sheet at either end. Make sure one end on both channel-pockets is sealed with Velcro.

Note: Make sure the bars are inserted into the two hanging loops sewn on both sides of the stretcher. This will ensure that the Stretcher Bars hang correctly from the hanging loops.

SET UP

Preparing the Scale for Use

1. Remove protective plastic film from keypad and display.
2. Place batteries in the battery holder (see “Replacing Batteries”).
3. Plug the scale’s AC adapter into the power source.
4. Press the ON/OFF button to turn the scale on. The display will show “Health o Meter Pro-Plus” and then “000Lb00oz”.^(a) If the Stretcher weight has previously been tared out, the display will show a negative weight equivalent to the weight of the Stretcher.
5. Hang the stretcher on the Stretcher Hoops to register a 0 weight on the display. Press the Tare key if necessary.
6. Place a weight [not exceeding 400 Lb (180 Kg)] on the stretcher. The display should read the correct weight.
7. Press the KG/LB button to select the weight mode (Lb/Kg).^{(b) *}
8. Press the Reweigh button; the scale should perform the weighing process again. The display should read “WEIGHING” until the weight calculation is complete and the weight is displayed.
9. Remove the weight from the stretcher, the scale will return to zero and the display will read “ZERO” on the left upper side of the screen along with “000Lb00oz”.
10. Disconnect the scale’s AC adapter from the power source, and the scale will switch to battery power.

Note: If the set up procedure failed, refer to the troubleshooting instructions. If the problem is not corrected, refer to qualified service personnel.

^(a) To adjust the display backlight and/or contrast, please refer to page 19.

^(b) To change the display mode in pounds to fractions or decimals, please refer to page 14.

Replacing Batteries




Figure 6. Replacing Batteries

(refer to the Parts List on page 25 for details on parts followed by (#))

1. Unplug the scale.
2. Remove the battery cover from the display assembly.
3. Disconnect battery holder cable connector from the scale-battery connector.
4. Carefully remove the battery holder by sliding it out of the display assembly.
5. Replace the batteries with new ones.**
6. Carefully slide the battery holder into the display assembly.
7. Plug the battery holder cable connector to the scale-battery connector.
8. Attach the battery cover to the display assembly and install the screw.

* To change the scale default units to kilograms, or to deactivate the LB/KG button, please contact Technical Support at 1 800 638-3722 or +1 708 598-9100.

** We recommend the use of EVEREADY Energizer® e2™ batteries.

OPERATION INSTRUCTIONS

Preweighing the Stretcher

1. Make sure the scale is free of load (nothing is hanging from it), including the stretcher.
2. Press the ON/OFF button to turn the scale on.
3. Wait until "000Lb00oz" and "ZERO" on the left side of the display will appear. The scale will remember the last amount that was tared, so if the weight of the stretcher had previously been tared, the scale will display a negative weight when it is turned on again.

Note: The 'ZERO' sign will show on the upper left side of the display only when the scale is free of load, including the stretcher.

4. If the weight of the stretcher is known, enter the TARE value according to instructions in the "Menu", para 02. If the stretcher's weight is unknown, place the stretcher only on the scale; press REWEIGH. Enter the Automatic TARE by pressing the TARE button. Display will read "000Lb00oz".

Note: Once the tare is set the "TARE" sign will appear in the upper right side of the display.

Placing Patient in Stretcher

1. After setting the stretcher's weight as tare, take the stretcher off the lifter. The display will now show a negative weight, equivalent to the weight of the stretcher.
2. Lay the stretcher on a bed or on any other stable surface you can use to lay the patient you wish to weigh.
3. Carefully move the patient to a lying position on the stretcher.
Note: lay the patient's center mass as close as you can to the center of the stretcher.
4. Wheel the lifter scale close to the stretcher and lower the boom with the hoops as far as 5-10" (12 cm -25 cm) from the patient's body.
5. Connect the four hanging loops of the stretcher to the four hooks at the end of each Stretcher Hoop, making sure that the loops are hanging freely from the four hooks.

Note: Make sure both hooks connecting the load cell assembly beam lie on the lifter boom. If only one of the hooks lies on the beam it can cause inaccurate readings of weight.

6. Tighten the hydraulic control valve knob by turning it clockwise. The control valve knob is located on the Pump near the Pump handle.
7. Use the Base Adjustment Lever to widen the Lifter Base. The Base must be spread to the widest possible position to maximize stability.
8. Pump the hydraulic Pump handle to lift the Stretcher and patient until the entire Stretcher is hanging from the Stretcher Hoops and is not touching the bed or any other object. Do not lock the Wheel on the Lifter Base. If needed, the Base will move slightly to adjust for balance.

Caution: Never move patient away from bed or leave patient unattended.

Taking Weight Measurements

1. The scale will calculate the patient's weight. Depending upon the movement of the stretcher, the scale may take several seconds to lock into the weight. We recommend reweighing after the patient is lying still on the stretcher by pressing the REWEIGH button.
2. After weight measurement is obtained and recorded, lower the Stretcher and the patient to the bed. Open the hydraulic pump control valve knob by turning it counterclockwise **very slowly** to control the descent.
3. Guide the patient and Stretcher down onto the bed.

OPERATION INSTRUCTIONS

4. Release the hanging loops from the lifter.

Note: The scale is programmed to save the last setting chosen as well as the weight of the stretcher as tare weight. This enables weighing a number of patients without having to re-set the scale.

If a weight was added or reduced from the stretcher scale it is required to reset the tare weight.

Note: To zero the scale by using the ZERO key, take the stretcher off the lifter before pressing ZERO key. Replace the stretcher after zeroing.

5. Press ON/OFF key to turn scale off.

OPERATION INSTRUCTIONS



Figure 7. Keypad

ITEM	DESCRIPTION	FUNCTION
	POWER	Turns scale ON and OFF.
	ZERO	Zeros the scale prior to weighing.
	HOLD/RELEASE	Holds the value of the weighed object on the display until the button is pressed again to clear the value. Also used to scroll down in the menu.
	REWEIGH	Allows repeated weighing of the patient without stepping off the scale.
	KG/LB button	Toggles between kilograms or pounds. Also used to scroll up in the menu.
	MENU	Enters the menu of the scale.
	BMI	Prompts entry of data to calculate the patient's Body Mass Index (BMI).
	PRINT	Prints patient's data (if printer is connected to the scale).
	ID	Prompts entry of patient's identification number (ID). This ID will be stored with all the weighing made until is cleared or a different ID is stored.
	TARE	Prompts entry of TARE value that will be deducted from the weight on the platform. Also releases tare weight (returns display to zero).
	EXIT	Reverts back one step when in the menu and data entry modes.
	ENTER	Used to enter commands and values into the scale.

OPERATION INSTRUCTIONS

MENU

In the menu screen the user can set preferences and/or instruct the scale how to handle stored data. The menu can be navigated using the up and down keys (▲▼) or by entering the associated menu position number with the keypad. The menu has a “roll-over” way of working: when the user scrolls to the bottom of the menu and presses the down button, it will return to the top of the menu.




Figure 8. Main Menu

01 WEIGHT DISPLAY MODE

(only applies to pound values, NOT the metric values)

The user can set the screen display value in either fraction of pound ($\frac{1}{4}$, $\frac{1}{2}$ or $\frac{3}{4}$ lb), in pounds and ounces (resolution of 4 ounces) or in decimals (resolution of 0.2 lb). When kg is selected as the units of use, these settings have no effect. The mode that is selected is used through all the screens.




Figure 9. Weight Display Mode

OPERATION INSTRUCTIONS

02 TARE WEIGHT

NOTE: Due to the scale's sensitivity, we recommend using the REWEIGH function prior to setting the TARE weight, in order to eliminate any operator interference with the item to be tared out.

There are two ways to manually enter a tare weight (such as the weight of a stretcher, shoes, etc): press MENU (1) and select option 02-Tare Weight or press TARE (9). If there is no weight on the scale platform (value displayed is zero and there is no tare) and the user presses the TARE button (9), the Tare Weight window appears and prompts the user to enter the TARE value and to press ENTER. NOTE: The tare weight must be entered using the following increments: 4oz, 0.2lb, 1/4lb. The value entered will display in minus (-). After the TARE has been entered, the scale goes back to normal operation. This TARE value is stored in memory until changed or cleared.

Important: The TARE weight cannot exceed 125 lb.




Figure 10. Tare Weight

Automatic Tare

The user can set a tare weight by pressing the TARE button (9) while there is a weight on the scale platform. The display will come to zero and TARE will be displayed to indicate there is a tare value in the memory (as displayed in Figure 10 above).

Removing the Tare

Additional pressing of the TARE button will delete the tare value from memory, TARE will disappear from the display and the scale will resume normal operation.

03 DATA MANAGEMENT

The scale manages patient data including patient ID, weight, height, tare and BMI. The value is stored in memory or transferred to PC. This function is performed by opening a new data file.

OPENING A NEW DATA FILE:

1. Press the ID button (7).
2. Using the keypad, type in the identification number.
3. Press ENTER.

This scale will also calculate Body Mass Index (BMI).


CALCULATING BMI:

1. Place the patient on the scale platform.
2. Press the BMI button (2).


OPERATION INSTRUCTIONS

- The scale will prompt you to enter the patient's height. Use the keypad to enter the height in 1/4" (1=1/4, 2=1/2, 3=3/4) or 1 cm increments (use inches for weight in lb, cm for kg). Press ENTER.
- Display will read the patient's BMI.

The scale offers you two options to manage your information: to transfer/download the values or to store them. The first option automatically downloads (transfers) the value to your PC. The second option stores the value in memory. The maximum capacity of the scale is 270 files of different data.




01 Auto Download




Automatic Download is the default option and will transfer the value to the PC as soon as the patient steps off the scale or when the user presses the HOLD/RELEASE button if it was kept in "HOLD". If no PC is connected, the value is not transferred and not stored and will be lost after the load is removed from the scale.


02 Store in memory



The value is stored in memory for later download to PC. If the memory is close to full the user will be warned and given the option to transfer all values to the PC or to clear the memory of all values.




03 Do not store data



All the data will be cleared.


OPERATION INSTRUCTIONS

04 Transfer now




All the values stored in memory are transferred to the PC and the scale memory is cleared of all values. If the transfer was unsuccessful, the values are kept in memory until successfully transferred or cleared.

05 Clear memory




All the values stored in memory will be cleared.

04 SCALE SETTINGS




01 Auto Hold Time



The user can determine how long to display the weight reading once it is determined, regardless of whether the patient remains on the platform. The scale defaults to no Auto Hold Time. The maximum setting is 20 seconds Hold Time.

02 Sleep Time



The user can set the time elapsed before the scale goes into the sleep mode. The default is 1 minute. When the scale goes into sleep mode, STANDBY is displayed on the screen.


OPERATION INSTRUCTIONS

03 Auto Off Time




The user can determine how long the scale will operate before turning off automatically due to inactivity. Default time is 10 minutes. If the value is set to zero, the auto off function is disabled.

04 Tone Volume




There is an option to adjust the beeping tone of the scale. This tone should sound when the scale has determined the weight on the platform, when a key is pressed, after the scale is turned on, at the end of self-test, or in the case of fault or warning.



Use the ▲ and ▼ keys on the keypad to adjust the volume. Whenever the user presses the key to change a volume, a beep will sound to indicate the set volume level.


05 Display Date and Time



This option will turn on or turn off the date and time display.


OPERATION INSTRUCTIONS

06 Display Backlight




The user can set the brightness of the backlight.

07 Display Contrast



The user can set the brightness of the LCD.


08 Live Weight



By selecting “Yes” the user can set the Live Weight mode to deactivate the motion-sensing mode. In the Live Weight mode the weight displayed will fluctuate with the patient’s movement; the scale will not lock on quickly to the weight as is the case in the motion-sensing mode.


Press the REWEIGH button to operate the motion-sensing mode and to determine the correct weight on the screen. To revert back to motion-sensing mode, change the Live Weight setting to “No”.

05 SYSTEM SETTINGS




01 Set Time & Date

The user can set the time and date using the keypad.




To set the time move between hours, minutes and seconds using the up and down keys ▲▼ and enter the values on the keypad. To jump to the AM/PM line press the ENTER button once.

OPERATION INSTRUCTIONS




Set the date using the up and down keys ▲▼ and enter the values on the keypad.

02 About




This screen displays the software version of the scale.

06 SYSTEM TEST



01 Battery Test



The scale will display the estimated amount of battery life remaining until the batteries will have to be replaced.

NOTE: In order to complete the battery test, the scale must be powered by batteries only. Unplug the scale from AC power source prior to battery test.


02 USB Connection




The scale will test the connection to the PC and will display a message "Connection is OK" or "NO Connection". If "NO Connection" is displayed, check your USB connections on the scale and on your PC and retest the connection. Refer to qualified service personnel if problem persists.

OPERATION INSTRUCTIONS


03 UI Test



The scale has a diagnostic routine where it tests the User Interface (UI) hardware functionality (LCD, keypad). In order to do this the user has to press all the keys according to the messages displayed on the screen.




If the requested command was not received or wrong button was pressed, the following message will be displayed.




If after 10 seconds the requested command was not received, the following message will be displayed. If "UI Error Failure" is displayed, refer to qualified service personnel.

07 USER SETTINGS



01 Retain Entered Values



This option allows the user to use the same values for ID, height and TARE between weighing. If this option is disabled, the user has to re-enter these values for each reading. If the values are not entered, only the weight is stored.

NOTE: These values cannot be retained by ID number.

OPERATION INSTRUCTIONS

02 Disp Height & ID

02-Disp Height & ID

Yes **No**

▲▼ to change Enter

01/14/2004 11:30pm
ZERO

5³/₄ lb

Height: 6ft00 0/0in
ID: 981723737123

01/14/2004 11:30pm
ZERO

5³/₄ lb

When the user selects to display the height and ID of the user, it will be displayed at the bottom of the screen. We recommend the use of this function to ensure that the patient's correct ID and height have been entered.

03 Prompt For Height

03-Prompt For Height

Yes **No**

▲▼ to change Enter

When this option is activated, the user will be asked to enter the patient's height after every weighing.

01/14/2004 11:30pm

153³/₄ lb

Please enter Height:
0ft00 0/0in

▲▼ to change Enter

04 Prompt For ID

04-Prompt For ID

Yes **No**

▲▼ to change Enter

When this option is activated, the user will be asked to enter the patient's ID number after every weighing.

01/14/2004 11:30pm

153³/₄ lb

Please enter ID:
2345872336

▲▼ to change Enter

MAINTENANCE

GENERAL

This section provides instructions for maintenance, cleaning, troubleshooting and operator replaceable parts for the ProPlus® Stretcher Scale Model 2000KL. Maintenance operations other than those described in this section should be performed by qualified service personnel.

MAINTENANCE

Before first use and after periods of non-use, check the scale for proper operation and function. If the scale does not operate correctly, refer to qualified service personnel.

1. Check overall appearance of the total scale for any obvious damage, wear and tear.
2. Inspect AC adapter for cord cracking or fraying or for broken or bent prongs.

CLEANING

Proper care and cleaning is essential to ensure a long life of accurate and effective operation.

Disconnect the scale from the AC power source.

1. Clean all external surfaces with a clean damp cloth or tissue. Mild soap and water solution may be used. Dry with a clean soft cloth.
2. Do not immerse the scale into cleaning or other liquid solution.
3. Do not use Isopropyl Alcohol or other solutions to clean the display surface.

TROUBLESHOOTING

Refer to the following instructions to check and correct any failure before contacting service personnel.


SYMPTOM	POSSIBLE CAUSE	CORRECTIVE ACTION
Scale does not turn on	<ol style="list-style-type: none"> 1. Dead Battery 2. Faulty electrical outlet 3. Bad power supply 	<ol style="list-style-type: none"> 1. Replace batteries 2. Use a different outlet 3. Replace adapter
Questionable weight or the scale does not zero	<ol style="list-style-type: none"> 1. External object interfering with the scale 2. The weight of the stretcher as tare was not set correctly 3. The patient's center mass is not at the center of the stretcher and one of the hanging hooks of the lifter boom does not lie on the lifter beam. 4. Scale is out of calibration 	<ol style="list-style-type: none"> 1. Remove interfering object from the scale 2. Take the patient and the stretcher off the scale, zero the scale, hang the stretcher on the scale and enter the tare weight. After the tare weight is set start weighing action. 3. Take the patient off the stretcher and re-position so that the center mass will lie between the two hoops. 4. Check weight with known weight value
Weighing is performed but the display shows "weigh" and "reweigh" every few seconds; the weighing process takes too long and no weight is displayed.	The patient is not still	Ask the patient to be still or change the scale to use the live weight setting (p.19).
The display shows "Overload" message	The load on the scale exceeds the capacity (400 Lbs)	Remove the excess weight and use the scale according to its limits
The display shows "LOW BAT" message	The batteries are empty	Replace batteries according to instructions
The display shows "Load Cell Error" message	There is a problem with one or more load cells or the load cell cable is disconnected.	Check load cell cable connection at the display and platform assembly ports. If the problem is not corrected, refer to qualified service personnel to replace the defective load cell

CALIBRATION PATH


The calibration is performed in Kg or Lb, according to the units of measure used upon entering into the calibration path.

Operator Action	Display
1. Press ON/OFF button and immediately press and hold HOLD/RELEASE button for 3-4 seconds	Calibration Process
2. Press ENTER	Enter load weight 400.0
3. Using ▲▼ keys and/or the keypad, set the calibration load to at least 200 lbs (for best results and greater accuracy, use a 400 lb weight). Press ENTER	Zero calibration Please clear the scale
4. Please clear the stretcher and press ENTER	Zero calibration Please wait
5. Please wait without touching or shaking the scale until the zero calibration process is finished	Weight calibration Put: XXX.X
6. Load the stretcher with the required weight and press ENTER	Weight calibration Please wait
7. Please wait without touching or shaking the scale until the calibration process is finished	Calibration factor: X.XXXXX
8. Press ENTER	Calibration Please clear the scale
9. Please remove weight from the stretcher	Rebooting Please wait
10. Please wait until the scale resumes normal operation	Health O Meter Pro Plus

EXPLODED VIEW OF SCALE



EXPLODED VIEW OF SCALE (continue)



PARTS LIST

Key No.	Part No.	Description	Qty.
1	1109403-0	DISPLAY ASSEMBLY	1
2	2268101-0	TILT MECH. ASSEMBLY	1
3	3848401-0	PLASTIC WASHER	2
4	401369	PLASTIC OVAL BAR COVER	3
5		HEX. HEAD SCREW NC 3/8-16*2 1/2	1
6		HEX. HEAD SCREW NC 3/8-16*1 3/4	1
7	2268001-0	LIFTER BOOM ASSEMBLY	1
8	2034701-0	MAIN BOARD TO CON. DIG CABLE	1
9	3809801-0	BEAM ASSEMBLY	1
10		PHIL.PAN HEAD SCREW NC6-32*1/4"	2
11	2138903-0	CON. DIG BOARD	1
12		SPACER WITH EX.THREAD STUD NC6*1/4"	2
13		SOCKET BUTTON HEAD CAP SCREW NC1/4-20*3/4"	2
14		SAFETY WASHER NC3/8"	2
15	401304	HEX. SOCKET HEAD SCREW NC3/8-16*1 1/2	2
16	3808401-0	STRETCHER HOOP	2
17	401338	STRETCHER	1
18	3835701-0	STRETCHER STORAGE TRAY	1
19	3833801-0	BAR ADAPTER	2
20	3808501-0	STRETCHER BAR	4
21	3816901-0	FORK JOINT	2
22		HEX. SOCKET HEAD SHOULDER SCREW (5/16) 3/8-1"	4
23	3816801-0	LOAD CELL ADAPTER	2
24	3828301-0	LOAD CELL CONNECTING BEAM	1
25	1304202-0	LOAD CELL	2
26		HEX. SOCKET HEAD CAP SCREW NC5/16-18*3/4"	4
27		NUT NC3/8"	2
28	401337	PLASTIC RHS50*50 COVER	2
29	400579	CABLE HOLDER P.G.-7	1
30	400066	SELF LOCK NUT NC3/8"	2
31		PHIL. PAN HEAD SCREW NC8*3/8"	6
32	3822701-0	CABLE COVER	1
33	400412	SELF LOCK NUT NC6	2
34		PHIL. PAN HEAD SCREW NC6-32*3/8"	2
35	3817001-0	ADAPTER HOLDER	1
36	420988	ADAPTER 6VDC 120VAC	1
37	411220	KEYPAD 400LB	1
38		INDICATOR COVER	1
39		LCD WINDOW	1
40		LCD BOARD	1
41		WN1412 CROSS HEAD SCREW K22L6	4
42		DISPLAY TO EAGLE P.C.B CABLE	1
43		WN1412 CROSS HEAD SCREW K30L6	4
44		EAGLE BOARD	1
45		INDICATOR BASE	1
46		PHIL. PAN HEAD SCREW NC6-32*3/8"	6
47		MODEL LABEL	1
48	3245801-0	INPUT/OUTPUT PORTS LABEL	1
49		PHIL. PAN HEAD SCREW NC4-40*1/4"	1
50	3822801-0	BATTERY COVER	1
51	400152	RUBBER O RING	1
52	2033801-0	BATTERY HOLDER	1
---	63855	PATIENT LIFT	1

WARRANTY

LIMITED WARRANTY

What does the Warranty Cover?

Pelstar LLC scales are warranted from date of purchase against defects of materials or in workmanship for a period of three (3) years. If product fails to function properly, return the product, freight prepaid and properly packed to Pelstar. See "To Get Warranty Service" below for instructions. If manufacturer determines that a defect of material or in workmanship exists, customers' sole remedy will be repair or replacement of scale at no charge. Replacement will be made with a new or remanufactured product or component. If the product is no longer available, replacement may be made with a similar product of equal or greater value. All parts including repaired and replaced parts are covered only for the original warranty period.

Who is Covered?

The original purchaser of the product must have proof of purchase to receive warranty service. Pelstar dealers or retail stores selling Pelstar products do not have the right to alter, or modify or any way change the terms and conditions of this warranty.

What is Excluded?

Your warranty does not cover normal wear of parts or damage resulting from any of the following: negligent use or misuse of the product, use on improper voltage or current, use contrary to the operating instruction, abuse including tampering, damage in transit, or unauthorized repair or alternations. Further, the warranty does not cover Acts of God, such as fire, flood, hurricanes and tornadoes. This warranty gives you specific legal rights, and you may also have other rights that vary from country to country, state to state, province to province or jurisdiction to jurisdiction.

To get Warranty Service

Make sure you keep your sales receipt or document showing proof of purchase.

Call 1 (800) 638-3722 or 1 (708) 598-9100 to receive a return authorization number. Attach proof of purchase to your defective product along with your name, address, daytime telephone number and description of the problem. Carefully package the product and send with shipping and insurance prepaid to:

Pelstar LLC
Attention R/A# _____
Repair Department
7400 W. 100th Place
Bridgeview, IL 60455

If your scale is not covered by warranty, or has been damaged, an estimate of repair costs or replacement costs will be provided to you for approval prior to servicing or replacing.

Pelstar LLC

7400 West 100th Place, Bridgeview IL 60455 • 1-800-638-3722 or 1-708-598-9100

www.healthometermedical.com

Health o meter ProPlus® products are manufactured, designed and owned by Pelstar LLC.
Health o meter® is a registered trademark of Sunbeam Products Inc., Boca Raton, FL 33431.
ProPlus® is a trademark of Pelstar LLC

Patents Pending

Free Manuals Download Website

<http://myh66.com>

<http://usermanuals.us>

<http://www.somanuals.com>

<http://www.4manuals.cc>

<http://www.manual-lib.com>

<http://www.404manual.com>

<http://www.luxmanual.com>

<http://aubethermostatmanual.com>

Golf course search by state

<http://golfingnear.com>

Email search by domain

<http://emailbydomain.com>

Auto manuals search

<http://auto.somanuals.com>

TV manuals search

<http://tv.somanuals.com>