TANITA

Monitoring Your Health

GUIDE

InnerScan™ *Body Composition Monitor*

Model : BC-550

- GB** A guide to your body composition monitor
- D** Ein Handbuch für Ihre Körperanalysewaage
- F** Guide de votre impédancemètre
- NL** Een gids bij uw lichaamssamenstellingsmonitor
- I** Una guida al vostro monitor della composizione corporea
- E** Guía para el uso de su monitor de composición corporal
- P** Um guia para o seu monitor de composição corporal

Tanita UK Ltd.

The Barn, Philpots Close, Yiewsley,
Middlesex, UB7 7RY UK
Tel: +44 (0)1895 438577
Fax: +44 (0)1895 438511
www.tanita.co.uk

Tanita Corporation

14-2, 1-Chome, Maeno-Cho
Itabashi-Ku, Tokyo, Japan 174-8630
Tel: +81(0)-3-3968-2123
Fax: +81(0)-3-3967-3766
www.tanita.co.jp
ISO 9001 Certified

Tanita Europe GmbH

Dresdener Strasse 25
D-71065 Sindelfingen
Germany
Tel: +49 (0)-7031-6189-6
Fax: +49 (0)-7031-6189-71
www.tanita.de

Tanita Corporation of America Inc.

2625 South Clearbrook Drive
Arlington Heights, IL 60005 USA
Tel: +1-847-640-9241
Fax: +1-847-640-9261
www.tanita.com

Tanita France S.A.

Villa Labrouste
68 Boulevard Bourdon
92200 Neuilly-Sur-Seine France
Tel: +33 (0)-1-55-24-99-99
Fax: +33(0)-1-55-24-98-68
www.tanita.fr

Tanita Health Equipment H.K. Ltd.

Unit 301-303 3/F Wing On Plaza
62 Mody Road, Tsimshatsui East
Kowloon, Hong Kong
Tel: +852 2838 7111
Fax: +852 2838 8667

©2005 Tanita Corporation

BC5507901

Download from www.Somanuals.com. All Manuals Search And Download.

WHO CAN USE A BODY COMPOSITION MONITOR

This Body Composition Monitor is intended for adults aged 18 – 99 years. Children aged 7 – 17 years can use the monitor for weight and body fat percentage readings only; the other features are not applicable to children.

Your monitor is also equipped with an athlete mode for adults with athletic body types. Tanita defines an athlete as a person involved in intense physical activity of approximately 10 hours per week and has a resting heart rate of 60 beats per minute or less. Individuals who have been an athlete for a number of years but currently exercises less than 10 hours per week can also use athlete mode.

Pregnant women should use the weight function only.
All other functions are not intended for pregnant women.

The Tanita Body Composition Monitor is intended for home use only. This unit is not intended for professional use including hospitals, medical or fitness facilities; it is not equipped with the adequate specifications required for heavy usage. Using the unit in a professional location will invalidate your warranty.

The Tanita InnerScan Body Composition Monitor product provides readings for informational purposes only. This product is not intended to diagnose or treat any disease or abnormalities. Please consult with your physician with any questions or concerns related to your health.

WHY MONITOR BODY COMPOSITION?

Body composition monitors are designed for healthier living by giving you an insight into key health indicators that will enable you to monitor the impact of changes to your lifestyle:

- See the impact of a change in diet on the body to make sure you are dieting the healthy way
- Fine-tune your fitness programs by monitoring progress of muscle mass and Basal Metabolic Rate
- Monitor visceral fat, a key indicator of the risk of Type 2 diabetes and heart disease.
- Set a target for your physique and monitor your progress towards it.

HOW DOES A BODY COMPOSITION MONITOR WORK?

Tanita Body Composition Monitors calculate your body composition using Bio-electrical Impedance Analysis (BIA). Safe, low-level electrical signals are passed through the body via the patented Tanita footpads on the monitor platform. It is easy for the signal to flow through fluids in the muscle and other body tissues but meets resistance as it passes through body fat, as it contains little fluid. This resistance is called impedance. The impedance readings are then entered into medically researched mathematical formulas to calculate your body composition.

WHEN IS THE BEST TIME TO USE MY BODY COMPOSITION MONITOR?

Your body water levels naturally fluctuate throughout the day and night. Any significant changes in body water may affect your body composition readings; for example, the body tends to be dehydrated after a long night sleep so if you take a reading first thing in the morning your weight will be lower and your body fat percentage higher. Eating large meals, drinking alcohol, menstruation, illness, exercising, and bathing may also cause variations in your hydration levels.


To get the most reliable reading it is important to use your Body Composition Monitor at a consistent time of day under consistent conditions. We suggest taking a reading before your evening meal.

WHAT IS BODY FAT PERCENTAGE?

Body fat percentage is the amount of body fat as a proportion of your body weight.

Reducing excess levels of body fat has shown to reduce the risk of certain conditions such as high blood pressure, heart disease, diabetes and cancer. The chart below shows the healthy ranges for body fat.

¹ Based on NIH/WHO BMI Guidelines.

² As reported by Gallagher, et al. at NY Obesity Research Center.

To determine the percentage of body fat that is appropriate for your body, consult your physician.

Healthy Range Indicator (Standard Adults only*)

Your Body Composition Monitor automatically compares your body fat percentage reading to the Healthy Body Fat Range chart. After your body fat percentage has been calculated, a black bar will flash along the bottom of the display, identifying where you fall within the Body Fat Ranges for your age and gender.

- (-) : Underfat; below the healthy body fat range. Increased risk for health problems.
- (0) : Healthy; within the healthy body fat percentage range for your age/gender.
- (+) : Overfat; above the healthy range. Increased risk for health problems.
- (++) : Obese; high above the healthy body fat range.
Greatly increased risk of obesity-related health problems.

*Note :

If you select Athlete mode, the unit will not display the Healthy Range Indicator. Athletes may have a lower body fat range depending on their particular sport or activity. (Range of user age is between 18-99)

WHAT IS TOTAL BODY WATER PERCENTAGE?

Total Body Water Percentage is the total amount of fluid in a person's body expressed as a percentage of their total weight. Water plays a vital role in many of the body's processes and is found in every cell, tissue and organ. Maintaining a healthy total body water percentage will ensure the body functions efficiently and will reduce the risk of developing associated health problems.

Your body water levels naturally fluctuate throughout the day and night. Your body tends to be dehydrated after a long night and there are differences in fluid distribution between day and night. Eating large meals, drinking alcohol, menstruation, illness, exercising, and bathing may cause variations in your hydration levels.

Your body water percentage reading should act as a guide and should not be used to specifically determine your absolute recommended total body water percentage. It is important to look for long-term changes in total body water percentage and maintain a consistent, healthy total body water percentage.

Drinking a large quantity of water in one sitting will not instantly change your water level. In fact, it will increase your body fat reading due to the additional weight gain. Please monitor all readings over time to track the relative change.

Every individual varies but as a guide the average total body water percentage ranges for a healthy adult are:

Female : 45 to 60%

Male : 50 to 65%

Source : Based on Tanita's Internal Research

Note: The total body water percentage will tend to decrease as the percentage of body fat increases. A person with a high percentage of body fat may fall below the average body water percentage. As you lose body fat the total body water percentage should gradually move towards the typical range given above.

WHAT IS DAILY CALORIE INTAKE(DCI)?

"Daily Calorie Intake (DCI)" is the sum of calories for basal metabolism, daily activity metabolism (activities including daily household chores), and diet-induced thermogenesis (energy used in connection with digestion, absorption, metabolism, and other eating activities). It is an estimate of how many calories you can consume within the next 24 hours to maintain your current weight.

HOW DOES A TANITA BODY COMPOSITION MONITOR CALCULATE DCI?

The basic way of calculating Basal Metabolic (BMR) Rate is a standard equation using weight and age. Tanita has conducted in-depth research into the relationship of BMR and body composition giving a much more accurate and personalized reading for the user based on the impedance measurement. This method has been medically validated using indirect calorimetry (measuring the breath composition).*

*Reliability on equation for Basal Metabolic Rate: At 2002 Nutrition Week : A Scientific and Clinical Forum and Exposition Title: International Comparison: Resting Energy Expenditure Prediction Models: The American Journal of Clinical Nutrition.

DCI=BMR × Activity Level
Activity Level

	1	2	3
Female	1.56	1.64	1.82
Male	1.55	1.78	2.10

Source: World Health Organization(WHO)

WHAT IS BONE MASS?

This feature indicates the amount of bone (bone mineral level, calcium or other minerals) in the body. Research has shown that exercise and the development of muscle tissue are related to stronger, healthier bones. While bone structure is unlikely to make noticeable changes in a short period, it is important that you develop and maintain healthy bones by having a balanced diet and plenty of exercise. People worried about bone disease should consult their physician. People who suffer from osteoporosis or low bone densities due to advanced age, young age, pregnancy, hormonal treatment or other causes, may not get accurate estimations of their bone mass.

Below is the result of estimated bone masses of persons aged 20 to 40, who are said to have the largest amounts of bone masses, by weight. (Source:Tanita Body Weight Science Institute)

Please use the below charts as a guide to compare your bone mass reading.

Women : Average of estimated bone mass(lb)

Weight		
Less than 110 lb	110lb - 165 lb	165 lb and up
4.3 lb	5.3 lb	6.5 lb

Men : Average of estimated bone mass(lb)

Weight		
Less than 143 lb	143 lb - 209 lb	209 lb and up
5.9 lb	7.3 lb	8.1 lb

Note: Persons described below may obtain varying readings and should take the values given for reference purposes only.
- Elderly persons
- Women during or after menopause
- People receiving hormone therapy

Note: "Estimated bone mass" is a value estimated statistically based on its correlation with the fat-free amount (tissues other than the fat). "Estimated bone mass" does not give a direct judgment on the hardness or strength of the bones or the risks of bone fractures. If you have concerns over your bones, you are recommended to consult a specialist physician.

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