

6RUNNER® 14 (888WNLL, 888WNLLHD)

SHOPRIDER®



The 6Runner® series powerchairs, all with 6 points of contact with the ground, provide the desired stability of traditional rear driven powerchairs, along with the maneuverability of modern mid wheel driven powerchairs. The 6Runner® series powerchairs deliver the best features available in today's powerchair technology all in one finished package. Tomorrow's products today, only from Shoprider.

Available Colors



Red



14" Wheel

6RUNNER® 14 (888WNLL,888WNLLHD)

- Luxurious leatherette seat
- 21" high back executive seat
- Stable four post seating system
- Reclining seat and adjustable seat height
- Adjustable length, height and angle foot plate
- Adjustable width, height and flip-up arm rest
- Midwheel drive maneuverability
- 14" wheels
- 400 lb. capacity (888WNLL)
450 lb. capacity (888WNLLHD)



Adjustable height



Many seating options available



Flip-up arm rest



Reclining seat

6RUNNER® 14 SPECIFICATIONS

		(888WNLL)	(888WNLLHD)
Overall Dimensions (L* ¹ x W x H)	in	41x 25 x 46	41x 25 x 46
Number, Size of Tires (front, middle, rear)	ea, in	2, 7"; 2, 14"; 2, 7"	2, 7"; 2, 14"; 2, 7"
Max. Single User Weight (on the level road)	lb	400	450
Weight of the heaviest part (lower frame)	lb	82	82
Total Unit Weight (w/battery, 75Ah*2)	lb	277	277
Each Battery Weight (75Ah)	lb	57	57
Battery Module Capacity	_V_Ah x Pcs	12V75Ah x 2	12V75Ah x 2
Charger	Amp	8 amp off board	8 amp off board
Output Power of Motor (for Reference only)	hp x Pcs	1.34 x 2	1.34 x 2
Turning Radius	in	25	25
Max. Speed* ²	mph	3.75	3.75
Range* ⁴ (Per Charge w/STD 75Ah battery/300lbs. User Weight) (After the battery and mechanical moving parts fully break in)	mile	26	20

All specification subject to change without prior notice.

WARNING! Please make sure that you are physically fit to operate this unit. Consult your dealer /doctor for advice.

*1. Include the Anti-Tip Wheel.

*2. Driver weight may exceed weight of the unit, speed must be reduced when turning.

*3. When proceeding up any incline, please move your seat to the most forward position or if you have the deluxe seat, make sure that it is in the 90 degree (upright) position.

*4. The actual driving range varies with the factors shown as below.

(A) The Weight of Occupant (E) Ambient Temperature
(B) Level Ground (F) The Way of Driving
(C) Battery Conditions (G) ETC.
(D) Type of Charger

*5. The Centers for Medicare and Medicaid Service (CMS).
The healthcare Common Procedure Coding System (HCPCS)

All specifications are subject to change without prior notice.

Shoprider Mobility Products, Inc. reserves the rights of any changes to the unit.

1) Include the anti-tip wheel or the rear castor.

2) Driver weight may exceed weight of the unit; speed must be reduced when turning.

3) The actual driving range varies with the factors shown as below:

a) The weight of occupant
b) Ground surface
c) Battery capacity and conditions
d) Type of charger
e) Ambient temperature
f) The way of driving
g) If the battery and mechanical moving parts fully break in
h) Etc.

SHOPRIDER

Toll-Free: 800.743.0772

www.shoprider.com

Free Manuals Download Website

<http://myh66.com>

<http://usermanuals.us>

<http://www.somanuals.com>

<http://www.4manuals.cc>

<http://www.manual-lib.com>

<http://www.404manual.com>

<http://www.luxmanual.com>

<http://aubethermostatmanual.com>

Golf course search by state

<http://golfingnear.com>

Email search by domain

<http://emailbydomain.com>

Auto manuals search

<http://auto.somanuals.com>

TV manuals search

<http://tv.somanuals.com>